

Operating Instructions

METTLER Analytical balance AE240 Dual Range Balance

METTLER TOLEDO



Leveling the balance  1. Turn the balance on. 2. Turn the leveling feet clockwise until the bubble is centered.	Weighing range selection  1. Turn the balance on. 2. Press the Range key. 3. Turn the selection wheel to choose the weighing range: - 0.1 mg (0.0001 g) - 0.01 mg (0.00001 g) - 0.001 mg (0.000001 g)	Calibrating the balance  1. Turn the balance on. 2. Turn the calibration adjustment knobs clockwise until the digital display shows "CAL".	Integration time Steps 1/2/3  1. Turn the balance on. 2. Press the Integ key. 3. Turn the selection wheel to choose the integration time: - 1 ms - 1 s - 10 s - 100 s	Stability detector Steps 1/2/4/Off  1. Turn the balance on. 2. Press the Stab key. 3. Turn the selection wheel to choose the stability detector: - Off - 1 ms - 1 s - 10 s - 100 s
Switching the display on  2. Turn the balance on. 3. Turn the display on/off switch clockwise. Digital display: - 0.0000 - 0.00000 - 0.000000 - 0.0000000	 2. Turn the balance on. 3. Turn the display on/off switch clockwise. Digital display: - TWS 200 - TWS 40 - 0.00000 - 0.000000	 2. Turn the balance on. 3. Turn the calibration adjustment knobs clockwise until the digital display shows: - CAL - 0.00000 - 0.000000 - 0.0000000	 2. Turn the balance on. 3. Turn the selection wheel to choose the integration time: - 1 ms - 1 s - 10 s - 100 s - TWS - TWS 200 - TWS 40 - 0.00000 - 0.000000 - 0.0000000	 2. Turn the balance on. 3. Turn the selection wheel to choose the stability detector: - Off - 1 ms - 1 s - 10 s - 100 s - TWS - TWS 200 - TWS 40 - 0.00000 - 0.000000 - 0.0000000
Taring  1. Turn the balance on. 2. Place a reference weight on the platform. 3. Press the Tare key. Digital display: - 0.0000 - 0.00000 - 0.000000 - 0.0000000	Weighing  1. Turn the balance on. 2. Place the sample on the platform. 3. Press the Start key. Digital display: - 0.0000 - 0.00000 - 0.000000 - 0.0000000	 4. Turn the balance on. 5. Place the sample on the platform. 6. Press the Start key. Digital display: - 0.0000 - 0.00000 - 0.000000 - 0.0000000	Short-term operating instructions  1. Turn the balance on. 2. Place the sample on the platform. 3. Press the Start key. Digital display: - 0.0000 - 0.00000 - 0.000000 - 0.0000000	

Mettler-Toledo AG, CH-8401 Greifensee, Switzerland,
Tel. (01) 944-22 00, Fax (01) 944-30 80, Telex 828150

A. Mettler-Toledo Gmbh & Co., Packerer Gewerbegebiet, D-7100 Stuttgart,
Tel. (0711) 502-01 00, Fax (0711) 502-00 00

B. Mettler-Toledo Ltd., Pines Way, Wetherby, LS23 7EE, UK,
Tel. (0193) 234-00 01, Fax (0193) 234-11 88

CH. Mettler-Toledo (Champlain) Ltd., 4th Floor, Building A, 4-471 Get Pregelstr,
Urbino 2000235, Tel. (011) 470 3332, Fax (011) 470 3331

CH. Mettler-Toledo (Zürich) AG, Zürichstrasse, CH-8006 Zürich, Switzerland,
Tel. (01) 944-45 45, Fax (01) 944-45 30

D. Mettler-Toledo GmbH, Postfach 100040, D-3300 Kassel, Germany,
Tel. (05241) 5027-0, Fax (05241) 5027-15

DK. Mettler-Toledo A/S, Søborgvej 172, DK-2100 Copenhagen, Denmark,
Tel. (03) 25 00 00, Fax (03) 25 04 29

E. Mettler-Toledo S.A.E., Paseo de la Castellana 77, E-28040 Madrid, Spain,
Tel. (03) 22 22 22, Fax (03) 22 22 02 71

F. Mettler-Toledo S.p.A., Via F. 14, 16-20, via de' Pepparini, I-20123 Milano, Italy,
Tel. (02) 30 17 17 17, Fax (02) 30 17 16 00

G. Mettler-Toledo S.p.A., Via Ugo, 42, I-20020 Novate Milanese,
Tel. (02) 22 02 02, Fax (02) 22 02 02

H. Mettler-Toledo S.p.A., Via Ugo, 42, I-20020 Novate Milanese,
Tel. (02) 22 02 02, Fax (02) 22 02 02

I. Mettler-Toledo AB, Västgötagatan 12, S-113 87 Stockholm, Sweden,
Tel. (08) 522 00 00, Fax (08) 524 40 02

JK. Mettler-Toledo Inc., 200 Ayer Street, Somerville, MA 02145, USA,
Tel. (617) 652 0773, Fax (617) 652 0773

TH. Mettler-Toledo (Thailand) Division, CH-Germany, 17/236 Rama 11, P.O. Box 1000,
Nonthaburi 10110, Tel. (031) 271 00 00, Fax (031) 271 00 79

UK. Mettler-Toledo Ltd., 14 Weston Road, Wetherby, LS23 7EE, UK,
Tel. (0193) 234-00 01, Fax (0193) 234-00 00


US. Mettler-Toledo Inc., 500 Morris Avenue, PO Box 21, Hightstown,
NJ 08520-0001, USA, Tel. (609) 441-3000, Fax (609) 441-0425

For all other countries: Analyse (France), AG, PO Box VI-100, CH-1400 Gossau,
Tel. (01) 844-22 01, Fax (01) 844-30 00

METTLER AE 240 Dual Range Balance

Operating elements and connections

- 1 Brief Operating Instructions
- 2 Single control bar
- 3 Level indicator
- 4 Calibration lever
- 5 Weighing pan / windshiedring
- 6 Leveling screw



- 10 Voltage selector
- 11 Microfuse
- 12 Power line connection socket
- 13 Socket for handkey / foot pedal



Preliminary steps

Checking the operating voltage

The operating voltage setting must agree with your local power-line voltage; please check this setting and, if needed, change it.

Admissible power-line voltage in the switch positions:
115 V...92 V...132 V
220 V...184 V...265 V

Location

- A stable location, as free from vibration as possible.
- Make sure there are no large temperature fluctuations.
- Avoid direct sunlight and drafts.
- Connect the power cable of the work station.

Installing the weighing pan and the windshiedring; leveling the balance

- Place weighing pan on balance; the control peg centers the pan in the opening in the base of the weighing chamber.
- The two leveling screws should be adjusted so that the bubble is in the middle of the circle.

Whenever the location of the balance is changed, the balance should be leveled.

Operation

Short-form operating instructions

Short-form operating instructions can be found on a card that swings out from underneath the balance housing.

Switching the display on/off

- Briefly press the single control bar; all display segments light up for several seconds.

8.8.8.8.8.8.8.

- Afterwards, the display automatically sets itself to zero.
- Lightly lift the control bar; the display is switched off.

Weighing range selection

- Press control bar until "mg" appears (weighing range).
- Releasing the control bar and pressing it again briefly, permits to alternate between the 40 g and the 200 g ranges.

After the desired range has been selected wait for "Display" — followed by zero. Balance is now in weighing mode.

Calibration

Make absolutely sure:

The balance must be left connected to the power supply for at least 60 minutes before "calibrating". For calibration it is of no consequence whichever weighing range has been selected. Calibration is being performed for both weighing ranges simultaneously.

- Press and hold the single control bar until "CAL" appears in the display, then release control bar. The display changes to "CAL" — then to CAL 100 (calibration).
- Move calibration lever all the way to the rear; the display changes to CAL — followed by 100.000, then to CAL 0 (pink).
- Move calibration lever all the way back towards the front of the balance; the display changes to — followed by zero.

Note: After selecting the weighing range, "calibration", "integration time" or "stability check" can be selected by extended pressing of the control bar.

Measuring cycle / measuring accuracy

By selecting a particular integration cycle, as well as a particular stability detection step, the balance can be configured according to your weighing location and needs.

Integration time

Step 1: Used for very stable/vibration-free weighing table (short measuring cycle).

Step 2: Normal setting.

Step 3: Used for unfavorable ambient conditions (long measuring cycle).

- Press the control bar and hold it until "int" appears in the display, then release the control bar.
- Immediately press the control bar briefly; the display will change to the next step.
- Step of the step you wish to use and wait for the display to return to the weighing mode (zero).

Stability detector

Step 1: Great sensitivity (long pause before data are released).

Step 2: Less sensitivity (short pause before data are released).

Normal setting.

off: The stability detector is switched off. Please note that even this is the case, DeltaDisplay is also switched off (described in Paragraph entitled "Weighing-in").

- Press the control bar and hold until "off" appears in the display, then release control bar.
- Immediately press control bar again briefly; the display changes to the next step.
- Step of the step you wish to use and wait for the display returns to the weighing mode (zero).

Note: After selecting the integration time, you can go directly to the selection of the stability detector setting by holding the control bar down.

Taring

- Open the sliding glass door.
- Place a tare container on the weighing pan.
- Close the sliding glass door.
- Press the control bar briefly; the display changes to zero.

Note: If it is possible to carry out external taring by using the handkey or foot pedal from the "accessories optional" (connection sockets on rear of balance).

The weight of the container is now tared out. To weigh in, the balance weighing range — minus the weight of the tare container — is now available.

Weighing-in (Aks 1 or 2)

- Open the sliding glass door.
- Fill in substance up to the desired target weight (to read the weight accurately, the door must be closed).

If different components are to be weighed, one after the other, into the same container, it is possible to tare after each weighing and start the next weighing from zero. This can be done until the tare container and all the components together reach the end of the weighing range.

DeltaDisplay:



The DeltaDisplay switches on automatically when weighing in substances quickly; the last two digits are blocked out and the display change sequence speeds up. This allows the increase in weight to be followed better. When weighing in slowly towards the target weight, the two digits switch back on. The display change sequence remains repeatable. Only when weight changes are very small does the balance switch back to the normal display change with the full number of decimal places.

Stability detector:


When stability is achieved (determined by the step selected for the stability detector), the green dot in the display goes out.

The result is then stable.

Note: When the green dot lights up in the display, the data interface is blocked; when the green dot goes off (stability), the data interface is unlocked.

Further capabilities

- 20 Cover
- 21 Screw
- 22 Hook



GD hanger (for weighing below the balance)

- Open all sliding glass doors.
- Remove the weighing pan.
- Place the balance on its back.
- Loosen the screw on the bottom of the balance.
- Swing the cover to one side.
- Refreshen the screw.

A hook is visible in the opening; the object or substance can be weighed by attaching a hanger from this hook.

- Place the balance back on its feet, place the weighing pan back on and level the balance.

- With the hanger attached to the hook, press here.

Note: The weighing pan does not have to be placed back on if the hanger is of least as heavy as the pan. The hanger is not available from METTLER.

Care and maintenance

Cleaning

A cloth with some soapy water is sufficient to clean the weighing pan and housing. Do not use any strong solvents. To remove residues from the weighing chamber, use the small metal brush that is included in the balance standard equipment (do not blow air through the chamber).

Replacing the microfuse

- Disconnect the power cable.
- Turn out the fuse holder (in the power-line connection socket) with a screwdriver.
- Replace the fuse (captive fuse-in the fuse holder).
- Push the fuse holder back on.
- Plug the power-line cable back in.

Accessories

Optional equipment

	Order No.
- Windshield ring, can be stocked, 1 unit	38598
- Tweezers, 210 mm long (with plastic tips)	70209
- Density (specific gravity) determination kit	33340
- Foot pedal	48278
- Handkey	42500
- Microfuses, 160 mA slow-blowing (set of 3)	55144
- Data interfaces:	
011 Option - CL/RB232C unidirectional	38750
012 Option - CL/RB232D bidirectional	38751
013 Option - IEEE 488	38752
040 Opt. Output (unidirectional mode)	38755

Standard equipment

- Power-line cable	neutral	80578
	Switzerland	87920
	Germany	87925
	USA	86668
- Weighing pan, 80 mm dia.	38590	
- Centering disk (for windshield ring)	38929	
- Hair-bristle brush	70114	
- Windshield ring	38680	

Specifications

AE240 Dual Range Balance	
40 g range	200 g range
Repeatability	0.01 mg
Weighting range	0...41 g
Tare range (subtractive)	0...41 g
Reproducibility (standard deviation)	0.02 mg
Linearity	± 0.03 mg
Stabilization time (typical)	8 sec
Integration time (adjustable)	3/6/12 sec
Display sequence	- Mettler DeltaDisplay off - Mettler DeltaDisplay on
Stability detector:	0.4 sec: 0.2/0.4 sec
- Sensitivity selectable in these steps:	1/2/ off
Sensitivity drift (10...30°C)	± 2 x 10⁻⁴ /°C
Calibration weight (built-in), adjusted to an apparent mass of 0.0 g/cm³ at an air density of 1200 mg/l	100 µg, adjusted to ± 0.1 mg
Dimensions:	Weighting pan (stainless steel) Open space above weighing pan Balance housing (W x D x H) Net weight
Power supply:	115 V/220 V 92...132 V, 184...265 V 50...60 Hz Power consumption
Admissible ambient conditions during operation	10...40°C 25...85%
Temperature	10...40°C
Relative humidity (non-condensing)	25...85%
What's wrong if...	then...
- the entire display does not light up?	- no power reaching the instrument. - the fuse is defective.
- the OFF display appears?	- a temporary power failure has taken place. (Press the control bar.)
- only the upper horizontal segments light up in the display?	- the weighing range has been exceeded. - the calibration weight has been activated. - there was weight on the pan when the instrument was switched on.
- only the lower horizontal segments light up in the display?	- the weighing pan is not installed. - there was weight on the pan when the instrument was switched on.
- the weighing result is unreliable?	- there are too many drafts. - the weighing table is unstable. - the integration time setting is too low. - the object being weighed is not at room temperature. - the balance must be calibrated or has been calibrated using the wrong external weight.
- the weighing result is obviously incorrect?	- a temporary malfunction has occurred (pull out power cord and plug it back in).
- only a portion of the display lights up?	- the weighing table or the load is too unsteady (close sliding glass doors, set a longer integration time and/or change the stability detection setting).
- the middle horizontal segments in the display are blurring (for more than 30 sec)?	- the weighing pan was not unloaded before calibrating the balance, or the wrong internal calibration weight was used (return to the weighing mode by pressing and holding the control bar).
- CAL Err appears in the display?	- a temporary malfunction has occurred (recalibrate balance).
- no CAL appears in the display?	- the weighing table or the load is too unsteady (close sliding glass doors, set a longer integration time and/or change the stability detection setting).
- a zero display does not appear after pressing zero?	- the weighing table or the load is too unsteady (close sliding glass doors, set a longer integration time and/or change the stability detection setting).