### GeneAmp<sup>®</sup> PCR System 9700 Base Module

User's Manual



#### © Copyright 2001, Applied Biosystems. All rights reserved.

#### For Research Use Only. Not for use in diagnostic procedures.

Information in this document is subject to change without notice. Applied Biosystems assumes no responsibility for any errors that may appear in this document. This document is believed to be complete and accurate at the time of publication. In no event shall Applied Biosystems be liable for incidental, special, multiple, or consequential damages in connection with or arising from the use of this document.

#### Authorized Thermal Cycler

This GeneAmp\* PCR System 9700 Base Unit, Serial No\_\_\_\_\_\_, in combination with its immediately attached sample block modules, comprise an Authorized Thermal Cycler. The purchase price of this Base Unit includes the up-front fee component of a license under United States Patent Nos. 4,683,195, 4,683,202 and 4,965,188, owned by Roche Molecular Systems, Inc., and under corresponding claims in patents outside the United States, owned by F. Hoffmann-La Roche Ltd, covering the Polymerase Chain Reaction ("PCR") process to practice the PCR process for internal research and development using this instrument. The running royalty component of that license may be purchased from Applied Biosystems or obtained by purchasing Authorized Reagents. This instrument is also an Authorized Thermal Cycler for use with applications licenses available from Applied Biosystems. Its use with Authorized Reagents also provides a limited PCR license in accordance with the label rights accompanying such reagents. Purchase of this product does not itself convey to the purchaser a complete license or right to perform the PCR process. Further information on purchasing licenses to practice the PCR process may be obtained by contacting the Director of Licensing at Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404.

ABI PRISM and its Design, Applied Biosystems, and MicroAmp are registered trademarks of Applera Corporation or its subsidiaries in the U.S. and certain other countries.

ABI, AB (Design), Applera, and BigDye are trademarks of Applera Corporation or its subsidiaries in the U.S. and certain other countries.

AmpErase, AmpliTaq, AmpliTaq Gold, EnviroAmp, GeneAmp, and TaqMan are registered trademarks of Roche Molecular Systems, Inc.

AppleScript and Macintosh are registered trademarks of Apple, Inc.

Microsoft, Windows, and Windows NT are trademarks of Microsoft Corporation.

All other trademarks are the sole property of their respective owners.

Printed in Singapore 11/2001 Part Number 4303481 Rev. E

#### *Contents*

#### 1 Introduction and Safety

Overview
About This Chapter 1-1
In This Chapter
About This Manual
Overview
Instrument Safety
Before Operating the Instrument
Documentation User Attention Words
Safety Alert Symbols1-4
Chemical Safety
Chemical Hazard Warning1-6
Chemical Waste Hazard Warning 1-6
About MSDSs 1-7
Ordering MSDSs 1-7
About Waste Disposal 1-8
GeneAmp PCR System 9700 Labels and Warnings 1-9
Instrument Safety Labels 1-9
Instrument Warnings Diagram 1-9
Danger of Burns
Electrical Shock Hazard1-10
Fire Hazard1-11
Grounding and Electrical Safety 1-11
Fuses
Laboratory Environmental Requirements
Introduction
Temperature, Humidity, and Environment1-12
Pollution
Emission/Immunity Statement
Routine Maintenance for Safe Operation1-12

#### 2 Product Overview

Overview	2-1
About This Chapter	2-1
In This Chapter	2-1
About the GeneAmp PCR System 9700	2-2
About the Instrument	2-2

Control Panel	2-2
Control Panel Diagram	2-2
Using the Keys	2-3
Selecting a Field	2-3
Entering Numeric Values	2-3
Transfer or Storage of Methods	2-4
Upgrading the Firmware	2-4
How the 9700 Differs from the GeneAmp PCR System 9600	2-5

#### 3 Instrument Setup

Overview
About This Chapter
In This Chapter
Unpacking the Instrument
Inspecting the Instrument
Additional Part Numbers
Setting Up your Laboratory
Choosing a Location
Operating Temperature
What to Avoid
Do Not Block Air Vents
Fuse Service Configuration    3-4
Input Voltage
Interchangeable Sample Block Modules
Converting from a Single to Double-Line Fuse
Introduction
Fuse Part Number   3-5
Power Entry Module Diagram
Converting Single to Double-line Fuse
Turning On the Instrument
Turning On the Instrument   3-7
Setting Custom Parameters
Overview
Displaying the Configuration Screen
Setting the Time
Setting the Date
Enabling or Disabling a Printer
Turning the Run-Time Beeper On or Off
Setting the Pause Time Out
Defining the Idle State Setpoint Temperature
Defining the Baud Rate for Your Printer Port
Configuring the Contrast and Screen Saver Options

Setting the IP Address	
Running Instrument Diagnostics	
Overview	
Viewing the Diagnostic Screen	
Testing the Electrical Components	
Testing the Display Screen	
Testing the Keypad	
Using the T <sub>m</sub> Calculator	
How to Use the T <sub>m</sub> Calculator	
Upgrading the System Firmware	
Introduction	
Upgrading Through the Serial Port	3-19
Downloading Firmware Using Windows 95 or Windows NT	
Downloading Firmware Using Windows 3.1	
Upgrading Firmware from a PCMCIA Flash Memory Card	
Copying Methods	
How to Copy Methods	
Connecting and Configuring a Printer	
Connecting a Printer	
Configuring the Printer	

#### 4 Running PCR Samples

Overview	. 4-1
About This Chapter	. 4-1
In This Chapter	. 4-1
Selecting a Method	. 4-2
Introduction	. 4-2
What Is a Method	. 4-2
Predefined Methods	. 4-2
Selecting a Method	. 4-3
Viewing Method Parameters	. 4-4
Searching for Methods	. 4-4
Sorting Methods	. 4-5
Running a Method	. 4-6
Starting a Run	. 4-6
About the Run Time Screen.	. 4-7
Viewing Method Information	. 4-8
Pausing a Run	. 4-8
Stopping a Run Before It Completes	. 4-9
Reviewing the History of a Run	4-10
How to Review the History of a Run	4-10
History Formats	4-10

]	History File Records	4-11
When a	a Run Completes	4-12
]	End of Run Screen	4-12
]	Method Not Stored Screen	4-12

#### 5 Creating and Editing Methods

Overview
About This Chapter
In This Chapter
Adding or Changing Users
Introduction
Adding a New User Name
Protecting Methods
Changing a User Name 5-5
Deleting a User Name
Creating Methods
About the Default Method
Basic Parameters. 5-6
Displaying the Create Methods Screen 5-7
Entering Temperature Control Parameters 5-8
Defining Pre-PCR Holds 5-9
Defining PCR Parameters 5-9
Defining Post-PCR Holds 5-10
Modifying Cycling
Introduction
Changing Temperature Control Parameters 5-11
Modifying Ramp Rates 5-12
Inserting Holds
Deleting Holds
Inserting Cycles
Inserting Programmed Pauses 5-14
Editing Programmed Pauses 5-16
Naming and Storing Methods 5-16
Printing a Method
Introduction
Printing a Method 5-17
Editing or Deleting Methods
Introduction
Editing a Method
Deleting a Method

#### 6 Converting Hold Times

Overview	. 6-1
About This Chapter	. 6-1
In This Chapter	. 6-1
About Setting Hold Times	. 6-2
Introduction	. 6-2
Hold Time Differences	. 6-2
Guidelines for Converting Hold Times	. 6-3
Introduction	. 6-3
What to Base the New Hold Times On	. 6-3
Setting Up Ramps	. 6-3
Table of Conversions	. 6-3
Up Ramp Example	. 6-4
Setting Down Ramps	. 6-4
Table of Conversions	. 6-5
Down Ramp Example	. 6-5

#### 7 Routine Maintenance

Overview	/-1
About This Chapter	7-1
In This Chapter	7-1
Changing the External Fuses	7-2
Introduction	7-2
Power Entry Module Diagram7	7-2
Changing Fuses for Single-Line Operation	7-3
Changing Fuses for Double-Line Operation7	7-4

#### 8 Troubleshooting

Overview	8-1
About This Chapter	8-1
In This Chapter	8-1
If There Is a Power Failure	8-2
Introduction	8-2
During a Power Failure	8-2
Display Screen Error Messages	8-3
Error Messages Table	8-3
Troubleshooting Information	8-6
Troubleshooting Table	8-6

#### A Instrument

#### **Specifications**

Overview
About This Appendix
In This Appendix
System Specifications
Dimensions
Power Configurations
Electrical Requirements
Control Panel Specifications
Display Screen
Keys
Sample Temperature Information
Temperature
Printer Specifications. A-3
Serial Interface Board Specifications
Cable Part Number

#### **B** Supplied Methods

Overview	B-1
About this Appendix	B-1
In This Appendix	B-1
About the Methods	B-2
Introduction	B-2
AmpliCycle Sequencing	B-2
AmpliTaq Gold	B-2
BigDye	B-3
General PCR	B-3
LMS2	B-4
Time Release PCR	B-4
Touchdown PCR	B-5
XL PCR	B-5

#### C Contacting Technical Support

Technical Support	. C-1
Contacting Technical Support	. C-1
To Contact Technical Support by E-Mail	. C-1
To Contact Technical Support by Telephone or Fax (North America)	. C-3
To Contact Technical Support by Telephone or Fax (Outside North America)	. C-4
To Reach Technical Support Through the Applied Biosystems Web Site	. C-6
To Obtain Technical Documents	. C-6

To Obtain Customer Training Information	C-7
---	-----

#### **D** Limited Warranty Statement

Index

# 1

# Introduction and Safety

#### **Overview**

About This Chapter This chapter provides information to help you safely operate the GeneAmp® System 9700.		
In This Chapter	The following topics are covered in this chapter:	
	Topics	See Page
	About This Manual	1-2
	Instrument Safety	1-3
	Chemical Safety	1-6
	GeneAmp PCR System 9700 Labels and Warnings	1-9
	Laboratory Environmental Requirements	1-12

#### **About This Manual**

	s manual describes how to use the GeneAmp <sup>®</sup> PCR System 9700. It includes the owing chapters and appendixes:
+	Chapter 1, "Introduction and Safety," contains safety information.
<b>+</b>	Chapter 2, "Product Overview," describes the instrument, its components, and requirements for installation.
+	Chapter 3, "Instrument Setup," provides information on how to set up the instrument and place it correctly in the laboratory, how to configure the instrument and how to connect and configure a printer.
+	Chapter 4, "Running PCR Samples," provides information on selecting a method starting and stopping a run, reviewing the history of a run, and what to do when the run is completed.
+	Chapter 5, "Creating and Editing Methods," describes how to create and edit PCF methods and how to work with stored methods.
+	Chapter 6, "Converting Hold Times," provides information about setting hold times for the GeneAmp <sup>®</sup> PCR System 9700 compared to the DNA Thermal Cycler or DNA Thermal Cycler 480.
<b>+</b>	Chapter 7, "Routine Maintenance," provides procedures for routine maintenance on the instrument.
+	Chapter 8, "Troubleshooting," lists error messages and suggestions for dealing with other problems you may encounter.
+	Appendix A, "Instrument Specifications," describes the dimensions, power, and electrical specifications of the GeneAmp® PCR System 9700 system, including the control panel, sample temperature information, and printer specifications.
+	Appendix B, "Supplied Methods," contains information about the methods that are supplied with the instrument.
<b>+</b>	Appendix C, "Contacting Technical Support," explains how to contact Applied Biosystems' Technical Support staff.
+	Appendix D, "Limited Warranty Statement," contains the instrument warranty statement.

#### **Instrument Safety**

<b>Before Operating the</b>	Ensure that everyone involved with the operation of the instrument has:
Instrument	<ul> <li>Received instruction in general safety practices for laboratories</li> </ul>
	<ul> <li>Received instruction in specific safety practices for the instrument</li> </ul>
	<ul> <li>Read and understood all related MSDSs</li> </ul>
	<b>ACAUTION</b> Avoid using this instrument in a manner not specified by Applied Biosystems. Although the instrument has been designed to protect the user, this protection can be impaired if the instrument is used improperly.
Documentation User Attention Words	Five user attention words appear in the text of all Applied Biosystems user documentation. Each word implies a particular level of observation or action as described below.
	Note Calls attention to useful information.
	<b>IMPORTANT</b> Indicates information that is necessary for proper instrument operation.
	<b>ACAUTION</b> Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury. It may also be used to alert against unsafe practices.
	<b>A WARNING</b> Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.
	<b>ADANGER</b> Indicates an imminently hazardous situation which, if not avoided, will result in death or serious injury. This signal word is to be limited to the most extreme situations.

# **Safety Alert Symbols** The following chart is an illustrated glossary of all electrical symbols that are used on Applied Biosystems instruments. Whenever such symbols appear on instruments, please observe appropriate safety procedures.

#### **Electrical Symbols**

	This symbol indicates the On position of the main power switch.
Ο	This symbol indicates the Off position of the main power switch.
Φ	This symbol indicates the On/Off position of a push-push main power switch.
Ŧ	This symbol indicates that a terminal may be connected to another instrument's signal ground reference. This is not a protected ground terminal.
	This symbol indicates that this is a protective grounding terminal that must be connected to earth ground before any other electrical connections are made to the instrument.
~	A terminal marked with this symbol either receives or delivers alternating current or voltage.
₹	A terminal marked with this symbol can receive or supply an alternating and a direct current or voltage.
A	This symbol indicates the presence of high voltage and warns the user to proceed with caution.
	This symbol alerts you to consult the manual for further information and to proceed with caution.

#### **Electrical Safety Testing**


Routine safety testing of analytical instruments (*e.g.*, high potential voltage testing) may be required by various safety agencies.



Testing should only be carried out by qualified personnel after seeking advice from the Applied Biosystems Service Department.

#### **Nonelectrical Symbols**

The following is an illustrated glossary of all nonelectrical safety alert symbols found on Applied Biosystems instruments.



#### **Chemical Safety**

Chemical Hazard Warning	• WARNING CHEMICAL HAZARD. Some of the chemicals used with Applied Biosystems instruments and protocols are potentially hazardous and can cause injury, illness, or death.
	<ul> <li>Read and understand the material safety data sheets (MSDSs) provided by the chemical manufacturer before you store, handle, or work with any chemicals or hazardous materials.</li> </ul>
	<ul> <li>Minimize contact with chemicals. Wear appropriate personal protective equipmer when handling chemicals (<i>e.g.</i>, safety glasses, gloves, or protective clothing). Fo additional safety guidelines, consult the MSDS.</li> </ul>
	<ul> <li>Minimize the inhalation of chemicals. Do not leave chemical containers open. Use only with adequate ventilation (<i>e.g.</i>, fume hood). For additional safety guidelines consult the MSDS.</li> </ul>
	<ul> <li>Check regularly for chemical leaks or spills. If a leak or spill occurs, follow the manufacturer's cleanup procedures as recommended on the MSDS.</li> </ul>
	<ul> <li>Comply with all local, state/provincial, or national laws and regulations related to chemical storage, handling, and disposal.</li> </ul>
Chemical Waste Hazard Warning	<b>WARNING</b> CHEMICAL WASTE HAZARD. Wastes produced by Applied Biosystems instruments are potentially hazardous and can cause injury, illness, or death.
	<ul> <li>Read and understand the material safety data sheets (MSDSs) provided by the manufacturers of the chemicals in the waste container before you store, handle, or dispose of chemical waste.</li> </ul>
	<ul> <li>Handle chemical wastes in a fume hood.</li> </ul>
	<ul> <li>Minimize contact with chemicals. Wear appropriate personal protective equipmer when handling chemicals (<i>e.g.</i>, safety glasses, gloves, or protective clothing). Fo additional safety guidelines, consult the MSDS.</li> </ul>
	<ul> <li>Minimize the inhalation of chemicals. Do not leave chemical containers open. Use only with adequate ventilation (<i>e.g.</i>, fume hood). For additional safety guidelines consult the MSDS.</li> </ul>
	<ul> <li>After emptying the waste container, seal it with the cap provided.</li> </ul>

About MSDSs Some of the chemicals used with this instrument may be listed as hazardous by their manufacturer. When hazards exist, warnings are prominently displayed on the labels of all chemicals.

Chemical manufacturers supply a current MSDS before or with shipments of hazardous chemicals to new customers and with the first shipment of a hazardous chemical after an MSDS update. MSDSs provide you with the safety information you need to store, handle, transport and dispose of the chemicals safely.

We strongly recommend that you replace the appropriate MSDS in your files each time you receive a new MSDS packaged with a hazardous chemical.

**A WARNING** CHEMICAL HAZARD. Be sure to familiarize yourself with the MSDSs before using reagents or solvents.

**Ordering MSDSs** You can order free additional copies of MSDSs for chemicals manufactured or distributed by Applied Biosystems using the contact information below.

To order documents by automated telephone service:

1	From the U.S. or Canada, dial <b>1.800.487.6809</b> , or from outside the U.S. and Canada, dial <b>1.858.712.0317</b> .	
2	Follow the voice instructions to order documents (for delivery by fax).	
	<b>Note</b> There is a limit of five documents per fax request.	

To order documents by telephone:

In the U.S.	Dial 1.800.345.5224, and press 1.	
	To order in English, dial 1.800.668.6913 and press 1, then 2, then 1	
In Canada	• To order in French, dial 1.800.668.6913 and press 2, then 2, then 1	
From any other country	See the specific region under "To Contact Technical Support by Telephone or Fax (Outside North America)".	

To view, download, or order documents through the Applied Biosystems web site:

Step	Action
1	Go to http://www.appliedbiosystems.com
2	Click <b>SERVICES &amp; SUPPORT</b> at the top of the page, click <b>Documents on Demand</b> , then click <b>MSDS</b> .
3	Click <b>MSDS Index</b> , search through the list for the chemical of interest to you, then click on the MSDS document number for that chemical to open a pdf of the MSDS.

For chemicals not manufactured or distributed by Applied Biosystems, call the chemical manufacturer.

About WasteAs the generator of potentially hazardous waste, it is your responsibility to perform the<br/>actions listed below.

- Characterize (by analysis if necessary) the waste generated by the particular applications, reagents, and substrates used in your laboratory.
- Ensure the health and safety of all personnel in your laboratory.
- Ensure that the instrument waste is stored, transferred, transported, and disposed of according to all local, state/provincial, or national regulations.

**Note** Radioactive or biohazardous materials may require special handling, and disposal limitations may apply.


#### **GeneAmp PCR System 9700 Labels and Warnings**

Labels


Instrument Safety Safety labels are located on the instrument. Each safety label has three parts:

- A signal word panel, which implies a particular level of observation or action (e.g., CAUTION or WARNING). If a safety label encompasses multiple hazards, the signal word corresponding to the greatest hazard is used.
- A message panel, which explains the hazard and any user action required.
- A safety alert symbol, which indicates a potential personal safety hazard.


The following diagram shows where the hazards and warnings labels are located on Instrument the GeneAmp® PCR System 9700. Warnings Diagram



**Danger of Burns AWARNING PHYSICAL INJURY HAZARD**. Hot Surface. Use care when working around this area to avoid being burned by hot components.



Electrical Shock Hazard AWARNING ELECTRICAL HAZARD. To reduce the chance of electrical shock do not remove covers that require tool access. No user serviceable parts are inside. Refer servicing to Applied Biosystems qualified service personnel.



The three-pronged power cord and receptacle contain the grounding connector.

**WARNING** ELECTRICAL HAZARD. Grounding circuit continuity is vital for safe operation of equipment. Never operate equipment with grounding conductor disconnected.




**A WARNING** ELECTRICAL SHOCK HAZARD. Severe electrical shock, which could cause physical injury or death, can result from working on an instrument when the high voltage power supply is operating. To avoid electrical shock, disconnect the power supply to the instrument, unplug the power cord, and wait at least 1 minute before working on the instrument.



Grounding and

**Electrical Safety** 

Fire Hazard AWARNING FIRE HAZARD. For continued protection against the risk of fire, replace fuses only with Listed and Certified fuses of the same type and rating as those currently in the instrument.



The system 9700 must be grounded for protection against electrical shock.

A CAUTION ELECTRICAL HAZARD. Do not use an adapter to a two-terminal outlet since this does not provide positive ground protection.

Fuses Improper fuses can damage the wiring system and cause a fire.

A WARNING ELECTRICAL FIRE HAZARD. Before turning on the instrument, verify that the fuses are properly installed.

#### Laboratory Environmental Requirements

Introduction	Take the precautions described in this section whenever you operate the system 9700. Read this section before you install the instrument.
	<b>CAUTION</b> The instrument should be used according to the instructions provided in this manual. If used otherwise, the protection provided by this instrument may be impaired.
Temperature,	<b>IMPORTANT</b> This instrument is designed for indoor use.
Humidity, and Environment	<b>IMPORTANT</b> Do not operate in a Cold Room or a refrigerated area. The system 9700 will operate safely when the ambient temperature is 5 °C to 40 °C (41 °F to 104 °F) and will meet performance specifications when the ambient temperature is 15 °C to 30 °C and the ambient relative humidity is 20 to 80%. These specifications have been calculated for altitudes between 0 and 2,000 meters.
	<b>ACAUTION</b> FIRE HAZARD. This instrument is not designed for operation in an explosive environment. Do not place the instrument close to potentially explosive materials or objects.
	<b>IMPORTANT</b> The instrument should be stored between $-20$ °C and $60$ °C ( $-4$ °F and $140$ °F) at altitudes between 0 and 12,000 meters.
	<b>Note</b> This instrument is able to withstand transient overvoltage according to Installation Category II as defined in IEC 1010-1.
Pollution	The installation category (overvoltage category) for this instrument is II, and it is classified as portable equipment. The instrument has a pollution degree rating of 2 and may be installed in an environment that has non-conductive pollutants only.
Emission/Immunity Statement	For our European customers, any product marked with the CE label meets the European EMC directive 89/336/EEC and the Low Voltage Directive 72/23/EEC. This product meets Class B emission limits.
Routine Maintenance for	If you use any cleaning or decontamination method, except those recommended in the manual, you may risk damaging the equipment.
Safe Operation	Maintain your instrument in good working order. In the event that the instrument has been subjected to adverse environmental conditions (such as fire, flood, earthquake, etc.), contact your local sales office for advice.


# 2

### **Product Overview**

#### Overview

About This Chapter This chapter describes the GeneAmp® PCR System 9700, its compon requirements for installation.		
In This Chapter	The following topics are contained in this chapter:	
	Торіс	See page
	About the GeneAmp PCR System 9700	2-2
	About the GeneAmp PCR System 9700	2-2

#### About the GeneAmp PCR System 9700



#### Using the Keys The following table describes the Control Panel keys in the previous figure.

Key	Use to
Soft keys (F1-F5)	Select the function specified above the key.
	The function of each key is defined on the display screen above the key, and is redefined as you view different screens.
Numeric keys	Enter numbers from left to right into a field you highlight.
Stop key	Stop a method while it is running.
Enter key	Enter information typed into a field and advance the highlight box to the next field on a screen.
Clear Entry key (CE)	Remove information from a field.
Circular Key Pad	Move the highlight box to different fields on the display screen in the direction of the arrow.

#### Selecting a Field There are two ways to select a field.

If you want to	Then
move the highlight box in one of four directions	Use the Circular Key Pad.
advance the box to the next field	Press the Enter key.

### Entering Numeric The following table lists how to enter numeric values for the Temperature and Hold Time parameters.

Parameter	Description
Temperature parameters	Enter values for temperature in decimal form. It is not necessary to type a decimal point.
	For example, to specify 89.0 °C, press 8 9 0, then press Enter.
Hold Time parameters	A hold time is the length of time the samples will be maintained at a specified temperature.
	Specify all hold times in minutes and seconds, then press Enter.
	For example, to specify one minute and five seconds, press 1 0 5, then press Enter.

#### Transfer or Storage You can transfer or store methods using a Methods Storage Card. The instrument can of Methods run a method from its own software or from the methods stored on a Methods Storage card.

Use	То	See
Methods Storage Card, Centennial 256kb SRAM (P/N 0940-1064)	transfer methods from the Methods Storage card to the instrument or from the instrument to the card for storage.	"Copying Methods" on page 3-25.

**Upgrading the** There are two ways to upgrade the firmware:

- Firmware 🗼 With a serial cable connection (PC communication cable P/N N805-1327 or Macintosh communication cable P/N N805-1328) from the instrument's RS485 ports
  - With a PCMCIA Flash Memory Card

Note Upgrade firmware can be ordered from the Applied Biosystems web site or by contacting PCR Technical Support.

For more information about upgrading the firmware, see the following.

If you are using	See
a serial cable connection (RS485 ports)	"Upgrading Through the Serial Port" on page 3-19.
a PCMCIA Flash Memory Card	"Upgrading Firmware from a PCMCIA Flash Memory Card" on page 3-23.
Microsoft Windows <sup>®</sup> 95 or Windows NT <sup>®</sup>	"Downloading Firmware Using Windows 95 or Windows NT" on page 3-21.
Windows <sup>®</sup> 3.1	"Downloading Firmware Using Windows 3.1" on page 3-22.

 How the 9700 Differs
 The GeneAmp PCR System 9700 differs from the GeneAmp® PCR System 9600 in the following ways:

 PCR System 9600
 Reaction volumes of up to 50 µL only can be run and in the "9600 Mode". as

- ✦ Reaction volumes of up to 50 µL only can be run and in the "9600 Mode", as opposed to 100 µL in the GeneAmp PCR System 9600.
- Methods are stored under a user name, allowing users to keep track of and protect their own methods.
- Hold, Cycle, and Auto programs no longer exist as menu items under the Create function. Instead of linking programs together to create a method, you create a single method on the GeneAmp PCR System 9700.
- If the allowed pause time elapses during a manual pause of a method, the method will continue running (instead of stopping).
- The last method run cannot be re-run or edited unless it is first saved.
- The history file no longer includes individual setpoint time and temperature information.
- An automated restart or incubate function is available for use after a power outage.

# 3

## Instrument Setup

#### **Overview**

About This Chapter This chapter provides information on how to set up the GeneAmp® PCR System 9700 and place it correctly in the laboratory, how to configure the instrument, and how to connect and configure a printer.

In This Chapter The following topics are covered in this chapter:

Торіс	See page
Unpacking the Instrument	3-2
Setting Up your Laboratory	3-3
Converting from a Single to Double-Line Fuse	3-5
Turning On the Instrument	3-7
Setting Custom Parameters	3-8
Running Instrument Diagnostics	3-14
Using the Tm Calculator	3-18
Upgrading the System Firmware	3-19
Copying Methods	3-25
Connecting and Configuring a Printer	3-27

#### Unpacking the Instrument

ry to ed in
-
ed in

**Note** See the *Interchangeable Sample Block Module Users Manual* for part numbers of disposable items.

#### Setting Up your Laboratory

Choosing a Location	To prevent vibration, place the GeneAmp PCR System 9700 on a solid, stable, level surface that allows free airflow overhead and around the sides and back. You should keep all ventilation slots in the instrument cover free of obstruction, for example from excess printer paper.
Operating Temperature	The instrument will meet performance specifications when the ambient temperature is 15 °C to 30 °C (59 ° to 86 °Fahrenheit) and the ambient relative humidity is 20% to 80%.
What to Avoid	When setting up the instrument, avoid:
	<ul> <li>Placing the instrument under overhanging shelves, especially when there is a wall behind the unit.</li> </ul>
	<ul> <li>Proximity to other instruments on the same bench or other heat-generating equipment.</li> </ul>
	<ul> <li>Locations subject to wide temperature fluctuation, such as direct sunlight, or air drafts.</li> </ul>
	♦ Damp areas.
	<ul> <li>Environments where there is an oil mist.</li> </ul>
Do Not Block Air Vents	Do not block the circulation of air to the vents located on the sides and bottom of the instrument
	air vents
	<b>Note</b> When operating multiple 9700 instruments, provide at least 8 in. of space in-between the instruments.

**Fuse Service** The instrument is shipped with single-line service configuration. If your installation requires double-line service, you must remove the jumper clip. See "Converting from a Single to Double-Line Fuse" on page 3-5.

**A WARNING** ELECTRICAL HAZARD. To protect yourself against shock hazards, use a properly wired three-terminal outlet. Do not use an adapter to a two-terminal outlet.

**A WARNING** ELECTRICAL FIRE HAZARD. Improper fuses can damage the wiring system and cause a fire. Before turning on the instrument, verify that the fuses are properly installed.

Input Voltage IMPORTANT You must be able to disconnect the main power supply to the instrument immediately if necessary.

The following table specifies the electrical operating range for the instrument in various parts of the world. Select the appropriate fuse configuration based on the voltage used.

Location	Voltage (VAC) <sup>a</sup>	Frequency	Amperage (A) Nominal
Japan	100 ±10%	50/60 Hz ±1%	3.16
USA/Canada	120 ±10%	50/60 Hz ±1%	4.20
EC	230 ±10%	50/60 Hz ±1%	3.14

a. Acceptable AC line voltage tolerances: 100, 120, 220, 230  $\pm$ 10%; 240 VAC +6%/-10%, 50/60 Hz  $\pm$  1%.

**Note** The Volt-Amp number for this instrument is 725 Volt Amps.

Interchangeable The instrument features an interchangeable sample block module allowing portability sample Block and versatility in sample configuration.


#### Modules

Sample Block Module	Part Number
96-Well Gold Sample Block Module	4314443
96-Well Aluminum Sample Block Module	4314445
60-Well 0.5 mL Sample Block Module	4309131
Dual 384-Well Sample Block Module	N805-0400
Auto-Lid Sample Block Module	4312904

**Note** The GeneAmp PCR System 9700 will not operate without an interchangeable sample block module installed.

For an example, refer to the *GeneAmp PCR System 9700 96-Well Sample Block Module User Guide*, for instructions on how to install a block and configure the samples.

#### **Converting from a Single to Double-Line Fuse**



# Converting Single to The following procedure describes how to convert from a single to a double-line fuse. Double-line Fuse Note Remove the jumper clip from the fuse compartment to convert the instrument from single-line to double-line service.


To convert from a single to a double-line fuse:

Step	Action
1	Disconnect the power cord from the power entry module.
	<b>WARNING</b> ELECTRICAL SHOCK HAZARD. Severe electrical shock, which could cause physical injury or death, can result from working on an instrument when the high voltage power supply is operating. To avoid electrical shock, disconnect the power supply to the instrument, unplug the power cord, and wait at least 1 minute before working on the instrument.
2	Insert a fine flat-tip screwdriver into the slot in the upper portion of the power entry module, and open the door.
3	Pull the fuse compartment out.
	<b>Note</b> Observe the orientation of the metal bars that run from the fuse compartment to the two tabs.
4	Remove the jumper clip from the end portion of the fuse compartment.
5	Insert one fuse into the left side of the compartment.
6	Place the fuse compartment back into the Power Entry Module and close the fuse compartment door.
	Press the fuse compartment door until it locks in place. Connect the instrument power cord.

#### **Turning On the Instrument**

Turning On the Instrument

**Turning On the** To turn on the instrument:



#### **Setting Custom Parameters**

**Overview** The instrument is shipped with default configuration values. The following procedure describes how to set customized values as well as how to enable or disable optional features.

Displaying the Configuration	To display the Configuration Screen:		
Screen	Step	Action	
	1	From the Main menu, press Util.	
		The Utilities screen appears.	
		Utilities	
		Diag – Instrument diagnostics TmCalc – Calculates melting temp Config – Instrument configuration	
		Diag TmCalc Config More Exit	
		F1 F2 F3 F4 F5	
		From the Utilities screen, you can take the following action:	
		If you want to	See page
		Delete a method	5-20
		Run instrument diagnostics	3-14
		Review the history of a run	4-10
		Access the T <sub>m</sub> calculator	3-18
		Copy methods to or from a methods transportability card	3-25
	2	Press <b>Config</b> . The first Configuration screen appears.	
		Instrument ConfigurationTime:11:30 AMDate:01/25/00 M/D/YRun Time Printer: OffRun Time Beep:OffAccept24HrPMMoreF1F2F3F4F5	

Setting the Time In the first Configuration screen, you can set the current time and date for file memory maintenance, and run-time displays. You can also enable or disable the run-time printer and the run-time beep.

To set the time:

Step	Action	
1	Use the Circular Key Pad to select the <b>Time</b> field.	
2	Press the <b>24 Hr</b> or <b>PM (AM)</b> soft keys until the format you want for the current time displays in the Time field.	
3	Use the numeric keys to type in the hours followed by minutes.	
4	Press Accept when your entries are complete.	
	Note CE clears an entry.	

Setting the Date You must set the instrument to the correct date.

To set the date:

Step	Action	
1	Use the circular key to select the Date field.	
2	The three fields to set in the Date field are the:	
	♦ Days field	
	♦ Month field	
	♦ Year field	
3	Press the <b>D/M/Y</b> or <b>Y/M/D</b> soft keys until the format you want for the current date displays in the Date field.	
4	Use the numeric keys and type in a number for each of the three fields.	
	The order of these three fields depends on the format you chose in step 3.	

**Disabling a Printer** 

Enabling or Enabling the printer allows you to print method parameters or records of run time events directly from the display screen. The default value for the optional printer is Off.

Note For connection and configuration see "Connecting and Configuring a Printer" on page 3-27.

To enable or disable the optional printer:

Step	Action		
1	Select the Run Time Printer field. This changes the functions of the soft keys.		
2	You can take the following action:		
	If you want to	Then	
	enable the printer	Press ON.	
	disable the printer	Press Off.	
3	Press Enter to accept yo	our entry.	
4	Press <b>Cancel</b> to cancel all entries and return to the previous screen.		

Turning the When turned on, the run-time beeper beeps during a pause and once at the Run-Time Beeper completion of a run. The default value is Off. On or Off

To turn the Run Time-Beeper on or off:

Step	Action		
1	Select the Run Time Beep field. This changes the value of the soft keys.		
2	You can take the following action:		
	If you want to	Then	
	turn the beeper on	Press ON.	
	turn the beeper off	Press Off.	
3	Press Enter to accept you	ur entry.	

### Setting the Pause The Pause Time Out field sets the time in minutes: seconds format for the length of a Time Out pause when you use the Pause soft key to pause a run from the Run Time screen. For more information on Pausing a Run, see "Pausing a Run" on page 4-8.

To set the pause time out:

Step	Action		
1	Select More to display the second configuration screen.		
	Instrument Configuration		
	Pause Time Out: 10:00 (00:01-99:59)		
	Idle State Setpoint: 25.0°C (4.0-99.9) Baud Rate: 9600		
	Accept More Cancel		
	F1 F2 F3 F4 F5		
2	Set the <b>Pause Time Out</b> field by using the numeric keys to enter a time in minutes and seconds.		
	From this screen, you can also define the:		
	<ul> <li>Idle state setpoint temperature.</li> </ul>		
	<ul> <li>Baud rate for your printer port.</li> </ul>		
3	Press Accept when all information on this screen is correct.		
	Note CE clears an entry.		

### Defining the Idle State Setpoint Temperature

The Idle State Set Point temperature is the temperature at which the instrument will remain when powered up, but idle.

e IMPORTANT After a run is completed or terminated, there is approximately a 30 second delay before the instrument attains the specified idle state temperature. This allows you to stop one method and start another before the instrument temperature changes.

To define the Idle State Set Point temperature:

Step	Action	
1	Select the Idle State Set Point field.	
2	Use the numeric keys and type in a temperature between 4.0° C and 99.9 °C.	
	Note The default is 25 °C.	
3	Press Enter to accept your entry.	

### Defining the Baud Rate for Your Printer Port

The following section contains a table that lists serial board specifications and a procedure that describes how to define the printer port value in the Baud Rate field.

### **Serial Board Specifications**

You can connect the GeneAmp PCR System 9700 to any printer with a serial interface board and the following specifications.

Baud Rate	9600
Parity	NONE
Data Bits	8
Stop Bits	1

### How to Set the Baud Rates

To set baud rates:

Step	Action			
1	Select the Baud Rate field.			
2	You can take the following action:			
	If you want to Then			
	increase the baud rate Press Up.		-	
	decrease the baud rate Press Down.			
	200, 600, and 300.			
3	Press Enter to accept your e	entry.		

Configuring the Contrast and Screen Saver Options You can use the third Configuration screen to set: Screen contrast from 1 to 20 (light to dark)

Screen saver mode

### **Displaying the Screen**

From the second Configuration screen, choose More to display the third Configuration screen.

### **Screen Saver Options**

The following table lists the screen saver options:

Use This Option	To Activate the Screen Saver
Smart	After 15 minutes unless a method is running
Always	After 15 minutes
Never	Never

Setting the You can use the fourth Configuration screen to set the IP address for the instrument. IP Address The instrument can be linked to a network via the RS485 ports located on the right-rear side of the instrument (see the figure on page 3-19). Once connected and configured, the instrument will be visible from the network under the IP address and instrument name you select in the following procedure.

**Note** See the System 9700 Networking Software User's Manual (P/N 4309575) for more information about networking GeneAmp PCR System 9700 instruments.

To set the IP Address for an instrument:

Step	Action			
1	Select More to display the second configuration screen.			
	Instrument Configuration Set IP: RS485 Name: Cycler 01			
	IP:	1		
	Accept	- + Cancel		
	Fl	F2 F3 F4 F5		
2	Set the Set I	P field using the + and – keys.		
	The following settings are available:			
	Setting Definition			
	RS485	Activates the RS 485 ports.		
		Must be used with the networking software		
	Off	Turns off the networking capability and turns on the RS232 ports.		
3	Use the Circ	ular Key Pad to select the <b>Name</b> field.		
4	Select Name	ł.		

To set the IP Address for an instrument: (continue
--

Step	Action			
5	Using the Circular Key Pad and <b>Enter</b> key, enter up to a 16 character name for identifying the instrument on the network.			
	abcdefghi jklmnopqr IP Name: Cycler 01 stuvwxyz			
	Image: A constraint of the const			
	F1 F2 F3 F4 F5			
6	Press Accept when finished entering a name.			
	The software returns to the previous screen and displays the new name in the Name field.			
7	Use the Circular Key Pad to select the last IP field.			
8	Set the IP field by using the numeric keys to enter an IP address for the instrument.			
	Instrument Configuration			
	Set IP: RS485 Name: Cycler 01 IP: 0. 0. 01			
	Accept More Cancel			
	F1 F2 F3 F4 F5			
9	Press Accept when all information on this screen is correct.			

### **Running Instrument Diagnostics**

**Overview** The instrument provides a number of internal diagnostic utilities, tests you can run to verify that the instrument hardware and software components meet performance specifications.

Viewing the The following Diagnostic Screen

Viewing the The following procedure describes how to view the Diagnostic screen.

To view the Diagnostic screen:

Step	Action					
1	From the Main menu, press	Util.				
	The Utilities screen appears.					
	Utilities Diag - Instrument diagnostics TmCalc - Calculates melting temp Config - Instrument configuration					
	Diag TmCalc (	Config More Exit				
	F1 F2	F3 F4 F5				
2	From the Utilities screen, pr	ress Diag.				
	The Diagnostics screen app	bears.				
	Diagnostics Hard - Hardware Diagnostics System - System Performance Tests TmpVer - Temperature Verification Upgrad - Firmware Upgrade Hard System TmpVer Upgrad Exit					
	F1 F2	F3 F4 F5				
3	Press the soft key that acce	ccesses the diagnostic utility you want to run:				
	Press this soft key	То				
	Hard	Access hardware diagnostic utilities that allow you to test hardware and electrical components such as the Display screen and Keypad.				
	System	Access the system performance tests.				
		You can run these two system tests to test the rate of system heating and cooling, and the performance of PCR cycling.				
	TmpVer	Access the temperature verification tests that allow you to verify sample block calibration and uniformity.				
		These two tests are discussed in the users manual for the interchangeable sample block.				
	Upgrade	Access a utility that allows you to update the 9700 instrument firmware through the RS485 serial port or using a PCMCIA Flash Memory card.				

Testing the<br/>ElectricalUse the hardware diagnostic utilities to test the electrical components of the 9700<br/>instrument.

Components

To test the electrical components:

Step	Action					
1	From the Diagnostics screer	n, press <b>Hard</b> .				
	The Hardware Diagnostics appears.					
	Disp – LCD Display Diagnostic Keypad – Keypad Diagnostic					
	Disp [	Keypad Exit				
	F1 F2	F3 F4 F5				
2	want to perform one of the fo					
	Press this soft key	То				
	Display	Visually determine if the Liquid Crystal Display (LCD) screen is functioning properly.				
		(LCD) screen is functioning properly.				

Testing the Display<br/>ScreenThe Display diagnostic test allows you to visually determine if the display screen is<br/>properly functioning by turning on and off all the LCD pixels.

To test the display screen:

Step	Action						
1	From the Hardware Diagnostics screen, press Display.						
	The Display Diagnostic screen appears.						
	Display Dia	agnostics					
	<ol> <li>Read all instructions first.</li> <li>Press Run to turn ON all pixels.</li> <li>Press STOP to turn OFF all pixels.</li> <li>Press STOP to exit.</li> </ol>						
	Run	Run Exit					
	F1 F2 F3	B F4 F5					
2	You can take the following action:	:					
	Press To						
	Run Turr	n ON all pixels.					
	Che	eck display.					
	Stop Turr	n OFF all pixels.					
	Che	eck display.					
	Stop, again Exit	t.					

**Testing the Keypad** Use the Keypad diagnostic test to verify that all 22 keys on the control panel are functioning properly.

To test the keypad:

Step	Action					
1	From the Hard	ware Diagno	stics screen	, press <b>Keyp</b>	oad.	
	The Keypad D	iagnostic sci	reen appear	S.		
	After pr key or p					
	Run F1	F2	F3	F4	EXIC F5	
2	Press <b>Run</b> . The Control Pa					
	STOP	F	F2     F3     F       1     2       4     5       7     8	F4 F5 3 6 9 2E		
	F1	F2	F3	F4	F5	
3	Press the indic	cated flashin	g key.			

### Using the T<sub>m</sub> Calculator

> Action Step 1 From the Main menu, press Util. The Utilities screen opens as shown below. Utilities - Instrument diagnostics Diag TmCalc - Calculates melting temp Config - Instrument configuration TmCalc Config More Exit Diag F1 F2 F3 F4 F5 2 Press TmCalc. The Tm Calculator appears. [Salt]: 50 mM [Primer] 0.20 uM P1: 5' P2: 5′ Tm of P1= Tm of P2= Press ENTER to calculate Tm's Enter the salt concentration. 3 Note The default is 50. Enter values 5 to 1000. 4 Enter the primer concentration. Note The default it 0.20. Enter values 0.01 to 10.00. 5 Enter primer sequence in P1. Enter primer sequence in P2 and press Enter to calculate the T<sub>m</sub>s. 6 The melting points are displayed. Use this information to program a run. For more information, see Chapter 5, "Creating and Editing Methods." 7 Press Return to display the Utilities screen.

To use the  $T_m$  Calculator:

## Upgrading the System Firmware

Introduction	You can update the GeneAmp PCR System 9700 firmware from a serial connection to the RS485 port or from a PCMCIA Flash Memory Card.				
	Contact Memory		esentative to obtain a PCMCIA Flash		
Upgrading Through the Serial Port	The follo serial po	•	pgrade the firmware through the RS485		
	To upgra	ade the firmware through the RS485	5 serial port:		
	Step	Action			
	1	Connect the appropriate communicati on the 9700 instrument Base Module	on cable from a computer to the RS485 port (see the figure below).		
		You can use the following cables:			
		Communication cable	Part number		
		PC communication cable	N805-1327		
		Macintosh communication cable	N805-1328		
	2	For information on connecting the print	RS485 Ports		
	2	The Firmware Upgrade screen appea			
		The Fillinware Opgrade Screen appea	15.		
		Firmware Upg PCcard - Upgrade from Serial - Upgrade from PCcard Serial	m PC card		
		F1 F2 F3	F4 F5		

To upgrade the firmware through the RS485 serial port: (continued)

From the Firmware Upgrade			Action					
3 From the Firmware Upgrade screen, press Serial. Serial Port Firmware Upgrade								
						Attach cable from RS485 port on 9700 to serial port on computer. Press Update, then begin download of file from the computer.		
F1 F2	F3 F	74	F5					
<b>IMPORTANT</b> Do not remo	ove the cable or tu	urn off the i	nstrume	nt until the upgrade				
From the Serial Port Firmware Upgrade screen, press Update.								
DO NOT TURN OFF THE INSTRUMENT								
Upgrading over the serial port will take approximately 10 minutes. Instrument will re-boot when complete.								
· · · · · · · · · · · · · · · · · · ·								
Take the following action:	<u> </u>	1	1.2					
If you are using See								
If you are using     See       Microsoft <sup>™</sup> Windows <sup>®</sup> "Downloading Firmware Using Windows 95 or       95, or Microsoft     Windows NT" on page 3-21.       Windows NT <sup>®</sup> Windows NT <sup>®</sup>								
Microsoft Windows® 3.1	"Downloading F page 3-22.	irmware Us	sing Win	dows 3.1" on				
	Attach cable from serial port on conserial port on conserial port on conserial port on conservation downloc computer.         F1       F2         IMPORTANT Do not remains completed.         From the Serial Port Firmward         DO NOT TURN OR REMOVE         Upgrading ov         will take approximation of the following action:         F1       F2         Take the following action:         If you are using         Microsoft <sup>™</sup> Windows <sup>®</sup> 95, or Microsoft         Windows NT <sup>®</sup> Microsoft Windows <sup>®</sup>	Attach cable from RS485 port serial port on computer. Pretthen begin download of file computer.         Update         F1       F2       F3       F         IMPORTANT       Do not remove the cable or to is completed.         From the Serial Port Firmware Upgrade scree         DO       NOT       TURN       OFF       THE       INST OR         DO       NOT       TURN       OFF       THE       INST OR         Upgrading       over       the serial will       take approximately       10         Instrument       will       re-boot       where         F1       F2       F3       F         Take the following action:       See       Microsoft <sup>™</sup> Windows <sup>®</sup> "Downloading F         95, or Microsoft       Windows NT <sup>®</sup> "Downloading F       Windows NT <sup>®</sup> Microsoft Windows <sup>®</sup> "Downloading F       "Downloading F	Attach cable from RS485 port on 970         serial port on computer. Press Upda         then begin download of file from th         computer.         IMPORTANT         Do not remove the cable or turn off the i         is completed.         From the Serial Port Firmware Upgrade screen, press Upgrading over the serial port         OR REMOVE SERIAL CABLE!         Upgrading over the serial port         will take approximately 10 minute         Instrument will re-boot when completed         F1       F2         F3       F4         Take the following action:       See         Microsoft <sup>™</sup> Windows <sup>®</sup> "Downloading Firmware Us         95, or Microsoft       "Downloading Firmware Us         Windows NT <sup>®</sup> "Downloading Firmware Us         Microsoft Windows <sup>®</sup> "Downloading Firmware Us	Attach cable from RS485 port on 9700 to serial port on computer. Press Update, then begin download of file from the computer.         Update       Exit         F1       F2       F3       F4       F5         IMPORTANT Do not remove the cable or turn off the instrument is completed.       From the Serial Port Firmware Upgrade screen, press Update.         From the Serial Port Firmware Upgrade screen, press Update.         DO NOT TURN OFF THE INSTRUMENT OR REMOVE SERIAL CABLE!         Upgrading over the serial port will take approximately 10 minutes.         Instrument will re-boot when complete.         F1       F2       F3       F4       F5         Take the following action:       See       Microsoft™ Windows®       "Downloading Firmware Using WindWindows NT" on page 3-21.         Microsoft Windows®       "Downloading Firmware Using WindWindows NT" on page 3-21.				

1

Downloading Firmware Using Windows 95 or Windows NT

**Downloading** To download the system firmware using Windows 95 or Windows NT:

1       Click Start > Programs > Accessories > HyperTerminal.         2       Click the HyperTerminal icon. The Connection Description dialog box appears.         3       In the Connection Description dialog box, enter the following information a OK. The Connect Using dialog box appears.         1       In this field         Name       Enter a name, for example TEC.         Icon       Choose one of the icons.         4       In the Connect Using dialog box, connect using either Direct to Com 1 or D Com 2, depending on the port to which the cable is attached.         Note       For Windows NT applications this dialog box is titled Connect To, a fields are COM1 and COM2.         5       Depending on whether you chose Com 1 or Com 2 in the previous step, in Properties dialog box, enter the following.         In this field       Enter         Bits Per Second       9600         Data Bits       8         Parity       None         Stop Bits       1         Flow Control       None         When you have competed entering the information, click OK.         6       From the Transfer menu, choose Send file. The Send File dialog box, enter the following and click Send.         7       In this field       Take this action Filename         Locate the upgrade file with the name under which y saved the file.       Protocol <th>tep</th> <th>Action</th> <th></th> <th></th> <th></th>	tep	Action					
The Connection Description dialog box appears.         In the Connect Using dialog box, enter the following information a OK. The Connect Using dialog box appears.         In this field       Take this action         Name       Enter a name, for example TEC.         Icon       Choose one of the icons.         In the Connect Using dialog box, connect using either Direct to Com 1 or D Com 2, depending on the port to which the cable is attached.         Note       For Windows NT applications this dialog box is titled Connect To, a fields are COM1and COM2.         Depending on whether you chose Com 1 or Com 2 in the previous step, in Properties dialog box, enter the following.         In this field       Enter         Bits Per Second       9600         Data Bits       8         Parity       None         Stop Bits       1         Flow Control       None         When you have competed entering the information, click OK.         From the Transfer menu, choose Send file.         The Send File dialog box, enter the following and click Send.         In this field       Take this action         Filename       Locate the upgrade file with the name under which y saved the file.	1	Click Start > Programs > Accessories > HyperTerminal.					
<ul> <li>In the Connection Description dialog box, enter the following information a OK. The Connect Using dialog box appears.</li> <li>In this field Take this action <ul> <li>Name</li> <li>Enter a name, for example TEC.</li> <li>Icon</li> <li>Choose one of the icons.</li> </ul> </li> <li>In the Connect Using dialog box, connect using either Direct to Com 1 or D Com 2, depending on the port to which the cable is attached.</li> <li>Note For Windows NT applications this dialog box is titled Connect To, a fields are COM1 and COM2.</li> <li>Depending on whether you chose Com 1 or Com 2 in the previous step, in Properties dialog box, enter the following.</li> <li>In this field Enter</li> <li>Bits Per Second</li> <li>9600</li> <li>Data Bits</li> <li>Parity</li> <li>None</li> <li>Stop Bits</li> <li>In Flow Control</li> <li>None</li> <li>When you have competed entering the information, click OK.</li> <li>From the Transfer menu, choose Send file.</li> <li>The Send File dialog box, appears.</li> <li>In this field Take this action</li> <li>Filename</li> <li>Locate the upgrade file with the name under which y saved the file.</li> </ul>	2	Click the HyperTern	ninal icon.				
OK. The Connect Using dialog box appears.         In this field       Take this action         Name       Enter a name, for example TEC.         Icon       Choose one of the icons.         In the Connect Using dialog box, connect using either Direct to Com 1 or D         Com 2, depending on the port to which the cable is attached.         Note       For Windows NT applications this dialog box is titled Connect To, a fields are COM1 and COM2.         Depending on whether you chose Com 1 or Com 2 in the previous step, in Properties dialog box, enter the following.         In this field       Enter         Bits Per Second       9600         Data Bits       8         Parity       None         Stop Bits       1         Flow Control       None         When you have competed entering the information, click OK.         From the Transfer menu, choose Send file.         The Send File dialog box, enter the following and click Send.         In this field       Take this action         Filename       Locate the upgrade file with the name under which y saved the file.		The Connection Description dialog box appears.					
Name       Enter a name, for example TEC.         Icon       Choose one of the icons.         4       In the Connect Using dialog box, connect using either Direct to Com 1 or D Com 2, depending on the port to which the cable is attached.         Note       For Windows NT applications this dialog box is titled Connect To, a fields are COM1 and COM2.         5       Depending on whether you chose Com 1 or Com 2 in the previous step, in Properties dialog box, enter the following.         In this field       Enter         Bits Per Second       9600         Data Bits       8         Parity       None         Stop Bits       1         Flow Control       None         When you have competed entering the information, click OK.         6       From the Transfer menu, choose Send file.         The Send File dialog box, enter the following and click Send.         In this field       Take this action         Filename       Locate the upgrade file with the name under which y saved the file.	3		-	-	-		
Icon       Choose one of the icons.         4       In the Connect Using dialog box, connect using either Direct to Com 1 or D Com 2, depending on the port to which the cable is attached.         Note       For Windows NT applications this dialog box is titled Connect To, a fields are COM1and COM2.         5       Depending on whether you chose Com 1 or Com 2 in the previous step, in Properties dialog box, enter the following.         In this field       Enter         Bits Per Second       9600         Data Bits       8         Parity       None         Stop Bits       1         Flow Control       None         When you have competed entering the information, click OK.         6       From the Transfer menu, choose Send file.         The Send File dialog box, enter the following and click Send.         In this field       Take this action         Filename       Locate the upgrade file with the name under which y saved the file.		In this field	Take t	his action			
<ul> <li>In the Connect Using dialog box, connect using either Direct to Com 1 or D Com 2, depending on the port to which the cable is attached.</li> <li>Note For Windows NT applications this dialog box is titled Connect To, a fields are COM1and COM2.</li> <li>Depending on whether you chose Com 1 or Com 2 in the previous step, in Properties dialog box, enter the following.</li> <li>In this field Enter</li> <li>Bits Per Second 9600</li> <li>Data Bits 8</li> <li>Parity None</li> <li>Stop Bits 1</li> <li>Flow Control None</li> <li>When you have competed entering the information, click OK.</li> <li>From the Transfer menu, choose Send file.</li> <li>The Send File dialog box, enter the following and click Send.</li> <li>In this field</li> <li>Take this action</li> <li>Filename Locate the upgrade file with the name under which y saved the file.</li> </ul>		Name	Enter	a name, fo	r example TEC.		
Com 2, depending on the port to which the cable is attached.         Note       For Windows NT applications this dialog box is titled Connect To, a fields are COM1and COM2.         5       Depending on whether you chose Com 1 or Com 2 in the previous step, in Properties dialog box, enter the following.         In this field       Enter         Bits Per Second       9600         Data Bits       8         Parity       None         Stop Bits       1         Flow Control       None         When you have competed entering the information, click OK.         6       From the Transfer menu, choose Send file.         The Send File dialog box, enter the following and click Send.         7       In this field         Take this action         Filename       Locate the upgrade file with the name under which y saved the file.		lcon	Choos	se one of th	ne icons.		
5       Depending on whether you chose Com 1 or Com 2 in the previous step, in Properties dialog box, enter the following.         In this field       Enter         Bits Per Second       9600         Data Bits       8         Parity       None         Stop Bits       1         Flow Control       None         When you have competed entering the information, click OK.         6       From the Transfer menu, choose Send file.         The Send File dialog box, enter the following and click Send.         In this field       Take this action         Filename       Locate the upgrade file with the name under which y saved the file.	4	Com 2, depending of Note For Window	on the port to vs NT applicat	which the	cable is attached.		
Properties dialog box, enter the following.         In this field       Enter         Bits Per Second       9600         Data Bits       8         Parity       None         Stop Bits       1         Flow Control       None         When you have competed entering the information, click OK.         6       From the Transfer menu, choose Send file.         The Send File dialog box appears.         7       In the Send File dialog box, enter the following and click Send.         In this field       Take this action         Filename       Locate the upgrade file with the name under which y saved the file.	5			e Com 1 oi	r Com 2 in the previous step, in the		
Bits Per Second       9600         Data Bits       8         Parity       None         Stop Bits       1         Flow Control       None         When you have competed entering the information, click OK.         6       From the Transfer menu, choose Send file.         The Send File dialog box appears.         7       In the Send File dialog box, enter the following and click Send.         In this field       Take this action         Filename       Locate the upgrade file with the name under which y saved the file.	•		-				
Bits Per Second       9600         Data Bits       8         Parity       None         Stop Bits       1         Flow Control       None         When you have competed entering the information, click OK.         6       From the Transfer menu, choose Send file.         The Send File dialog box appears.         7       In the Send File dialog box, enter the following and click Send.         In this field       Take this action         Filename       Locate the upgrade file with the name under which y saved the file.			<b>–</b>		1		
Data Bits       8         Parity       None         Stop Bits       1         Flow Control       None         When you have competed entering the information, click OK.         6       From the Transfer menu, choose Send file.         The Send File dialog box appears.         7       In the Send File dialog box, enter the following and click Send.         In this field       Take this action         Filename       Locate the upgrade file with the name under which y saved the file.					-		
Parity       None         Stop Bits       1         Flow Control       None         When you have competed entering the information, click OK.         6       From the Transfer menu, choose Send file.         The Send File dialog box appears.         7       In the Send File dialog box, enter the following and click Send.         In this field       Take this action         Filename       Locate the upgrade file with the name under which y saved the file.					-		
Stop Bits       1         Flow Control       None         When you have competed entering the information, click OK.         6       From the Transfer menu, choose Send file.         The Send File dialog box appears.         7       In the Send File dialog box, enter the following and click Send.         In this field         Take this action         Filename       Locate the upgrade file with the name under which y saved the file.					-		
Flow Control       None         When you have competed entering the information, click OK.         6       From the Transfer menu, choose Send file.         The Send File dialog box appears.         7       In the Send File dialog box, enter the following and click Send.         In this field       Take this action         Filename       Locate the upgrade file with the name under which y saved the file.			None		4		
When you have competed entering the information, click OK.         6       From the Transfer menu, choose Send file.         The Send File dialog box appears.         7       In the Send File dialog box, enter the following and click Send.         In this field         Take this action         Filename       Locate the upgrade file with the name under which y saved the file.		Stop Bits	1				
<ul> <li>From the Transfer menu, choose Send file.</li> <li>The Send File dialog box appears.</li> <li>In the Send File dialog box, enter the following and click Send.</li> <li>In this field</li> <li>Take this action</li> <li>Filename</li> <li>Locate the upgrade file with the name under which y saved the file.</li> </ul>		Flow Control	None				
The Send File dialog box appears.         7       In the Send File dialog box, enter the following and click Send.         In this field       Take this action         Filename       Locate the upgrade file with the name under which y saved the file.		When you have competed entering the information, click <b>OK</b> .					
7       In the Send File dialog box, enter the following and click Send.         In this field       Take this action         Filename       Locate the upgrade file with the name under which y saved the file.	6	From the Transfer menu, choose Send file.					
In the Send File dialog box, enter the following and click Send.         In this field       Take this action         Filename       Locate the upgrade file with the name under which y saved the file.		The Send File dialo	g box appear	S.			
Filename     Locate the upgrade file with the name under which y saved the file.	7						
saved the file.		In this field	Take this a	action			
Protocol Enter Kermit.		Filename	Locate the upgrade file with the name under which you saved the file.				
		A dialog box appears that displays the current packet number and the number					
automatically resets.							

	-	
Downloading Firmware Using	To dowr	load the system firmware using Windows 3.1:
Windows 3.1	Step	Action


Step	Action							
1	From the Program Manager, click the Terminal program icon.							
2	From the Settings r	nenu, choose <b>Communication</b> . A dialog box appears.						
3	Enter the following	information in the dialog box and click <b>OK</b> .						
	In this field Enter							
	Baud Rate	9600						
	Data Bits	8						
	Stop Bits	1						
	Parity	None						
	Flow Control	None						
	Connector Com1 or Com2 (depending on the port to which the cable is attached).							
4	From the Binary Tra	ansfers menu, choose Settings. A dialog box appears.						
5	Select Kermit and o	click <b>OK</b> .						
6	From the Send Binary File menu, choose Transfers.							
	A directory dialog b	oox appears.						
7	From the directory dialog box, find and select the <b>9700 upgrade file</b> , and click <b>OK</b> .							
	A message at the bottom of the window indicates that the program is sending the file and a progress bar appears.							
	When the download is complete the GeneAmp PCR System 9700 automatically resets.							
8	When the download parameters.	d is complete, from the File menu, choose Save to save these						

Upgrading Firmware from a PCMCIA Flash Memory Card

The following procedure describes how to upgrade the firmware using a PCMCIA Flash Memory Card.

**Note** Contact your local Applied Biosystems representative to obtain a PCMCIA Flash Memory Card.

To upgrade using a PCMCIA Flash Memory Card:




To upgrade using a PCMCIA Flash Memory Card: (continued)

Step	Action						
6	If the PCMCIA Flash Memory Card contains a valid firmware upgrade, the <b>Upgrade</b> screen appears.						
	DO NOT TURN OFF THE INSTRUMENT OR REMOVE THE PC CARD!						
	Upgrading the instrument firmware will take approximately 30 seconds. Instrument will re-boot when complete.						
	F1 F2 F3 F4 F5						
	If successful, the <b>Firmware Upgrade Successful</b> screen displays the firmware version to which you have upgraded. Pressing the <b>Run</b> key in the lower right corner brings up the <b>Main</b> menu with the new						
	revision number of the software listed on the screen.						
7	Remove the PCMCIA Flash Memory Card after you have successfully upgraded the firmware.						

### **Copying Methods**

### **How to Copy Methods** You can copy methods from a Methods Transportability card to the instrument or from the instrument to a Methods Transportability card (P/N 940-1064). Runs can be made directly from methods on a Methods Transportability card.

To copy a method:



### To copy a method: (continued)

Step	Action							
4	In the Copy Method dialog box, you can take the following action:							
	Press	То						
	ToInst	Transfer a instrumen		om a Methods	s Transportat	ility card to the		
	ToCard		Transfer a method from the instrument to a Methods Transportability card.					
	A screen appears with a list of the methods on the instrument or the Methods Transportability card, depending on whether you selected ToInst or ToCard.							
	Methods on InstUserSizeStoredXL PCR< <pe>&gt;119/9/96Touchdown PCR&lt;<pe>&gt;69/9/96AmpliTaq Gold™&lt;<pe>&gt;99/9/96AmpliCycle Seq&lt;<pe>&lt;89/9/96</pe></pe></pe></pe>							
	Сору	View	User	Sort	Cancel			
	F1	F2	F3	F4	F5			
5	Use the Circu	Ilar Key Pad to	o select a n	nethod and p	ress <b>Copy</b> .			
6	Confirmation	appears wher	1 the metho	od has been s	successfully	copied.		
	Copy Method							
	Method XXXXX was successfully copied to the PC card.							
	Copy an Yes	other met	hod?					
	F1	F2	F3	F4	F5			

### **Connecting and Configuring a Printer**

## **Connecting a Printer** If you connect an optional printer to your 9700 instrument, you can print out a hard copy of the time and temperature parameters for the PCR methods you create.

Connect one end of your printer cable (N805-1326) to the RS-485 serial port on the side panel of the 9700 instrument and connect the other end to the RS-232C interface serial adapter on the rear panel of the printer.

## Configuring the After you have connected the printer cable to the printer, you must configure the instrument for the printer. You can connect the 9700 instrument to any printer with a serial board and the following specifications:

Baud Rate	9600
Parity	NONE
Data Bits	8
Stop Bits	1

See your printer manual for instructions on how to complete any other necessary installation steps.

# 4

# **Running PCR Samples**

### Overview

About This Chapter	This chapter provides information on selecting a me reviewing the history of a run, and what to do when	
In This Chapter	The following topics are covered in this chapter:	
	Торіс	See page
	Selecting a Method	4-2
	Running a Method	4-6
	Reviewing the History of a Run	4-10
	When a Run Completes	4-12

**Note** For information on loading and unloading samples, see the instructions accompanying your interchangeable sample block module.

### Selecting a Method

Introduction	After you have prepared your samples and loaded them in the sample block, you can run a PCR amplification with a new or a stored method.
What Is a Method	A method is a set of instructions in which you specify how the instrument should heat and cool your samples in a PCR thermal profile.
	Methods are stored in the instrument software.
Predefined Methods	The GeneAmp <sup>®</sup> PCR System 9700 supplies eight predefined methods that you can run:
	♦ AmpliCycle <sup>®</sup> Sequencing
	♦ AmpliTaq Gold <sup>®</sup>
	♦ BigDye <sup>™</sup>
	♦ General PCR
	♦ LSM2
	♦ Time Release PCR
	Touchdown PCR
	★ XL PCR
	Each of these methods is stored under the user name < <pe>&gt;. You can edit these methods and store them under a different name, a different user name, or select any one and run it.</pe>
	For a detailed description of each of these pre-coded methods, and how you can use them, see Appendix B, "Supplied Methods."

Selecting a Method If the method you want to run has already been created and stored, you can select it from a list. If the method you want to run has not been created, see Chapter 5, "Creating and Editing Methods."

To select a method:

Step	Action
1	From the Main menu, press Run.
	<b>Note</b> If a PC card is inserted, choose whether to run the method from the PC card or the instrument.
	The Stored Methods screen appears.
	Methods on InstUser Size Storedexp001lisa 10 06/23/96exp002lisa 15 06/25/96
	Start View User Sort Cancel
	F1 F2 F3 F4 F5
	<b>Note</b> Stored represents the date the method was last saved. In the appropriate case, this column designates the date last used.
2	The units for the <b>Size</b> field are based on a calculation of the complexity and length of a method relative to a maximum size of 1102 size segments for the storage capacity of the instrument.
	If you need help deciding which method to select you can:
	<ul> <li>View method parameters.</li> </ul>
	<ul> <li>Sort methods by different categories.</li> </ul>
	<ul> <li>Search for a method by user name.</li> </ul>
3	Select a method by using the Circular Key Pad to move the highlight box to a method listed on the <b>Stored Methods</b> screen.
	<b>Note</b> You can use the up and down keys as repeat keys for quick scrolling.
4	The top line of the display continuously cycles between the following three lines:
	<ul> <li>Methods on Inst User Size Stored.</li> </ul>
	<ul> <li>Used Mem: xxx methods xxx segments.</li> </ul>
	<b>Note</b> The <b>Used Mem</b> field displays the number of size segments used by all stored methods.
	♦ Free Mem: xxx methods xxx segments.
	<b>Note</b> The <b>Free Mem</b> field displays the number of size segments available to store created methods.
5	Press Start and start running your samples (see "Running a Method" on page 4-6).

Viewing Method Parameters	To view	parameters of a method before running:
	Step	Action
	1	From the Stored Methods screen, press View.
		The <b>View Method</b> screen appears. The screen shows all the parameters of the method you selected.
		2 Hld 3 Tmp 25 Cycles 2 Holds 94.0 94.0 94.0 55.0 10:00 0:30 55.0 72.0 0:30 5:00 4.0 0:30 $\sim$
		Start Method: exp 001 Return
		F1 F2 F3 F4 F5
	2	After reviewing PCR and post-PCR parameters of a stored method, you can:
		Press Start to start the method.
		Press Return and return to the Stored Methods screen.
		Note You cannot edit parameters from the View Method screen.

Searching for You can find any method that has been stored under a user name. Methods

To search for a method:

Step	Action			
1	From the Stored Methods screen, press User.			
	Select which user's methods to view by taking the follow	ing action.		
	If you want to	Then		
	list all the methods currently stored on the instrument	Press All.		
	display the methods stored under that user's name	Press Accept.		
	Note You cannot add, delete, or modify a user name f	rom this screen.		
2	Making a selection returns you to the <b>Stored Methods</b> so the methods of the user you selected.	reen which now disp		

Sorting Methods If you have a large number of stored methods, you can sort them by name, date last used, date stored and size.

To sort methods:

Step	Action	
1	From the Stored Methods screen, press S	ort.
	The sorting criteria screen appears.	
	Sort Methods By: Method name	
	Date last u	
	Date stored Method size	
	Accept	Cancel
	F1 F2 F3	F4 F5
2	Use the up and down Circular Key Pad to	select the type of sort.
	The following table describes the sort me	hods:
	Choose this item	To sort
	Method name	Methods alphabetically.
	Date last used	Methods chronologically in descending
		order by date of use.
		The last method which ran or was stored is listed first.
	Date stored	Methods chronologically by date
	<b>Note</b> Uses the most recent title and	stored.
	date, between date last used and the data stored.	The last method stored is listed first.
	Method size	Methods in increasing order by the
		amount of memory used to store each method.
		The largest size method is listed first.
3	Press Accept to accept a selection.	
	This returns you to the Stored Methods so	
	sorted according to your selection in step	2.

### **Running a Method**

 $Starting \ a \ Run \quad \mbox{After choosing a method, follow these steps to start a run.}$ 

To start a run:


Step	Action	
1	From the Main me	enu, press <b>Run</b> to display stored methods.
2	Use the up and down Circular Key Pad to select the method.	
3	a. Press Start.	
	The Select M	ethod Options screen appears.
	b. Choose from t	he following options:
	96-Well Gold Sample Block	Select <b>9600</b> or <b>MAX</b> mode first, and then select <b>volume</b> .
	Module	Select Method Options Reaction Volume: 50 µL Ramp speed: 9600
		Enter a value from 5 to 50 µL Start 9600 MAX Cancel
		F1 F2 F3 F4 F5
	96-Well Aluminum	Select <b>9600</b> or <b>STD</b> mode first, and then select <b>volume</b> .
	Sample Block Module	Select Method Options Reaction Volume: 50 µL Ramp speed: 9600
		Enter a value from 5 to 50 µL Start 9600 STD Cancel
		F1 F2 F3 F4 F5
4	of your reactions	action volume differs from your reaction volume, enter the volume in the <b>Reaction Volume</b> field. me range is dependent on the mode selected.
	For	The range is
	9600 mode	5 to 50 µL.
	MAX/STD mod	e 5 to 100 μL.
	<b>Note</b> Each met Press the <b>CE</b> key	hod remembers the reaction volume. to clear an entry.

To start a run: (continued)

Step	Action					
5	Press Start to	start a run				
	If the tempera heating" appe		heated cover	is less than <i>'</i>	103 °C, this	message, "Cover is
	Ple	ease wai	t. Cover	is heatin	ng	
	Tł	ne run w		re: 65°C h when the hes 103°C		
					Cancel	
	F1	F2	F3	F4	F5	
6	When the hea method you s			C, the Run Ti	me screen o	lisplays and the
	For a descrip	tion of this s	creen, see "/	About the Rui	n Time Scre	en" on page 4-7.
	Sample 5		cle 1 of	25	HOT	
	55.0 10	4.0 94.0 :00 0:30	$\neg \dots \dots \dots$	72.0 72.0	_	
	Pause			Info		
	Fl	F2	F3	F4	F5	

## Screen

About the Run Time You can use the Run Time screen to chart progress at any time during the run. The Run Time screen displays the executing segment, and the next segment to execute.



Use the Run Time screen for	For more information, see page
Viewing Method Information	4-8
Pausing a Run	4-8
Stopping a Run Before It Completes	4-9

Viewing MethodYou can view the Method Information screen during a run by pressing Info. PressInformationReturn to return to the Run Time screen.

02:32 PM	In	formati	on	55.2°C
	ed at 0	1:32:30	PCR PM, 01/01 PM, 01/0	
Reaction	vol: 50	µL Ramp	Speed: 9 Return	600
Fl	F2	F3	F4	F5

Pausing a Run You can manually pause a run for a ten minute period of time during a run by pressing Pause (Figure 4-1). If you want to specify a different period of time for a pause, see "Setting the Pause Time Out" on page 3-10.

Press Resume to resume running a method before a pause expires.

**Note** You can pause a run in order to add a reagent. Do not to touch the sample block or the bottom of the heated cover during a pause.




Figure 4-1 Pause screen

Your samples remain at the temperature of the instrument when you pause the run. The time remaining in a pause appears at the bottom of the screen in minutes:seconds format. It decrements to zero, and the paused run resumes at the point where you paused it.

### Stopping a Run Before It Completes

Stopping a Run The following procedure describes how to stop a run before it completes.

To stop a run before it completes:

Step	Action
1	Press the <b>Stop</b> key.
	The Stop confirmation screen appears.
	Sample 50.1 Confirm Stop HOT
	Press STOP to abort.
	Press Resume to continue.
	Resume
	F1 F2 F3 F4 F5
	The run pauses for a pre-programmed period of time. When the pause time expires, the run aborts.
	Press <b>Resume</b> to resume.
	<b>Note</b> The stop and pause times are configured by setting the <b>Pause Time Out</b> . See "Setting the Pause Time Out" on page 3-10.
2	Press the <b>Stop</b> key again.
	This stops the run and the End of Run screen appears.
	11:30 AM End of Run 25.1°C
	Method:exp001
	Run aborted at 11:30:05 AM 01/01/00. Length of run is 01:34:25.
	Hist
	F1 F2 F3 F4 F5
	If any errors occur during a run the following message appears, "Exception
	occurred, check history file." Press <b>HIST</b> to review the history file.
3	Press Exit to return to the Main menu.

### Reviewing the History of a Run

### How to Review the History of a Run Vou can read a record of the events and errors that occurred during a run by reviewing the history file. The instrument stores the history file until it is overwritten by the next method used.

To review the history:

Step	Action			
1	To display the History File screen:			
	From the Utilities screen, press More, and then press Hist.			
	From the End of Run screen, press Hist.			
	History of method exp002 User: lisa Reaction volume: 50 µL Run started at 02:30:45 PM, 01/01/00. Run aborted at 02:50:42 PM, 01/01/00. Length of run 00:19:57 Ramp speed: 9600 No exceptions PageDn Print Return			
	F1 F2 F3 F4 F5			
2	Press PageUp to move up through the record, or PageDn to move down.			
3	Press Print to print the record.			

History Formats The following table lists the history line formats.

Pre-PCR hold	<exception> in Pre-PCR xx Setpt xx</exception>
PCR segment	<exception> in Cycle xx Setpt xx Repxx</exception>
Any other hold	<exception> in Hold xx Setpt xx</exception>

History File Records The following table lists the history file records.

Record	Description	Туре
History of method xxxxxxxxxxxx User xxxxx Reaction volume xxx µL Run started at hh:mm:ss am mm/dd/yy Run ended at hh:mm:ss am mm/dd/yy Length of run hh:mm:ss RampSpeed: 9600	This header record is always created. <b>Note</b> If you stopped a method before it completed running, then the message, Run ended at, will be Run aborted at	Report
Power failure in Cycle xx at Setpt xx. Power failed at hh:mm:ss am for hh:mm:ss. Run resumed at hh:mm:ss am	There was a power failure during a specified point in a cycle. The message, for >18, indicates that the power was off for more than 18 hours.	Report
Drift error in Cycle xx Setpt xx Repxx. Temperature drifted x.x°c from setpt	Block drift error. The block has drifted ± 2 °C from set point during the hold segment of a run.	Report
Cover error in Cycle xx Setpt xx Repxx. Heated cover at xx.x°c	Heated cover drift error. The cover has drifted ± 5 °C from 105 °C anytime during the run.	Report
Sensor error in Cycle xx Setpt xx Repxx. Block sensor failure.	Block sensor failure.	Fatal error. Call for service.
Sensor error in Cycle xx Setpt xx Repxx. Cover sensor failure.	Heated cover sensor failure.	Fatal error. Call for service.
Setpt error in Cycle xx Setpt xx Repxx. Could not reach xx.x in hh:mm:ss	This setpoint error is only logged for setpoints above 15 °C. The limit is 5 times the normal ramping time.	Fatal error. Call for service.
Program pause in Cycle xx Setpt xx Rep xxx Method paused at xx °C for hh:mm:ss	A programmed pause was encountered.	Report
Manual pause in Cycle xx Setpt xx Rep xxx Method paused at xx °C for hh:mm:ss	You paused the run.	Report


### When a Run Completes

End of Run Screen At completion of a run, the instrument beeps and the End of Run screen appears. From the End of Run screen you can perform the same functions as you can from the Stop Run screen.

11:30 A	М	End of R	un	25.1°C
Method: Run com Length	pleted a	nt 11:30:0 .s 01:34:2	)5 AM, 25.	01/01/00.
Hist				Exit
F1	F2	F3	F4	F5

If you have not yet stored the method, you must store it before exiting, or you will lose the settings. The Store soft key appears if the method has not yet been stored.

Method Not Stored If you attempt to exit the End of Run screen before storing the method, the Method Not Screen Stored screen displays.



The following table lists the actions you can take.

lf you	Then press
want to store the method	Store.
do not want to store the method	Exit.
	<b>Note</b> The <b>Method Not Stored</b> screen appears for a few seconds before the Exit key is recognized. This prevents you from losing a newly created method.
	You return to the top level screen and the method you created is not saved.

# Creating and Editing Methods

### Overview

About This Chapter	how to work with stored methods.		
In This Chapter			
	Торіс	See page	
	Adding or Changing Users	5-2	
	Creating Methods	5-6	
	Modifying Cycling	5-11	
	Printing a Method	5-17	
	Editing or Deleting Methods	5-18	

### **Adding or Changing Users**

**Introduction** The GeneAmp<sup>®</sup> PCR System 9700 stores methods by user's names. You can add up to 19 different user names to the instrument. Once you've added your name to a list of users, and stored a method under that name, you can run the method at any time by selecting it from the Stored Methods screen (See "Selecting a Method" on page 4-3).

Use the User function on the Main menu to add new users or edit existing user names. The name you add or the name you select from a list of existing user names becomes the current user name. All new methods that you create are stored by default under the current user name.

Adding a New User You add a new user name by entering an alphanumeric name on the User Name screen.

To add a new user name:

Step	Action			
1	From the Main menu, press User.			
	The <b>Select User Name</b> screen appears. A list of names of all users who have been added to the instrument displays in a $4 \times 5$ matrix.			
	Select User Name			
	< <pre>&lt;<pre>&lt; description of the second seco</pre></pre>			
	$\begin{array}{c ccccccccccccccccccccccccccccccccccc$			
2	Press <b>New</b> to add a new name to the above list.			
	The <b>User Name</b> screen appears.			
	<b>Note</b> The blank space after the z is used to insert the blank spaces.			
	abcdefghi       jklmnopqrj       User Name       .,-+/():=			
	Use ENTER key to select a character.			
	AcceptBackspCancel			
	F1 F2 F3 F4 F5			

To add a new user name: (continued)

Step	Action			
3	In the <b>User Name</b> field, enter an alphanumeric name up to six characters in length. You can take the following action:			
	If you want to	Then		
	select a character in the list shown in the upper right portions of the screen	use the Circular Key Pads.		
	put the alphabetic character in the User Name field	press Enter.		
	enter the numbers directly into the User Name field	use the numeric keys.		
	go back one space and remove a	hold down the Soft key and go back multiple spaces.		
	single character			
4	single character Press Accept to accept a name. The Sec You can take the following action:			
4	Press <b>Accept</b> to accept a name. The <b>Sec</b> You can take the following action:	surity Code screen appears.		
4	Press Accept to accept a name. The Sec You can take the following action:	Then		

Protecting Methods You can protect methods and prevent other users from accidentally overwriting or deleting them by entering a Personal Identification Number (PIN#) on the Security Code screen.

User Nam	me: hank		umber: ction:	None Unlocked
Press PIN # to create a #. Then you set protection to Locked to prevent methods from being overwritten or deleted.				
Accept	Name	PIN#		Cancel
F1	F2	F3	F4	F5

The following table lists the two levels of protection.

If a	Then other users cannot
user has entered a PIN #	edit that user's name without knowing the PIN #.
method is locked	delete/overwrite the method.

Follow the procedure below to protect a method.

To protect a method:

Step	Action				
1	Press PIN #.				
	The New PIN Number screen appears.				
	Create a PIN Number				
	Your PIN number protects the access to your user name and protection level Enter a PIN number. New PIN #: XXXX				
	[Accept] [Cancel]				
	F1 F2 F3 F4 F5				
2	In the New PIN # field, use the numeric keys and type in a four-digit PIN.				
3	Press Enter.				
	The PIN Confirmation screen appears.				
	Configure DTM Number				
	Confirm PIN Number Your PIN number protects the access to				
	your user name and protection level				
	Enter a PIN number again. PIN #: XXXX Press Accept to confirm your PIN #.				
	Accept Co confirm your PIN #.				
	F1 F2 F3 F4 F5				
4	Confirm your PIN by typing your four-digit PIN in the Confirm PIN # field.				
5	Press Enter.				
	The Protection Status screen appears.				
	Username: hank PIN number: XXXX Protection: Unlocked				
	Press PIN # to create a #. Then you set				
	protection to Locked to prevent methods				
	from being overwritten or deleted. Accept Name PIN# Lock Cancel				
	F1 F2 F3 F4 F5				
6	Press Lock to lock your method.				
	This toggles between a Locked and Unlocked state. The <b>Protection</b> field displays the status of the method.				
7	Press Accept when you have entered a PIN you want to keep. The Select User Name screen appears.				
	The new name you entered should now display on the screen.				
	Press <b>Cancel</b> to cancel your entry and return to the previous screen.				
L	· · ·				

### Changing a UserIf you know the personal identification number for a user name, you can use the<br/>Circular Key Pad to select that name, and change it.

To change a user name:

Step	Action	
1	From the Main menu, press User.	
	The Select User Name screen appears.	
2	Use the Circular Key Pad to select the name you want to change.	
3	Press Edit.	
	If a PIN has been previously entered, the Security Check screen appears.	
4	Type in the four-digit PIN of the user name you selected.	
5	Press Name to enter a new user name.	
6	Press the CE key to clear the previous name.	
7	Enter a new user name.	
8	Press Accept.	
	The Security Code screen appears again.	

## Deleting a UserIf there aren't any methods stored under a user name, you can delete that name from<br/>the Select User Name screen.


To delete a user name:

Step	Action	
1	Use the Circular Key Pad from the Select User Name screen to select a user name.	
2	Press Delete to delete the name.	
	This removes the name from the Select User Name screen and allows you to add a new name to the instrument.	

### **Creating Methods**

## About the Default The GeneAmp PCR System 9700 comes with a default PCR thermal profile called a method. The create screen displays this default method. For information on displaying the create screen, see "Displaying the Create Methods Screen" on page 5-7.

You can run the default method shown above, or use it as a template to create a customized method.



**Basic Parameters** To create a method, you need to define the following four basic parameters:

- Temperature Control Parameters
- Pre-PCR Holds
- PCR Parameters
- Post-PCR Holds

For detailed information about these parameters refer to the table below:

For information about	See page
Entering Temperature Control Parameters	5-8
Defining Pre-PCR Holds	5-9
Defining PCR Parameters	5-9
Defining Post-PCR Holds	5-10

Displaying the Create Methods Screen

**Displaying the** Follow the steps below to display the Create Methods screen.

To display the Create Methods screen

Step	Action	
1	From the Main menu, press Create.	The create screen appears.
	94.0 94.0	rcles 2 Holds 72.0 72.0 0:30 7:00 $4.0More CancelF4 F5$
2	From the <b>Create</b> screen, you can take the following action:	
	If you want to	Then
	start running the default method	press Start.
	store the method under a user's name	press Store.
	print a record of the method parameters	press <b>Print</b> .
	enter temperature control parameters	enter the information on the create screen and create a new method.
	display more functions for modifying methods	press More. Note The More function only displays when you select a time or temperature parameter. For detailed information on using the More function to modify cycles, see "Modifying Cycling" on page 5-11.
	return to the previous screen	press Cancel.

Entering When you enter temperature control parameters, you define values for parameters in each of the three segments of a method: pre-PCR, PCR, and post-PCR.

**Control Parameters** 

To enter temperature control parameters:

Step	Action
1	On the Create screen, select a field.
	When you first display the Create screen, the HId field is highlighted.
2	Use the numeric keys to enter values.
3	Press <b>Enter</b> to accept a value. The next field is then selected in the order shown in Figure 5-1.




Figure 5-1 Order of advancement of the highlight box

Number	Description
1	Number of pre- PCR holds
2	Pre-PCR temperature parameter
3	Pre-PCR hold time parameter
4	Number of PCR segment temperatures
5	Number of cycles for the PCR segment
6	PCR segment temperature parameter
7	PCR segment time parameter
8	PCR segment temperature parameter
9	PCR segment time parameter
10	PCR segment temperature parameter
11	PCR segment time parameter
12	Number of post-PCR holds
13	Post-PCR temperature parameter
14	Post-PCR hold time parameter
15	Post-PCR temperature parameter
16	Post-PCR hold time parameter

#### Defining Pre-PCR Holds The Hld field on the Create screen defines the number of holds for the pre-PCR segment of your method. One (1) is the typical setting for most PCR amplifications and is the default value for the Hld field.

Generally, pre-PCR holds define a temperature and hold-time sufficient to denature any endonucleases or exonucleases that may contaminate your prepared samples.

To define pre-PCR holds:

Step	Action	
1	On the <b>Create</b> screen, in the <b>HId</b> field, type in the number of pre-PCR holds for your method.	
	<b>Note</b> You can enter 0 in this field to delete the pre-PCR hold segment from your method.	
2	Create the first temperature parameter:	
	a. Press Enter.	
	b. Type in a temperature value between 4.0 °C and 99.9 °C.	
3	Create the first hold-time parameter:	
	a. Press Enter.	
	b. Type in a hold-time value between 00:00 and 98:59 (minutes:seconds).	
4	Enter information for the next pre-PCR hold:	
	a. Press Enter.	
	<ul> <li>Repeat step 2 and step 3 until you have hold-time and temperature values for each of the pre-PCR hold parameters you defined in step 1.</li> </ul>	

## Defining PCR<br/>ParametersThe Tmp field on the Create screen defines the number of temperature control<br/>parameters in the PCR cycling segment of your method. Three temperature PCR is<br/>the typical setting for many PCR amplifications:

- Template denaturation
- Primer annealing
- Primer extension

To define the PCR parameters:

Step	Action	
1	On the <b>Create</b> screen, in the <b>Tmp</b> field, type in the number of temperature control parameters (2–6) you want for the PCR segment of your method.	
	<b>Note</b> Specifying only the minimum number of PCR cycles required for analysis will minimize the chance that unwanted targets will amplify competitively.	
2	Enter the number of cycles you want the method to run:	
	a. Press Enter to select the Cycles field.	
	b. In the <b>Cycles</b> field, type in the number of cycles (from 2–99).	
	Note Twenty-five cycles is the default setting.	
3	Create the first temperature parameter:	
	a. Press Enter.	
	b. Type in a temperature value between 4.0 °C and 99.9 °C.	

To define the PCR parameters: (continued)

Step	Action	
4	Create the first hold-time parameter:	
	a. Press Enter.	
	b. Type in a hold-time value between 00:00 and 98:59 (minutes:seconds).	
5	Repeat step 3 and step 4 until you have time and temperature values for each of the segment temperature control parameters you defined in step 1.	

Defining Post-PCR<br/>HoldsOn the Create screen, the Holds field defines the number of temperature control<br/>parameters in the post-PCR segment of your method.

The post-PCR incubation temperature and hold time parameters define how to hold your samples at a specified temperature until you are ready to analyze them.

Note If the idle state setpoint, or the last hold of the Method are below 15 °C, then the heated cover will automatically set to 50 °C.

### **Post-PCR Parameter Settings**

Typical Post-PCR parameter settings:

Temperature	Time (min:sec)	Use For
72 °C	7:00	Complete extension of all amplicons
72 °C	99:59 (×)	AmpErase <sup>™</sup> applications
4 °C	99:59 (×)	General storage

Follow the procedure below to define Post-PCR Holds

To define post-PCR holds:

Step	Action	
1	Select the Holds field.	
2	In the Holds field, type in the number of post-PCR steps for your method.	
3	Press Enter to select the first post-PCR temperature parameter.	
4	Type in a temperature value between 4.0 °C and 99.9 °C.	
5	Press Enter to select the first post-PCR hold time parameter.	
6	Type in a hold time value between 00:00 and 98:59 ( <i>min:sec</i> ).	
	<b>Note</b> The hold time $\times$ indicates a hold that lasts indefinitely. You can enter an $\times$ hold time, by typing a hold time value of 99:00 or greater.	
7	Press Enter. This selects the next temperature parameter.	
8	Repeat step 4 through step 7 until you have time and temperature values for each of the post-PCR hold parameters you defined in step 2.	

### **Modifying Cycling**

**Introduction** In addition to customizing values for PCR temperature control parameters, you can use the More function on the create screen and access cycle modification functions that allow you to:

- Auto-increment/decrement time and temperature parameters.
- Modify up-ramp and down-ramp rates in the cycling segment of a method.
- Insert holds, cycles, and programmed pauses.
- Delete temperature control parameters.

The time or temperature parameter you select on the create screen, determines which modification function you can access when you press More. Different modification functions are available depending on whether you select a temperature control parameter in the pre-PCR segment, the PCR segment, or the post-PCR segment of a method.

### Changing Temperature Control Parameters

Using the AutoX function, you can automatically increase or decrease the value for any PCR segment parameter by a fixed amount every cycle.

**Note** This feature is particularly useful towards the end of the amplification process since the amount of PCR product, available to be extended, increases with the number of cycles while the amount of available enzyme remains constant.

To automatically change temperature control parameters:

Step	Action	
1	From the <b>Create</b> screen, use the Circular Key Pad to select a time or temperature parameter in the PCR segment.	
2	Press More.	
	The Modify screen appears.	
3	Press Modify.	
	The Select Modification screen appears.	
4	Press AutoX. The AutoX screen appears.	
	2 Pre-PCR 3 Tmp 25 Cycles 2 Holds +0.0 +0.0 +0.0 +1.0 -0:30 * Cancel F1 F2 F3 F4 F5 Note If you have inserted a programmed pause, the AutoX screen displays the	
	pause, but you cannot modify it from the AutoX screen.	
5	Select the PCR time or temperature parameter that you want modified when you run your method.	
	<b>Note</b> From the AutoX screen, you cannot modify the number of parameters in each segment or the number of cycles.	


To automatically change temperature control parameters: (continued)

Step	Action				
6	You can take the following action:				
	If you want to	Then press			
	increase the value every cycle	+ (plus sign).			
	(a plus sign displays in the current field)				
	decrease the value every cycle	- (minus sign).			
	(a minus sign displays in the current field)				
	An asterisk * appears on method screens for parameter	ers that have been m	nodified.		
7	Press Accept to accept all entries on the AutoX screen	٦.			

### Modifying Ramp Rates

The ramp time is the time it takes the instrument to change from one temperature to another. Using functions accessible from the Modify screen, you can modify the up-ramp and down-ramp rates of the instrument by defining it as a percentage of the temperature's maximum rate of increase. The default maximum up-ramp and down-ramp rates is 100%.

To modify ramp rates:



To modify ramp rates: (continued)

Step	Action		
5	Modify the up-ramp and down-ramp rates you selected by definin percentage of the maximum of 100%:	g them as a	
	If you want to	Then press	
	increase the ramp rate by 10% up to a maximum of 100%	the Faster key.	
	decrease the ramp rate by 10% from 100% to 10%, and by 5% from 10% to 5%	the Slower key	
	<b>Note</b> You can also use the numeric keys to enter a value that defines the percentage by which you want to decrease the ramp rate for each cycle of the method. You can enter values between 5 and 95, or 100.		
	For ramp rates less than 100%, an asterisk * appears next to modified ramps. The asterisk remains beneath the modified temperature parameter to remind you that the method has been modified.		
	<b>Note</b> The following message appears if you enter a numeric va range of acceptable values, "Valid range is 5 to 95 and 100."	lue outside the	


### **Inserting Holds** Use the Insert function to insert holds and cycles into your method, and program pauses that the instrument automatically inserts into your method as it runs.

To insert holds:

Step	Action
1	From the <b>Create</b> screen, use the Circular Key Pad to select a time or temperature parameter to the right of where you want to insert a hold.
2	Press More.
	Depending on the parameter you select in step 1, one of three screens displays from which you can access the insert function.
3	Press Insert. The Insert screen appears.
	$\begin{array}{ c c c c c c c c c c c c c c c c c c c$
4	Press Hold to insert a hold of 4.0 °C for 30 seconds to the left of the parameter you selected in step 1.
5	Type in a value for the hold temperature.
6	Type in a value for the hold time.
7	Press More to return to the Create screen. Your modified method now appears.

### Deleting Holds Follow the procedure below to delete holds.

To delete holds:



Inserting Cycles Follow the procedure below to insert cycles.

To insert cycles:

Step	Action		
1	From the <b>Create</b> screen, use the Circular Key Pad to select a time or temperature parameter to the right of where you want to insert a cycle.		
2	Press More.		
	Depending on the parameter you select in step 1, one of three screens displays from which you can access the insert function.		
3	Press Insert. The Insert screen appears.		
	$ \begin{array}{ c c c c c c c c c c c c c c c c c c c$		
	F1 F2 F3 F4 F5		
4	Press <b>Cycle</b> to insert a cycle to the left of the segment you selected in step 2.		
	<b>Note</b> You can delete a cycle by entering 0 in the <b>Tmp</b> field.		

Inserting Follow the procedure below to insert a programmed pause.

Programmed Pauses

To insert a programmed pause:

Step	Action
1	From the <b>Create</b> screen, use the Circular Key Pad to select a PCR segment time or temperature parameter where you want to insert a programmed pause.
2	Press More.
3	Press Insert. The Insert screen appears.

To insert a programmed pause: (continued)

Step	Action			
4	From the Insert screen, press Pause . The Programmed Pause screen appears.			
	Define Programmed Pause			
	Start first pause at cycle 1 of 25. Pause every 25 cycles for 00:30. Beep during the pause? No			
	Accept			
	F1 F2 F3 F4 F5			
5	In the <b>Start first Pause at Cycle</b> field, type in the cycle number where you want the method to first pause.			
6	In the <b>Pause Every</b> field, type in the pause frequency in cycles.			
	The pause frequency specifies the number of cycles that will run between each pause.			
7	In the <b>Cycles For</b> field, type in the length of the pause in minutes:seconds (00:01–98:59) format.			
8	In the Beep During The Pause? field, press Yes or No.			
9	Press Accept to accept the pause information on the screen.			
	The word <b>Pause</b> now displays to the right of the incubation step where you programmed the pause.			
	Note You can only insert one pause in each cycle.			

EditingIf you have inserted a programmed pause in your method, you can edit the<br/>parameters for the pause at any time.

To edit programmed pauses:

Step	Action
1	From the <b>Create</b> screen, use the Circular Key Pad to highlight the word <b>Pause</b> . The <b>Edit</b> soft key appears.
2	Press Edit to access the programmed pause screen.
	From this screen, you can
	<ul> <li>Change any of the pause parameters, or</li> </ul>
	<ul> <li>Use the Circular Key Pad to select the pause time parameter on the screen, and edit it by entering a different time.</li> </ul>
	After you have entered all modifications to the customized method you are creating, you should store the method before running it (see "Naming and Storing Methods" below).

Naming and Storing Naming and storing completes the creation of the method. Methods

To name and	store	methods:
-------------	-------	----------

Step	Action		
1	From the Create screen, press Store.		
2	You can take the following action:		
	If you want to The	en	
		ss the e <b>pt</b> key.	
	The default method name is expxxx where xxx= a number from 0 to 999.		
	rename the method go t	o step 3.	
3	From the Store screen, press Method. The Method Name screen appe	ars.	
	<b>Note</b> The blank space after the letter Z is used to insert blank space:	6.	
	abcdefghi jklmnopqrj		
	Method Name exp001 stuvwxyz .,-+/():=		
	Use ENTER key to select a character.		
	AcceptBackspCancel		
	F1 F2 F3 F4 F5		

To name and store methods: (continued)

Step	Action			
4	In the <b>Me</b> name.	thod Name field, follow these s	teps to enter a 1–16 character alphanumeric	
	Step	Action		
	a.	Use the Circular Key Pad to select a character in the list shown in the upper right portion of the screen.		
	b.	After selecting a character, press Enter to place the character in the <b>Method Name</b> field.		
	C.	If you want to	Then	
		enter a number	press the appropriate number key.	
		go back one space	press the Backsp key.	
		clear the method name	press the CE key.	
5	Press Accept after you have entered a method name.			
		e screen appears again. The m ethod Name field.	ethod name you entered should now display	
6	From the Store screen, press Accept.			
	This stores the method under the name you entered.			

### Printing a Method

Introduction	parame	ave configured your instrument for a printer, you can print a record of the ters in a method. For more information on configuring a printer, see "Setting Parameters" on page 3-8.
Printing a Method	To print	a method:
	Step	Action
	1	Access the <b>Create</b> or <b>Edit</b> screen. From the <b>Create</b> or <b>Edit</b> screen you can print a copy of the parameters for the method displayed on the screen.
	2	Press Print.

### **Editing or Deleting Methods**

**Introduction** After you create a method, you can edit its parameters, and store the method by the same name, or change its name. At some time, you may also want to delete a method if you are no longer using it. You can access all editing functions from the Main menu. The delete screen is accessed through the utility menu.

Editing a Method The following procedure describes how to edit a method.

To edit a method:

Step	Action				
1	From the Main menu, press Edit.				
	<b>Note</b> If an (SRAM) PC C editing a method on the PC	ard is detected in the Card slot, you v C Card or the instrument.	will get a choice of		
	The top line of the display	continuously cycles between the follo	wing three lines:		
	Methods on Inst Use [or on PC card]	er Size Stored [or Last Used]			
	Used Mem: xxx methods xxx segments				
	Free Mem: xxx methods xxx segments				
	The following table describes these fields.				
	Field	Description			
	units for the Size field	Based on a calculation of the complexity and length of a method relative to a maximum size of 1102 size segments for the storage capacity of the instrument.			
	Used Mem field	Displays the number of segments used by all stored methods.			
	Free Mem field	Displays the number of segments a created methods.	available to store		
2	Select one of the methods follows:	displayed on the screen, or select an	other method as		
	If you want to		Then press		
	view the parameters of a method before making a selection		the View key.		
	(Refer to "Viewing Method				
	search for a method by user name		the User key.		
	(Refer to "Searching for Methods" on page 4-4.)				
	sort methods by different criteria		the Sort key.		
	(Refer to "Sorting Methods" on page 4-5.)				

To edit a method: (continued)

Step	Action	
3	Press Edit after selecting a method. The Edit screen appears. $ \begin{array}{c ccccccccccccccccccccccccccccccccccc$	
	F1 F2 F3 F4 F5	
4	Choose a temperature or time parameter within a PCR segment.	
5	Edit temperature control or time parameters. <b>Note</b> Editing parameters on the Create screen involves the same tasks and uses the same key combinations as you use when creating a method. The same functions for modifying methods are also available.	
6	From the Edit screen, press Store to store the method.	

### **Deleting a Method** The following procedure describes how to delete a method.

To delete a method:

Step	Action				
1	From the <b>Utilities</b> screen, screen appears.	press <b>More</b> and then press <b>Delete</b> . The I	Delete Method		
	Methods on Ins	ete Method st User Size Stored lisa 15 01/01/00			
	exp001 Press Yes to del Yes				
	F1 F2	F3 F4 F5			
	<b>Note</b> If a PC Card is demethod on the PC Card o	tected in the Card slot, you will get a cho r the instrument.	pice of deleting a		
		continuously cycles between the followi ser Size Stored [or Last Used]	ng three lines:		
	Used Mem: xxx methods xxx segments				
	Free Mem: xxx methods xxx segments				
	The following table describes these fields.				
	Field	Description			
	Units for the Size field	Based on a calculation of the complex a method relative to a maximum size of segments for the storage capacity of the	of 1102 size		
	Used Mem Field	Displays the number of segments used by all stored methods.			
	Free Mem Field	Displays the number of segments available to store created methods.			
2	Select one of the methods follows:	s displayed on the screen, or select anot	her method as		
	If you want to		Then press		
	view method parameters		the View key.		
	(Refer to "Viewing Method Information" on page 4-8.)				
	sort methods by different criteria (Refer to "Sorting Methods" on the <b>Sort</b> key. page 4-5.)				
3	Press <b>Delete</b> . The Delete Confirmation screen appears.				
	If the method is protected, is correct.	, enter a four-digit PIN and press Accept	when the number		

### To delete a method: (continued)

Step	Action		
4	Press <b>Yes</b> to confirm the deletion. This deletes the method and returns you to the Delete screen.		
	<b>Note</b> Even after you delete the last method stored under a User name, the name is removed from the instrument. To delete the name, see "Deleting a User Name" on page 5-5.		



## **Converting Hold Times**


### Overview

About This Chapter	This chapter provides information about setting hold times for the GeneAmp <sup>®</sup> PCR System 9700 compared to the DNA Thermal Cycler or DNA Thermal Cycler 480. This information is based upon differences in how the instruments heat and cool samples. <b>Note</b> Protocols using sample volumes between 5 $\mu$ L and 50 $\mu$ L, developed on either the GeneAmp <sup>®</sup> PCR System 2400 or 9600, may be transported to the GeneAmp PCR System 9700 without change when using the 9600 emulation mode.		
In This Chapter	The following topics are covered in this chapter:		
	Торіс	See page	
	About Setting Hold Times	6-2	
	Guidelines for Converting Hold Times	6-3	
	Guidelines for Converting Hold Times	6-3	

### **About Setting Hold Times**

**Introduction** The hold times specified for the GeneAmp PCR System 9700 are shorter than those used for the Step-Cycle or ThermoCycle files on the DNA Thermal Cycler and the DNA Thermal Cycler 480. This is because the DNA Thermal Cycler and the DNA Thermal Cycler 480 starts counting the hold time when the block reaches a temperature one degree before the target temperature, while the GeneAmp PCR System 9700 starts counting the hold time when the samples reach a temperature one degree before the target temperature.

Hold Time As shown in the figures below, since the block reaches the target temperature before the sample, the programmed hold time on the DNA Thermal Cycler and the DNA Thermal Cycler 480 must include enough time for the samples to reach the target temperature. A hold time of one minute or greater is required on the DNA Thermal Cycler and the DNA Thermal Cycler and the DNA Thermal Cycler 480 for samples to reach the target temperature. On the GeneAmp PCR System 9700, hold times of less than one minute are generally used.



### **Guidelines for Converting Hold Times**

Introduction	This section describes guidelines on how to convert hold times for the DNA Thermal
	Cycler or the DNA Thermal Cycler 480 to hold times for the GeneAmp PCR
	System 9700 when using the 9600 mode.

New Hold Times On or down ramp.

What to Base the The following table lists what to base the new hold time on if you are using an up ramp

If you are using	Then base the new hold time	For more information see
an up ramp	on the change in the temperature required to reach the next target temperature.	Table 6-1 on page 6-3.
a down ramp	on the starting temperature of the ramp and the change in temperature required to reach the next temperature.	Table 6-2 on page 6-5.

Setting Up Ramps The following table lists the process of setting the up ramp temperature.

Step	Action	
1	Determine the change in temperature required to reach the next target temperature and round this value off to the closest value found in Table 6-1.	
2	Subtract the number of seconds indicated from the hold time used for the DNA Thermal Cycler or the DNA Thermal Cycler 480.	
	The result is the hold time to use for the GeneAmp PCR System 9700.	
	<b>Note</b> The typical hold time is 10 to 15 seconds for denaturation.	


### Table of Conversions The following table lists the times for converting up ramp hold times greater than one minute from the DNA Thermal Cycler and the DNA Thermal Cycler 480 to the GeneAmp PCR System 9700.

 Table 6-1
 Converting up ramp hold times

$\Delta$ Temp (°C)	Seconds to subtract from DNA Thermal Cycler or DNA Thermal Cycler 480 hold times (>1 min.)
10°	38 sec.
15°	45 sec.
20°	49 sec.
30°	54 sec.
40°	55 sec,
50°	57 sec.
60°	57 sec.

**Up Ramp Example** In this example, the temperature was increased by 17 °C. This value was rounded to 15 °C. According to Table 6-1, subtract 45 seconds from the hold time on the DNA Thermal Cycler or the DNA Thermal Cycler 480, resulting in a new hold time of 35 seconds (see the figure below).

**Note** If methods developed on the GeneAmp PCR System 9700 will be used on the DNA Thermal Cycler or the DNA Thermal Cycler 480, you can also use Table 6-1 to convert the hold times. Add the indicated times instead of subtracting them.



Setting Down Ramps The following table lists the process of setting the down ramp temperature.

Step	Action	
1	Determine the change in temperature required to reach the next target temperature and round this value off to the closest value found in Table 6-2.	
2	Based on the starting temperature of the ramp, determine the number of seconds to subtract from the DNA Thermal Cycler or DNA Thermal Cycler 480 hold time to arrive at the new hold time.	
	<b>Note</b> The typical hold time is 10 to 15 seconds for annealing.	

## Table of ConversionsThe following table lists the times for converting down ramp hold times greater than<br/>one minute from the DNA Thermal Cycler or the DNA Thermal Cycler 480 to the<br/>GeneAmp PCR System 9700.


	Starting	g temperature of ra	mp (°C)
A Temp (°C)	35-55°	56-75°	<b>76-96</b> °
10°	39 sec.	39 sec.	40 sec.
15°	45 sec.	47 sec.	48 sec.
20°	49 sec.	52 sec.	53 sec.
30°	—	56 sec.	59 sec.
40°	_	57 sec.	62 sec.
50°	_	—	62 sec.
60°	_	_	60 sec.

Table 6-2 Converting down ramp hold times

### Down Ramp Example

In the example below, the temperature was decreased by 22 °C. This value was rounded to 20 °C. According to Table 6-2, with a starting temperature of 94 °C, we should subtract 53 seconds from the hold time on the DNA Thermal Cycler or the DNA Thermal Cycler 480, resulting in a new hold time of 27 seconds.

**Note** If the methods developed on the GeneAmp PCR System 9700 will be used on the DNA Thermal Cycler or the DNA Thermal Cycler 480, use Table 6-2 to convert the hold times. Add the indicated times instead of subtracting them.



# 7

## **Routine Maintenance**

### **Overview**

About This Chapter This chapter describes how to perform routine maintenance on the GeneAmp® PCR System 9700.

**CAUTION** Do not remove the instrument cover. There are no components inside the GeneAmp PCR System 9700 that you can safely service yourself. If you suspect a problem, contact an Applied Biosystems Technical Support Representative.

### In This Chapter The following topics are covered in this chapter:


Торіс	See page
Changing the External Fuses	7-2
Changing Fuses for Single-Line Operation	7-3
Changing Fuses for Double-Line Operation	7-4

### **Changing the External Fuses**

Introduction **CAUTION** FIRE HAZARD. For continued protection against the risk of fire, replace fuses only with listed and certified fuses of the same type and rating as those currently in the instrument.

All instruments have factory installed fuses. However, if you should ever need to change a fuse, or if you need to convert from the factory configured single-line service to double-line service, follow the instructions that follow in this chapter.

Power Entry Module The Power Entry Module is located at the rear of the left side of the instrument.



### Single-Line Operation

Changing Fuses for There is one Type T 250 V fuse accessible from the side of the instrument. You need a flat-tip screwdriver to open the door of the Power Entry Module and change the fuse.

To change the fuse for single-line operation:

Step	Action
1	Turn off the instrument and disconnect the power cord from the side of the instrument.
	<b>ADANGER</b> ELECTRICAL SHOCK HAZARD. Severe electrical shock, which could cause physical injury or death, can result from working on an instrument when the high voltage power supply is operating. To avoid electrical shock, disconnect the power supply to the instrument, unplug the power cord, and wait at least 1 minute before working on the instrument.
	The fuse is located in the Power Entry Module.
2	Insert a small flat-tip screwdriver into the slot in the upper portion of the power entry module, and open the door.
3	Pull the fuse compartment out. There is one fuse in the fuse compartment.
4	Pull out the fuse from the back of the fuse compartment and replace it with either:
	♦ an 8 Amp Type T (P/N 0999-1683) or,
	♦ a Slow-Blow 250 V fuse (P/N 0998-1643).
5	Place the fuse compartment back into the Power Entry Module and close the door.
	Press the door until it locks in place.
6	Connect the instrument power cord.

Changing Fuses for<br/>Double-LineInstruments configured for double-line operation have two 8 amp Type T 250 V<br/>5x20 mm fuses (P/N 0999-1683).OperationTo change the fuses for double-line instruments:

Step	Action
1	Turn off the instrument and disconnect the power cord from the side of the instrument.
	<b>ADANGER</b> ELECTRICAL SHOCK HAZARD. Severe electrical shock, which could cause physical injury or death, can result from working on an instrument when the high voltage power supply is operating. To avoid electrical shock, disconnect the power supply to the instrument, unplug the power cord, and wait at least 1 minute before working on the instrument.
	The fuses are located in the Power Entry Module.
2	Insert a small flat-tip screwdriver into the slot in the upper portion of the power entry module, and open the door.
3	Pull the fuse compartment out. There are two fuses in the fuse compartment.
4	Pull out the fuse from the back of the fuse compartment and replace the blown fuse with one 8 amp Type T 250 V fuse.
5	Place the fuse compartment back into the Power Entry Module and close the door.
	Press the door until it locks in place.
6	Connect the instrument power cord.

# 8

## Troubleshooting

### Overview

About This Chapter	This chapter describes instrument problems you may these problems, and any display screen messages y the GeneAmp <sup>®</sup> PCR System 9700.	
In This Chapter	The following topics are covered in this chapter:	
	Торіс	See page
	If There Is a Power Failure	8-2
	Display Screen Error Messages	8-3
	Troubleshooting Information	8-6

### If There Is a Power Failure

Introduction	An automated restart function allows for power failures and safe continuation of a PCR run after resumption of power.			
During a Power Failure				
	<b>Note</b> If the power is off for 15 seconds or longer and fails during execution of a cycle then the cycle currently running will restart. If the power fails while executing a hold, or approaching a hold, then that hold temperature will restart from the beginning.			
	Note If the power failure lasts longer than 18 hours, the Resume will not occur.			
	Do the following in a power failure:			
	Step Action			
	1	Restart or continue the PCR experiment.		
		The instrument determines what temperature was being approached, or was holding.		
	Upon resumption of power, it will go to that temperature and countdown the tim remaining in the hold as soon as the temperature is within the specified clock s limits.			
	2	Incubate the samples until you can continue the experiment.		
	3	Enter a record for any power outage in the history file.		

### **Display Screen Error Messages**

Error MessagesRefer to the following table for a description of error messages, and recommended<br/>actions that you should take.

### Table 8-1Error Messages

Message	Description	<b>Recommended Action</b>	
Battery RAM version number lost	This error is generated when the battery RAM has been lost and re-initialized.	Call Technical Support.	
Block Calibration initialized	ock Calibration initialized Software or hardware failure.		
Block Calibration reset to	System error.	Call Technical Support.	
default	Block data reset to defaults.		
Block isn't configured	The instrument has defaulted to a 96-well configuration.	Call Technical Support.	
Block initialized	Block module has been re-initialized.	Call Technical Support.	
Block version unknown, update firmware	The calibration data in the block is not recognized by the firmware.	Upgrade the firmware.	
Bus Error	System error.	Call Technical Support.	
Calibration battery RAM initialized	Calibration lost. Instrument may not perform to specification.	Call Technical Support.	
Can only enter an infinity hold at end	A method can only have a HOLD segment with an infinity hold as the last segment in a method.	Assign finite time segments to holds within a method.	
	This message occurs when you try to enter an infinity hold segment in the middle of a method.		
Can't allocate timer	System error.	Call Technical Support.	
Delete your methods first	User tried to delete a user name that has methods stored under it.	Delete or transfer the associated methods before deleting a user name.	
Enter a name or CANCEL	You did not enter at least one character on the User Name screen before pressing the Accept key.	Enter the user name to which the desired method is assigned.	
Enter oligo sequence	Incomplete TmCalc data.	Enter a value in the P1P2 fields of the TmCalc.	
Enter user and method names or CANCEL	You did not enter a user name and a method name before storing a method.	Specify the method name and choose a user to store a method.	
FATAL – Block shut off by hardware	Block thermal runaway.	Call Technical Support.	
FATAL – Block thermal runway	Fatal error.	Turn off system.	
		Call Technical Support.	
FATAL – Cover shut off by hardware	Heated cover thermal runaway.	Call Technical Support.	
FATAL – Heat sink is too hot	Ambient conditions may be too warm.	Call Technical Support.	
FATAL – Heat sink sensor failure	System error.	Call Technical Support.	
FATAL – Heated cover thermal	Fatal error.	Turn off system.	
runaway		Call Technical Support.	

### Table 8-1 Error Messages (continued)

Message	Description	Recommended Action
FATAL – Sample block sensor	Fatal error.	Turn off system.
failure		Call Technical Support.
Fatal – Stack Overflow	A warning or error message that displays which task stack overflowed.	Call Technical Support.
	The warning message is issued when the stack has reached within 10% of overflowing.	
Field is full	You tried to enter more data in a field than the field size allows.	Reenter data within the specifications of the field.
Heated cover sensor failure	The heated cover sensor failed.	Call Technical Support.
Infinity hold not allowed in cycle	A method can have an infinity hold segment as the last segment in the method.	Use finite values for cycle segments within the method.
	This message occurs when you tried to enter an infinity time in a CYCLE segment.	method.
Invalid password/pin#	You entered an incorrect PIN#.	Enter the correct PIN#.
LCD screen timed-out	System error.	Call a Technical Support.
	Display screen and firmware have a faulty connection.	
List of user names is full	The maximum number of users has been entered into the system.	Delete unused user names.
Maximum of 6 segments allowed	You tried to insert more than six temperature control parameters into a hold or cycle.	Do not assign more than six hold or cycle parameters to a method.
Method battery RAM initialized	Stored methods have been reset due to hardware or software failure. Not all methods may be lost.	Check method directory. Call Technical Support.
Method requires at least one segment	You deleted all temperature control parameters in a method.	Review and correct the method to include the
	A method must have at least one time and temperature parameter.	temperature parameter(s).
No seconds in time field	You did not include seconds in the time field.	Include seconds when entering the time.
Not enough method memory	This error occurs:	Determine how much
left	<ul> <li>When you attempt to exceed the limit of 137 methods.</li> </ul>	storage memory is available on the instrument or PC card.
	<ul> <li>When you attempt to store or create a new method which is larger than the available storage space.</li> </ul>	<ul> <li>Delete or store rarely used methods elsewhere.</li> </ul>
Not implemented yet	The feature is not implemented in the current firmware version.	Upgrade firmware when the new version is available from Applied Biosystems.
PC card and Flash do not verify	The PC (upgrade) card and instrument memory do not match.	Call Technical Support.
	Firmware upgrade unsuccessful.	
PC card does not contain valid data	The PC card being used to upgrade the instrument does not contain a valid program.	Call Technical Support.

### Table 8-1 Error Messages (continued)

Message	Description	<b>Recommended Action</b>
Preferences battery RAM initialized	User configuration has been reset due to software error.	Call Technical Support.
Printer not responding	The printer has been disconnected or is off line.	Check printer connections and power switch.
Remove infinity hold first	A method can have an infinity hold segment as the last segment in the method.	Add segments prior to the post-PCR infinity hold.
	This message occurs when a user tries to add a segment after one which contains an infinity hold.	
Setpoint could not be reached	<ul> <li>The instrument could not reach a temperature parameter set by the user.</li> </ul>	Call Technical Support.
	The unit has a Peltier or power amplifier failure.	
	<ul> <li>Ambient conditions may be out of recommended range.</li> </ul>	
SYSTEM ERROR invalid pointer	System error.	Call Technical Support.
Tm temperature out of range	Tm out of range.	Check input value and retry.
		Call Technical Support.
User name already defined	You entered a user name that already exists.	Do not duplicate user names.
WARNING: Block version	Some data in the block is unrecognized.	Upgrade firmware.
unknown	Instrument operation will not be effected.	
Watchdog timeout	Software failure.	Call Technical Support.
Write to block failed	Information written to the Interchangeable Sample Block Module has failed.	Call Technical Support.
Write to default block failed	System error.	Call Technical Support.
	Write to memory in block failed.	
Write to default Xicor failed	System error.	Call Technical Support.
	Write to memory in the block failed.	
Write to Xicor failed	Information written to the Interchangeable Sample Block Module has failed.	Call Technical Support.
Valid range is	You entered a number out of range.	Reenter a value within the
	The message include the valid range limits.	parameters of the field.

### **Troubleshooting Information**

Troubleshooting<br/>TableRefer to the following table for a description of potential problems, possible causes,<br/>and recommended actions that you should take.

 Table 8-2
 Troubleshooting Information

Problem	Possible Causes	Check and/or Remedy
Control panel not responding	Keypad failure.	Run keypad diagnostic.
		Call Technical Support.
Cooling rate too slow	<ul> <li>Ambient temperature is too warm.</li> <li>Peltier failure.</li> </ul>	<ul> <li>Move instrument to well-ventilated location(15-30°C).</li> </ul>
		<ul> <li>Run rate test diagnostic.</li> </ul>
		Call Technical Support.
Cycling time too long	Peltier failure.	Run cycle test diagnostic
Displayed temperature does not match specified temperature	Instrument may require calibration.	Run the Temperature Verification test.
Heated cover not responsive	Heated cover failure.	Call Technical Support.
Heating rate too slow	<ul> <li>Peltier failure.</li> </ul>	Run Rate Test diagnostic.
		Call Applied Biosystems Technical Support.
Instrument can't reach high or low	Ambient temperature is too warm.	♦ Run Rate Test diagnostic.
temperature range	<ul> <li>Peltier failure.</li> </ul>	<ul> <li>Run Cycle Test diagnostic.</li> </ul>
		Call Technical Support.
Instrument making too much noise	Fan failure.	Check for sidevent obstructions.
No beep	<ul> <li>Run time beeper disabled.</li> <li>Beeper failure.</li> </ul>	Check Run-Time Beep configuration. Call Technical Support.
No screen display	♦ Fuse blown.	Is power switch ON?
No response when you turn the	<ul> <li>Not connected to power source.</li> </ul>	Is power cord connected?
instrument on	<ul> <li>Interchangeable module not installed correctly.</li> </ul>	<ul> <li>Check fuses.</li> </ul>
Printer fails	<ul> <li>Incorrect printer configuration.</li> <li>Incorrect printer cable.</li> </ul>	<ul> <li>Check printer settings: baud rate = 9600, no parity, one stop bit, eight data bits.</li> </ul>
		<ul> <li>Purchase Applied Biosystems printer cable.</li> </ul>
Instrument cooling fan does not make	♦ Fuse blown.	Is power switch ON?
whirring sound	<ul> <li>Not connected to power source.</li> </ul>	Is power cord connected?
	<ul> <li>Interchangeable module not installed correctly.</li> </ul>	<ul> <li>Check fuses.</li> </ul>

# Instrument **Specifications**



### **Overview**

About This This appendix describes the dimensions, power, and electrical specifications for the Appendix GeneAmp® PCR System 9700 system, including the control panel, sample temperature information, and printer specifications.

In This Appendix The following topics are covered in this appendix:

Торіс	See page
System Specifications	A-2
Control Panel Specifications	A-3
Sample Temperature Information	A-3
Printer Specifications	A-3

### **System Specifications**

**Dimensions** The following tables list the footprint and the weight of the instrument with the 96–Well or 60–Well sample block module.

#### Footprint (With Sample Block Module Installed)

Height	26 cm (10 in)
Width	28 cm (11 in)
Depth	41 cm (16 in)

Note You must provide sufficient space around the instrument for unrestricted air circulation.

### Weight

Base Module	8.6 kg (19 lbs)
96-Well Gold/Aluminum Sample Block Modules	3.2 kg (7 lbs)

**Note** See the sample block module user's manual for physical information on a module.

**Power** There is one version of the instrument. The power requirements of the instrument under various power configurations are:

VAC ~100/120	8 AMP T (5x20 mm) or	50/60 Hz
	8 AMP Slow Blow (3 AB)	Use 250 V fuses
VAC ~220/230/240	8 AMP T (5x20 mm)	Max Power 725 VA

**Electrical IMPORTANT** You must be able to disconnect the main power supply to the instrument immediately if necessary.

The following table specifies the electrical operating range for the instrument in various parts of the world. Select appropriate fuse configuration based on the voltage used.

**WARNING** In Japan, the unit must have a dedicated 220-volt outlet! The unit will not operate properly with a 100-volt outlet.

Location	Voltage (VAC) <sup>a</sup>	Frequency	Amperage (A) Nominal
Japan	220 ±10%	50/60 Hz ±1%	3.16
USA/Canada	120 ±10%	50/60 Hz ±1%	4.20
Europe (pre-1992)	220 ±10%	50/60 Hz ±1%	3.16
EC	230 ±10%	50/60 Hz ±1%	3.14
UK (pre-1992)	240 +6%/–10%	50/60 Hz ±1%	3.12
Australia	240 +6%/–10%	50/60 Hz ±1%	3.12

a. Acceptable AC line voltage tolerances: 100, 120, 220, 230  $\pm$ 10%; 240 VAC +6%/–10%, 50/60 Hz  $\pm$  1%.

**Note** The Volt-Amp number for this instrument is 725Volt Amps.

#### **Control Panel Specifications**

Display Screen	The display screen is a 7 x 40 character display with a graphics mode of 60 x 240 pix resolution.				
Keys	The instrument control panel consists of a display screen and 22 keys. The keys are:				
	♦ Function keys				
	♦ Arrow keys				
	♦ Stop key				
	♦ Enter key				
	♦ 10 number keypad				

#### **Sample Temperature Information**

**Temperature** The following table lists sample temperature information.

Note Sample temperatures are displayed in degrees Celsius to the nearest 0.1 °C.

Sample Temperature Range	4.0 to 99.9 °C.
Temperature Calibration	Traceable to National Institute of Standards and Technology (NIST).

#### **Printer Specifications**

Serial Interface The instrument can use any printer with a serial (RS-232C) interface board with the Board Specifications following parameters.

Baud Rate	9600
Parity	NONE
Data Bits	8
Stop Bits	1
	Parity Data Bits

Cable Part Number Connect the printer to the Instrument port with printer cable part number N805-1326.

## B

### **Supplied Methods**

#### Overview

About this Appendix	GeneAmp <sup>®</sup> PCR System 9700.	
In This Appendix		
	Торіс	See page
	About the Methods	B-2
	AmpliCycle Sequencing	B-2
	AmpliTaq Gold	B-2
	BigDye	B-3
	General PCR	B-3
	LMS2	B-4
	Time Release PCR	B-4
	Touchdown PCR	B-5
	XL PCR	B-5

#### **About the Methods**


**Introduction** The instrument supplies you with eight pre-coded methods stored under the user name <<pre>repe>>.

The eight pre-coded methods are:

- AmpliCycle<sup>®</sup> Sequencing
- ♦ AmpliTaq Gold<sup>®</sup>
- ♦ BigDye<sup>™</sup>
- General PCR
- LSM2
- Time Release PCR
- Touchdown PCR
- ♦ XL PCR


See the following sections for detailed information about the methods.

AmpliCycleCycle sequencing has revolutionized Sanger sequencing of PCR products and otherSequencingDNA templates. The AmpliCycle Sequencing protocol is for 3-temperature cycling and<br/>for achieving clean sequence ladders from femtomole amounts of template.



This cycle sequencing process and the benefits of AmpliTaq<sup>®</sup> DNA Polymerase, CS, are described in the product insert for the AmpliCycle<sup>®</sup> Sequencing Kit (P/N N808-0175).


AmpliTaq Gold The AmpliTaq Gold protocol specifies a 5-minute pre-PCR heat step, required for the activation of AmpliTaq Gold<sup>®</sup> DNA Polymerase. This additional step provides seamless "hot start" PCR and replaces labor intensive methods such as manual hot start or wax bead-mediated hot start techniques.



Utilizing hot start techniques helps to minimize the formation of primer-dimers or non-specific products, thus increasing specificity and sensitivity of PCR.


You can find further information on AmpliTaq Gold DNA Polymerase in the product insert (P/N N808-0241) or at the Applied Biosystems website.

BigDye The BigDye<sup>™</sup> method consists of cycle sequencing parameters for dideoxy (Sanger) terminator sequencing using ABI PRISM<sup>®</sup> BigDye<sup>™</sup> Terminator Cycle Sequencing Ready Reaction Kits (available from Applied Biosystems). It consists of 25-cycle, three-temperature cycle sequencing followed by an infinite hold at 4 °C.




The pre-coded conditions in the BigDye<sup>™</sup> method are optimized for AmpliTaq<sup>®</sup> DNA Polymerase, FS and the cycle sequencing reagents supplied with the BigDye terminator kits. This process is further described in the *ABI PRISM BigDye Terminator Cycle Sequencing Ready Reaction Kits Protocol* (P/N 4303237).

General PCR The General PCR method is a basic one and can be easily modified with both preand post-PCR holds.



LMS2 The GeneAmp PCR System 9700 software includes a pre-coded LMS2 method for ABI PRISM<sup>®</sup> Linkage Mapping Set Version 2 (LMS2) thermal cycling. Linkage Mapping Set 2 employs over 400 fluorescent-labeled PCR primer pairs for analysis of select microsatalite loci from the Généthon human linkage map.<sup>1,2,3</sup> The following figures illustrate the thermal cycling profile for the LMS2 method.



This supplied method consists of an initial hold at 95.0 °C, two sets of three temperature cycles, followed by two additional holds at 72.0 °C and 4.0 °C.

See the ABI PRISM Linkage Mapping Set Version 2 User's Manual (P/N 904999) for comprehensive information on PCR amplification conditions, electrophoresis conditions, detection, and data analysis.

**Time Release PCR** The Time Release PCR method is designed for use with AmpliTaq Gold DNA Polymerase. The enzyme is activated more slowly than with the AmpliTaq Gold method. Here the pre-PCR hold is only 1 minute, and the number of cycles is increased to 40.


1 Hld 3	Tmp 40 (	Cycles 2	Holds	
95.0 9	5.0 55.0	72.0 72	. 0	
1.00 0	:15 0:15	72.0 72 0:30 7:0	<u>, , , , , , , , , , , , , , , , , , , </u>	
		1	$\sqrt{\frac{4.0}{\infty}}$	
Start	Mathad	Timo Pol	ease PCR	Return
Deare	Methou:	TIME VET	ease FCK	Incearing
F1	F2	F3	F4	F5

<sup>1.</sup> Weissenbach, J. *et al.* 1992. A second-generation linkage map of the human genome. *Nature* 359:794–801.

<sup>2.</sup> Gyapay, G., et al. 1994. Généthon Human Genetic Linkage Map. Nature Genet. 7:246-339.


<sup>3.</sup> Dib, C., et al. 1996. Nature 380: 152–154.

**Touchdown PCR** When the optimal annealing temperature is unknown, one strategy, touchdown PCR, incrementally decreases the annealing temperature in early cycles in order to maximize the yield of specific products.



This supplied method has an initial annealing temperature (65 °C) that incrementally decreases by an additional 0.5 °C in each of the first 20 cycles, followed by 10 cycles at 55 °C.

**XL PCR** XL PCR is the protocol specified for amplification of 5 kb-40 kb PCR products, using r*Tth* DNA Polymerase, XL, and unique reaction conditions.



This protocol uses two-temperature cycling (94 °C for 15 seconds; 68 °C for 10 minutes) and invokes a 15 second AutoX (automatic segment extension) for the anneal/extend step in the last 12 cycles.

By providing longer templates, XL PCR complements technologies for rapid, long-range PCR. More complete genes can be amplified in one reaction from known expressed sequences, thus more introns can be crossed. You can use XL PCR for the amplification of the control target, a 20.8 kb product from Lambda DNA, supplied in the kit.

This process is further described in the product insert for the GeneAmp<sup>®</sup> XL PCR Kit (P/N N808-0192).

## Contacting Technical Support

Synthesis

# C

#### **Technical Support**

Contacting	You can contact Applied Biosystems for technical support:				
<b>Technical Support</b>	◆ By e-mail				
	<ul> <li>By telephone or fax</li> </ul>				
	<ul> <li>Through the Applied Biosystems web s</li> </ul>	site			
To Contact Technical Support by E-Mail					
<b>T</b>	Product/Product Area E-mail address				
	Genetic Analysis (DNA Sequencing)     galab@appliedbiosystems.com				
	Sequence Detection Systems (Real-Time PCR) and PCR	pcrlab@appliedbiosystems.com			
	Protein Sequencing, Peptide, and DNA	corelab@appliedbiosystems.com			

Product/Product Area	E-mail address
<ul> <li>Biochromatography (BioCAD<sup>®</sup>, SPRINT<sup>™</sup>, VISION<sup>™</sup>, and INTEGRAL<sup>®</sup> Workstations and POROS<sup>®</sup> Perfusion Chromatography Products)</li> </ul>	tsupport@appliedbiosystems.com
<ul> <li>Expedite<sup>™</sup> 8900 Nucleic Acid Synthesis Systems</li> </ul>	
<ul> <li>MassGenotyping Solution 1<sup>™</sup> (MGS1) Systems</li> </ul>	
<ul> <li>PNA Custom and Synthesis</li> </ul>	
♦ Pioneer <sup>™</sup> Peptide Synthesizers	
<ul> <li>Proteomics Solution 1<sup>™</sup> (PS1) Systems</li> </ul>	
♦ ICAT™ Reagent	
♦ FMAT <sup>™</sup> 8100 HTS Systems	
<ul> <li>Mariner<sup>™</sup> ESI-TOF Mass Spectrometry Workstations</li> </ul>	
<ul> <li>Voyager<sup>™</sup> MALDI-TOF Biospectrometry Workstations</li> </ul>	
<ul> <li>CytoFluor<sup>®</sup> 4000 Fluorescence Plate Reader</li> </ul>	
LC/MS (Applied Biosystems/MDS Sciex)	support@sciex.com
Chemiluminescence (Tropix)	tropix@appliedbiosystems.com

**Telephone or Fax** 

To Contact Technical To contact Applied Biosystems Technical Support in North America, use the telephone Support by or fax numbers in the table below.

Note To schedule a service call for other support needs, or in case of an emergency, dial (North America) 1.800.831.6844, then press 1.

Product/Product Area	Telephone	Fax
ABI PRISM <sup>®</sup> 3700 DNA Analyzer	<b>1.800.831.6844</b> , then press <b>8</b> ª	1.650.638.5981
DNA Synthesis	<b>1.800.831.6844</b> , press <b>2</b> , then press <b>1</b> ª	1.650.638.5981
Fluorescent DNA Sequencing	<b>1.800.831.6844</b> , press <b>2</b> , then press <b>2</b> <sup>a</sup>	1.650.638.5981
Fluorescent Fragment Analysis (including GeneScan <sup>®</sup> applications)	<b>1.800.831.6844</b> , press <b>2</b> , then press <b>3</b> <sup>a</sup>	1.650.638.5981
Integrated Thermal Cyclers (ABI PRISM <sup>®</sup> 877 and Catalyst 800 instruments)	<b>1.800.831.6844</b> , press <b>2</b> , then press <b>4</b> <sup>a</sup>	1.650.638.5981
ABI PRISM <sup>®</sup> 3100 Genetic Analyzer	<b>1.800.831.6844</b> , press <b>2</b> , then press <b>6</b> <sup>a</sup>	1.650.638.5981
Peptide Synthesis (433 and 43x Systems)	<b>1.800.831.6844</b> , press <b>3</b> , then press <b>1</b> <sup>a</sup>	1.650.638.5981
Protein Sequencing (Procise <sup>®</sup> Protein Sequencing Systems)	<b>1.800.831.6844</b> , press <b>3</b> , then press <b>2</b> <sup>a</sup>	1.650.638.5981
Sequence Detection Systems (Real-Time PCR) and PCR	1.800.762.4001, then press:	1.240.453.4613
	1 for PCR <sup>a</sup>	
	2 for TaqMan <sup>®</sup> applications and Sequence Detection Systems including ABI Prism <sup>,</sup> 7700, 7900, and 5700 <sup>a</sup>	
	<b>6</b> for the 6700 Automated Sample Prep System <sup>a</sup> or	
	<b>1.800.831.6844</b> , then press <b>5</b> ª	
<ul> <li>Mariner<sup>™</sup> ESI-TOF Mass Spectrometry Workstations</li> </ul>	<b>1.800.899.5858</b> , press <b>1</b> , then press <b>3</b> <sup>b</sup>	1.508.383.7855
<ul> <li>♦ Voyager<sup>™</sup> MALDI-TOF Biospectrometry Workstations</li> </ul>		
<ul> <li>MassGenotyping Solution 1<sup>™</sup> (MGS1) Systems</li> </ul>		
<ul> <li>Proteomics Solution 1<sup>™</sup> (PS1)</li> <li>Systems</li> </ul>		
<ul> <li>ICAT<sup>™</sup> Reagent</li> </ul>		

Product/Product Area	Telephone	Fax
Biochromatography (BioCAD <sup>®</sup> , SPRINT <sup>™</sup> , VISION <sup>™</sup> , and INTEGRAL <sup>®</sup> Workstations and POROS <sup>®</sup> Perfusion Chromatography Products)	1.800.899.5858, press 1, then press 4 <sup>b</sup>	1.508.383.7855
Expedite™ 8900 Nucleic Acid Synthesis Systems	1.800.899.5858, press 1, then press 5⁵	1.508.383.7855
Pioneer™ Peptide Synthesizers	1.800.899.5858, press 1, then press 5 <sup>b</sup>	1.508.383.7855
PNA Custom and Synthesis	1.800.899.5858, press 1, then press 5 <sup>b</sup>	1.508.383.7855
<ul> <li>◆ FMAT<sup>™</sup> 8100 HTS Systems</li> <li>◆ CytoFluor<sup>®</sup> 4000 Fluorescence Plate Reader</li> </ul>	1.800.899.5858, press 1, then press 6 <sup>b</sup>	1.508.383.7855
Chemiluminescence (Tropix)	<b>1.800.542.2369</b> (U.S. only), or <b>1.781.271.0045</b> °	1.781.275.8581
LC/MS (Applied Biosystems/MDS Sciex)	1.800.952.4716	1.508.383.7899

a. 5:30 AM to 5:00 PM Pacific time.

b. 8:00 AM to 6:00 PM Eastern time.

c. 9:00 AM to 5:00 PM Eastern time.

#### To Contact Technical Support by Telephone or Fax (Outside North America)

To Contact Technical<br/>Support byTo contact Applied Biosystems Technical Support or Field Service outside North<br/>America, use the telephone or fax numbers below.

Region	Telephone	Fax		
Eastern Asia, China, Oceania				
Australia (Scoresby, Victoria)	61 3 9730 8600	61 3 9730 8799		
China (Beijing)	86 10 64106608 or 86 800 8100497	86 10 64106617		
Hong Kong	852 2756 6928	852 2756 6968		
India (New Delhi)	91 11 653 3743/3744	91 11 653 3138		
Korea (Seoul)	82 2 593 6470/6471	82 2 593 6472		
Malaysia (Petaling Jaya)	60 3 79588268	60 3 79549043		
Singapore	65 896 2168	65 896 2147		
Taiwan (Taipei Hsien)	886 2 2358 2838	886 2 2358 2839		
Thailand (Bangkok)	66 2 719 6405	66 2 319 9788		
	Europe	ŀ		
Austria (Wien)	43 (0)1 867 35 75 0	43 (0)1 867 35 75 11		
Belgium	32 (0)2 532 4484	32 (0)2 582 1886		
Denmark (Naerum)	45 45 58 60 00	45 45 58 60 01		
Finland (Espoo)	358 (0)9 251 24 250	358 (0)9 251 24 243		
France (Paris)	33 (0)1 69 59 85 85	33 (0)1 69 59 85 00		
Germany (Weiterstadt)	49 (0)6150 101 0	49 (0)6150 101 101		
Italy (Milano)	39 (0)39 83891	39 (0)39 838 9492		

Region	Telephone	Fax
Norway (Oslo)	47 23 12 06 05	47 23 12 05 75
Portugal (Lisboa)	351.(0)22.605.33.14	351.(0)22.605.33.15
Spain (Tres Cantos)	34.(0)91.806.1210	34.(0)91.806.12.06
Sweden (Stockholm)	46 (0)8 619 4400	46 (0)8 619 4401
Switzerland (Rotkreuz)	41 (0)41 799 7777	41 (0)41 790 0676
The Netherlands (Nieuwerkerk a/d IJssel)	31 (0)180 392400	<b>31 (0)180 392409</b> or <b>31 (0)180 392499</b>
United Kingdom (Warrington, Cheshire)	44 (0)1925 825650	44 (0)1925 282502
European Mar	aged Territories (EMT)	
Africa, English speaking (Johannesburg, South Africa)	27 11 478 0411	27 11 478 0349
Africa, French speaking (Paris, France)	33 1 69 59 85 11	33 1 69 59 85 00
India (New Delhi)	91 11 653 3743	91 11 653 3138
	91 11 653 3744	
Poland, Lithuania, Latvia, and Estonia (Warszawa)	48 22 866 40 10	48 22 866 40 20
For all other EMT countries not listed (Central and southeast Europe, CIS, Middle East, and West Asia)	44 1925 282481	44 1925 282509
	Japan	
Japan (Hacchobori, Chuo-Ku, Tokyo)	81 3 5566 6230	81 3 5566 6507
La	tin America	
Caribbean countries, Mexico, and Central America	52 55 35 3610	52 55 66 2308
Brazil	0 800 704 9004 or 55 11 5070 9654	55 11 5070 9694/95
Argentina	800 666 0096	55 11 5070 9694/95
Chile	1230 020 9102	55 11 5070 9694/95
Uruguay	0004 055 654	55 11 5070 9694/95

Support Through the Applied **Biosystems Web Site** 

To Reach Technical At the Applied Biosystems web site, you can search through frequently asked questions (FAQs) or a solution database, or you can submit a question directly to Technical Support.

#### Search FAQs

To search for FAQs:

Step	Action
1	Go to http://www.appliedbiosystems.com
2	Click SERVICES & SUPPORT at the top of the page, then click Frequently Asked Questions.
3	Click you geographic region for the product area of interest.
4	Follow the instructions under the <b>Frequently Asked Questions</b> section (1) to display a list of FAQs for your area of interest.

#### Search the Solution Database

To search for solutions to problems using the Solution Database:

Step	Action
1	Go to http://www.appliedbiosystems.com
2	Click SERVICES & SUPPORT at the top of the page, then click Frequently Asked Questions.
3	Follow the instructions under the <b>Search the Solution Database</b> section (2) to find a solution to your problem.

#### Submit a Question

To submit a question directly to Technical Support:

1	Go to http://www.appliedbiosystems.com
2	Click SERVICES & SUPPORT at the top of the page, then click Frequently Asked Questions.
3	In the Personal Assistance – E-Mail Support section (3), click Ask Us RIGHT NOW.
4	In the displayed form, enter the requested information and your question, then click <b>Ask Us RIGHT NOW</b> .
	Within 24 to 48 hours, you will receive an e-mail reply to your question from an Applied Biosystems technical expert.

#### To Obtain Technical You can obtain technical documents, such as Applied Biosystems user documents, Documents MSDSs, certificates of analysis, and other related documents for free, 24 hours a day. You can obtain documents:

- By telephone ٠
- Through the Applied Biosystems web site ŧ

#### **Ordering Documents by Telephone**

To order documents by telephone:

1	From the U.S. or Canada, dial <b>1.800.487.6809</b> , or from outside the U.S. and Canada, dial <b>1.858.712.0317</b> .
2	Follow the voice instructions to order documents (for delivery by fax).
	<b>Note</b> There is a limit of five documents per fax request.

#### **Obtaining Documents Through the Web Site**

To view, download, or order documents through the Applied Biosystems web site:

Step	Action
1	Go to http://www.appliedbiosystems.com
2	Click <b>SERVICES &amp; SUPPORT</b> at the top of the page, then click <b>Documents on Demand</b> .
3	In the search form, enter and select search criteria, then click <b>Search</b> at the bottom of the page.
4	In the results screen, do any of the following:
	<ul> <li>Click the pdf icon to view a PDF version of the document.</li> </ul>
	<ul> <li>Right-click the pdf icon, then select Save Target As to download a copy of the PDF file.</li> </ul>
	<ul> <li>Select the Fax check box, then click Deliver Selected Documents Now to have the document faxed to you.</li> </ul>
	<ul> <li>Select the Email check box, then click Deliver Selected Documents Now to have the document (PDF format) e-mailed to you.</li> </ul>
	<b>Note</b> There is a limit of five documents per fax request, but no limit on the number of documents per e-mail request.

To Obtain Customer Training Information

Step	Action
1	Go to http://www.appliedbiosystems.com
2	Click SERVICES & SUPPORT at the top of the page, then click Training.

## Limited Warranty Statement



PE Corporation (NY), through its Applied Biosystems Group ("Applied Biosystems") warrants to the customer that, for the Warranty Period ("the Warranty Period is as stated by Applied Biosystems in the accompanying documentation and begins on the date of shipment to the customer"), the GeneAmp® PCR Instrument System purchased by the customer ("the Instrument") will be free from defects in material and workmanship, and will perform in accordance with the specifications set out in the Instrument's User's Manual. If the Instrument fails to perform in accordance with the specifications during the Warranty Period, Applied Biosystems will replace or repair the Instrument, at a site designated by Applied Biosystems, at no charge to the customer.

The Warranty does not extend to any Instrument or part thereof that has been the subject of misuse, neglect or accident, that has been modified or repaired by any party other than Applied Biosystems or its designees, or that has been used in a manner not in accordance with the instructions contained in the Instrument User's Manual.

Applied Biosystems's obligation under this Warranty is limited to repairs or replacements that Applied Biosystems deems necessary to correct covered defects or failures of which Applied Biosystems is notified prior to expiration of the Warranty Period. All repairs and replacements under this Warranty shall be performed by Applied Biosystems's service engineers or designees.

No agent, employee, or representative of Applied Biosystems has any authority to bind Applied Biosystems to any affirmation, representation, or warranty concerning the Instrument that is not contained in the printed product literature or this Warranty Statement. Any such affirmation, representation, or warranty made by any agent, employee, or representative of Applied Biosystems shall not be binding on Applied Biosystems.

Applied Biosystems shall not be liable for any incidental, special, or consequential loss, damage or expense directly or indirectly arising from the purchase or use of the Instrument. Applied Biosystems makes no warranty whatsoever with regard to products or parts furnished by third parties; such products or parts will be subject to the warranties, if any, of their respective manufacturers.

This Warranty is limited to the original customer and is not transferable.

THIS WARRANTY IS THE SOLE AND EXCLUSIVE WARRANTY AS TO THE INSTRUMENT AND IS IN LIEU OF ANY OTHER EXPRESS OR IMPLIED WARRANTIES, INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE AND OF ANY OTHER OBLIGATION ON THE PART OF APPLIED BIOSYSTEMS.

### Index

Numerics 480 DNA Thermal Cycler converting hold times 6-2 to 6-5

#### A

alphanumeric characters, entering 5-3 AmpliCycle Sequencing 4-2, B-2 about supplied method B-1, B-2, B-3 AmpliTaq Gold 4-2, B-2 about supplied method B-2 AutoX function 5-11 to 5-12

#### B

baud rate, setting 3-11 beeper, run-time, turning on/off 3-10 BigDye 4-2, B-2 about supplied method B-3

#### С

calculator, T<sub>M</sub>, how to use 3-18 Changing 7-2 changing temperature control parameters 5-11 to 5-12 chemical hazard warning 1-6 chemical waste hazard warning 1-6 Circular Key Pad, definition of 2-3 Config, function of 3-8 Configuration Screen, displaying 3-8 Contacting Technical Support C-1 **Control Panel** diagram 2-2 entering numeric values, how to 2-3 selecting a field, how to 2-3 specifications A-3 using the keys 2-3 counting cycles. See cycles CPU board testing 3-16 Create Methods screen, displaying 5-7 custom parameters, setting 3-8 to 3-12 baud rate, setting 3-11 date, setting 3-9 Idle State Setpoint temperature, setting 3-11 Pause Time Out, setting 3-10 printer, enabling/disabling 3-9 screen contrast/saver mode, setting 3-12 time, setting 3-9 custom parameters, settings Configuration Screen, displaying 3-8 run-time beeper, turning on/off 3-10 customer support. See technical support C-1 cycles counting of 4-7 inserting 5-14

cycling, modifying 5-11 to 5-17 AutoX function, using 5-11 to 5-12 cycles, inserting 5-14 deleting holds 5-14 holds inserting 5-13 ramp rates, modifying 5-12 to 5-13

#### D

data bit, serial board specifications 3-11, 3-27 date last used, sorting methods by 4-5 setting 3-9 deleting methods 5-18 to 5-21 user name 5-5 diagnostic (tests), running 3-14 to 3-17 CPU board, testing 3-16 Diagnostic screen, viewing 3-14 Display screen, testing 3-16 electrical components, testing 3-15 keypad, testing 3-17 Diagnostic screen, viewing 3-14 dimensions, instrument A-2 display screen. See diagnostic (tests), running Display softkey, function of 3-15 DNA Thermal Cycler, converting hold times 6-2 to 6-5 Documents on Demand C-6

#### E

editing methods 5-18 to 5-21 electrical components testing 3-15 electrical safety 1-3 to 1-11 electrical shock hazard 1-10 e-mail, address for technical support C-1 End of Run screen, when a run completes 4-12 environment 1-12 error messages, display screen 8-3 to 8-5

#### F

Field Service in North America, contacting C-3 fields, selecting 2-3 fire hazard 1-11 firmware, upgrading 3-19 to 3-24 through serial port 3-19, 3-22 using PCMCIA card 3-23, 3-24 fuses changing for double-line operation 7-4 changing for single-line operation 7-3 converting from single to double-line 3-5 to 3-6 precautions about 1-11 service configuration 3-4

#### G

GeneAmp PCR System 9700 Control Panel about 2-2 specifications A-3 copying methods 3-25 to 3-26 custom parameters, setting 3-8 to 3-12 baud rate, setting 3-11 Configuration Screen, displaying 3-8 date, setting 3-9 Idle State Setpoint temperature, setting 3-11 Pause Time Out, setting 3-10 printer, enabling/disabling 3-9 run-time beeper, tuning on/off 3-10 screen contrast/saver mode, setting 3-12 time, setting 3-9 diagnostic (tests), running CPU board, testing 3-16 Diagnostic screen, viewing 3-14 Display screen, testing 3-16 electrical components, testing 3-15 keypad, testing 3-17 fuses, converting from single to double-line fuse 3-5 to 3-6 laboratory, setting up 3-3 to 3-4 air vents, air circulation 3-3 choosing a location 3-3 fuse service configuration 3-4 input voltage, table of 3-4 interchangeable block module, about 3-4 operating temperature 3-3 what to avoid 3-3 overview 2-2 to 2-5 about the instrument 2-2 control panel keys, using 2-3 control panel, diagram 2-2 entering numeric values, how to 2-3 how differs from 9600 System 2-5 methods, transferring or storing 2-4 selecting a field, how to 2-3 upgrading the firmware 2-4 printer connecting and configuring 3-27 specifications A-3 system firmware upgrading 3-19 to 3-24 system specifications A-2 T<sub>M</sub> calculator, how to use 3-18 turning on instrument 3-7 unpacking 3-2 General PCR 4-2, B-2 General PCR, about supplied method B-3 grounding 1-11

#### Η

Hard soft key, function of 3-14 hazards and warnings 1-3 to 1-8, 1-9 to 1-11 burns, danger of 1-10

electrical shock hazard 1-10 fire hazard 1-11 instrument warnings 1-9 heated cover screen 4-7 heating, of samples 5-8 history of run. reviewing history file records 4-11 history formats 4-10 procedure 4-10 hold time converting from DNA Thermal Cycler 6-2 to 6-5 parameters, defined 2-3 holds deleting 5-14 inserting 5-13 humidity 1-12

#### Ι

Idle State Setpoint temperature, setting 3-11 input voltage, table of 3-4 inserting holds 5-13 Internet address customer training information C-7 Documents on Demand C-6

#### K

keyboard softkey, function of 3-15 keypad, testing 3-17

#### L

labels, safety 1-9
laboratory environmental requirements 1-12
laboratory, setting up 3-3 to 3-4
air vents, air circulation 3-3
choosing a location 3-3
fuse service configuration 3-4
input voltage, table of 3-4
interchangeable block module, about 3-4
operating temperature 3-3
what to avoid 3-3
LCD pixels, turning on and off 3-16
locking/unlocking methods, about 5-4
LSM2 4-2, B-2

#### М

maintenance, routine 1-12 Memory Used field, definition of 4-3, 5-20 Method Not Stored screen, actions if displayed 4-12 methods adding or changing users 5-2 to 5-5 adding users 5-2 to 5-3 changing user name 5-5 deleting user name 5-5 protecting methods 5-3 to 5-4 copying using PCMCIA card 3-25 to 3-26 creating 5-6 to 5-10 basic parameters, to create method 5-6

Create Methods screen, displaying 5-7 default method, about 5-6 PCR holds, defining 5-9 to 5-10 post-PCR holds, defining 5-10 pre-PCR holds, defining 5-9 temperature control parameters, entering 5-8 editing or deleting methods 5-18 to 5-21 modifying cycling 5-11 to 5-17 AutoX function, using 5-11 to 5-12 cycles, inserting 5-14 deleting holds 5-14 naming and storing methods 5-16 to 5-17 programmed pauses 5-14, 5-16 ramp rates, modifying 5-12 to 5-13 printing a method 5-17 reviewing history of run 4-10 to 4-11 history file records 4-11 history formats 4-10 procedure 4-10 running 4-6 to 4-9 method information, viewing 4-8 pausing a run 4-8 Run Time screen, about 4-7 starting a run 4-6 to 4-7 stopping run before completion 4-9 selecting 4-2 to 4-4 method defined 4-2 method parameters, viewing 4-4 methods, searching for 4-4 methods, sorting 4-5 predefined methods, list of 4-2 procedure 4-3 supplied B-2 to B-5 AmpliCycle Sequencing B-2 AmpliTag Gold B-2 BigDye B-1, B-3 BigDye<sup>™</sup> B-1 LMS2 B-1 Touchdown PCR B-5 XLPCR B-5 transferring or storing, about 2-4 when a run completes 4-12 MSDSs 1-7

#### N

naming methods 5-16 to 5-17 numeric keys, definition of 2-3

#### 0

operating precautions 1-12

#### Р

parameters Hold Time 2-3 temperature control 5-11 temperature parameters 2-3 *See Also* custom parameters, setting parity bit 3-11, 3-27 part numbers for ordering accessories 3-2 fuse part number 3-5 Pause Time Out, setting 3-10 pauses, programmed editing 5-16 inserting 5-14 PCMCIA card using to copy methods 3-25 to 3-26 using to upgrade firmware 3-23, 3-24 PCR holds, defining 5-9 to 5-10 PCR samples, running methods running 4-6 to 4-9 selecting 4-2 to 4-4 PCR thermal profile, method defined 4-2 PIN#, entering 5-4 pollution 1-12 post-PCR holds, defining 5-10 power configurations A-2 power failure, what to do 8-2 precautions, operating 1-12 pre-PCR holds, defining 5-9 printer baud rate, setting 3-11 connecting and configuring 3-27 enabling/disabling 3-9 printer port, defining. See baud rate printing method 5-17 printer specifications A-3 programmed pauses editing 5-16 inserting 5-14 protecting methods 5-3 to 5-4

#### R

ramp rates, modifying 5-12 to 5-13 reaction volume entering 4-6 RS-232C, serial interface board 3-27 Run softkey, function of 3-16 run-time beeper, turning on/off 3-10

#### S

safety 1-1 to 1-8 sample temperature range A-3 Sanger sequencing B-2 screen contrast, setting 3-12 screen saver mode, setting 3-12 security code, entering 5-4 selecting fields 2-3 serial interface board specifications 3-11 serial port, upgrading firmware 3-19, 3-22 Size field, definition of 5-20 soft keys definition of 2-3 Sort, function of 4-5, 5-20 SRAM card, using to copy methods 3-25 to 3-26 starting instrument 3-7 stop bit 3-11, 3-27 Stop key, definition of 2-3 storing methods 5-16 to 5-17 system firmware, upgrading 3-19 to 3-24 through serial port 3-19, 3-22 using PCMCIA card 3-23, 3-24 System soft key, function of 3-14

#### T

technical support C-1 to C-7 e-mail address C-1 Internet address C-6 regional sales offices C-4 to C-5 telephone/fax (North America) C-3, C-4 temperature entering values 2-3 Idle State Setpoint, setting 3-11 operating precautions 1-12 operating temperature 3-3 sample temperature information A-3 temperature control parameters entering 5-8 temperature control parameters entering 5-8 using AutoX function 5-11 to 5-12 **Temperature Verification test** part number 3-2 Time Release PCR 4-2, B-2 Time Release PCR, about supplied method B-4 time, setting 3-9 T<sub>M</sub> calculator, how to use 3-18 Tmp field, using to define temperature control parameters 5-9 TmpVer soft key, function of 3-14 ToCard softkey, function of 3-26 ToInst softkey, function of 3-26 Touchdown PCR 4-2, B-2 about supplied method B-5 training obtaining information C-7 troubleshooting 8-1 to 8-6 display screen error messages 8-3 to 8-5 power failure, what to do 8-2 table of information 8-6 turning on instrument 3-7

#### U

unpacking the instrument 3-2 Upgrade soft key, function of 3-14 User, function of 4-4 users 5-2 to 5-5 adding users 5-2 to 5-3 changing user name 5-5 deleting user name 5-5 protecting methods 5-3 to 5-4 Util, function of 3-8

#### W

warning chemical hazard 1-6 chemical waste hazard 1-6 waste disposal 1-8 Win 3.1, downloading firmware using 3-22 Windows 95/Windows NT downloading firmware using 3-21

#### X

XL PCR 4-2, B-2 about supplied method B-5

#### Headquarters

850 Lincoln Centre Drive Foster City, CA 94404 USA Phone: +1 650.638.5800 Toll Free (In North America): +1 800.345.5224 Fax: +1 650.638.5884

#### Worldwide Sales and Support

Applied Biosystems vast distribution and service network, composed of highly trained support and applications personnel, reaches into 150 countries on six continents. For sales office locations and technical support, please call our local office or refer to our web site at www.appliedbiosystems.com or to the Technical Support and Training appendix in this document.

#### www.appliedbiosystems.com



Applera Corporation is committed to providing the world's leading technology and information for life scientists. Applera Corporation consists of the Applied Biosystems and Celera Genomics businesses.

Printed in Singapore 11/2001 Part Number 4303481 Rev. E

an Applera business