
# INSTRUCTIONS FOR USING THE TYPE 50 Ti ROTOR

In Beckman Class B, C, D, F, G, H, Q, and R Preparative Ultracentrifuges





#### **SPECIFICATIONS**

Maximum speed
At $r_{\text{max}}$ (80.8 mm)
·
At $r_{av}$ (59.1 mm)
At $r_{\min}$ (37.4 mm)
k factor at maximum speed
Number of tube cavities
Available tubes see Tables 1 and 2
Nominal dimensions of largest tube <sup>5</sup> / <sub>8</sub> x 3 in.
$(16 \times 76 \text{ mm})$
Nominal tube capacity
Nominal rotor capacity
Approximate acceleration time to maximum speed
(rotor fully loaded) in an L8M ultracentrifuge 5 <sup>1</sup> / <sub>2</sub> min
Approximate deceleration time from maximum speed
(rotor fully loaded) in an L8M ultracentrifuge 4 <sup>1</sup> / <sub>2</sub> min
Weight of fully loaded rotor
Rotor material titanium
Conditions requiring speed reduction see Run Speeds
Conditions requiring speed reduction see Run speeds

<sup>\*</sup> Relative Centrifugal Field (RCF) is the ratio of the centrifugal acceleration at a specified radius and speed ( $r\omega^2$ ) to the standard acceleration of gravity (g) according to the following formula:

$$RCF = \frac{r\omega^2}{g}$$

where r is the radius in millimeters,  $\omega$  is the angular velocity in radians per second  $(2\pi RPM/60)$ , and g is the standard acceleration of gravity (9807mm/s<sup>2</sup>). After substitution:

$$RCF = 1.12 r \left(\frac{RPM}{1000}\right)^2$$

### **DESCRIPTION**

The Type 50 Ti titanium fixed angle rotor, rated for 50 000 rpm in Beckman class B, C, D, F, G, H, Q, and R preparative ultracentrifuges, holds up to twelve 5/8-in diameter tubes at a 26-degree angle from the axis of rotation. This general purpose rotor is useful for many applications: concentrating proteins and subcellular particles and isopycnic banding of DNA. Up to 162 mL of gradient and sample can be centrifuged per run.

The exterior surface of the rotor is finished with black urethane paint for optimal temperature control. The aluminum handle and lid are red-anodized for corrosion resistance. The O-rings in the rotor handle and body are made of Buna N and maintain atmospheric pressure in the rotor during centrifugation. This rotor is protected from exceeding 50 000 rpm in ultracentrifuges with the photoelectric or the mechanical overspeed detection systems; it is equipped with an overspeed disk (see the Supply List) and two mechanical cartridges for this purpose. The Type 50 Ti rotor is warranted at maximum speed for 5000 runs, 10 000 hours of centrifugation, or 5 years, whichever occurs first (see the Warranty).

#### **OPERATION**

NOTE: Specific information about the Type 50 Ti rotor is given here. Information common to this and other rotors is contained in the Rotors and Tubes Manual, LR-IM, which should be used together with this bulletin for complete rotor and accessory operation.

#### **TUBES AND BOTTLES**

Tubes and bottles available for the Type 50 Ti rotor are detailed in Tables 1 and 2. Be sure to observe the maximum speed limits shown, and use *only* the tubes and bottles listed.

#### **Temperature Limits**

Polyallomer and polycarbonate containers should not be centrifuged below 2°C. Do not freeze polyallomer tubes before centrifugation, as they may become brittle and crack. Polyallomer and polycarbonate containers may be centrifuged at temperatures above 20°C, but they should be pretested under anticipated run conditions. Ultra-Clear TM tubes have been tested for use at temperatures between 4 and 20°C. For centrifugation at other temperatures, pretest these tubes as well. Stainless steel tubes can be centrifuged at any temperature.

## **Tube Preparation**

Thinwall polyallomer and Ultra-Clear tubes must be filled as full as possible and capped to prevent tube collapse during centrifugation. Thickwall plastic tubes may be used with or without tube caps. Do not overfill capless tubes. Stainless steel tubes can be filled to any level, but all opposing tubes for a run must be filled to the same level.

Table 1. Available 5/8-in. (16-mm) Diameter Tubes\* for the Type 50 Ti Rotor

	Dod	Required Accessory		Nominal Fill	Maximum
Description	Part Number	Description	Part Number	Volume (mL)	Speed (rpm)
Thinwall Ultra-Clear polyallomer	344085 326814	aluminum cap aluminum cap	330860 330860	13.5 13.5	50 000 50 000
Quick-Seal Ultra-Clear polyallomer polyallomer polyallomer polyallomer	344322 342413 344622 345830 346562	spacer spacer spacer floating spacer floating spacer	342696 342696 344634 345828† 345828†	13.5 13.5 10 6.3 4.2	50 000 50 000 50 000 50 000 50 000
Thickwall polyallomer	355640	none aluminum cap	none 338907	8 <sup>†</sup> 13.5	30 000 50 000
polycarbonate	355630	none aluminum cap	none 338907	8 <sup>†</sup> 13.5	50 000 . 50 000
Polycarbonate bottle	355603 (assembly) 355651 (bottle only)	Noryl cap	355604	10.4 <sup>‡</sup>	50 000
Stainless Steel	301108	aluminum cap	303319	13.5	**

<sup>\*</sup> All tubes are  ${}^{5}/_{8} \times 3$  in (16 x 76 mm) except Quick-Seal tubes 344622 ( ${}^{5}/_{8} \times 2$   ${}^{5}/_{8}$  in.; 16 x 67 mm), 345830 ( ${}^{5}/_{8} \times 1$   ${}^{9}/_{4}$  in.; 16 x 45 mm) and 356562 ( ${}^{6}/_{8} \times 1$   ${}^{1}/_{2}$  in.; 16 x 38 mm). k factors for these shorter tubes are 72, 50 and 44, respectively.

Table 2. Small Open-Topped Tubes Used with Delrin Adapters for the Type 50 Ti Rotor

	Tube			7	Required Maximum	
Dimensions and Volume	Description	Part Number	Required Cap	Adapter	Required Coood I	
<sup>1</sup> / <sub>2</sub> x 2 <sup>1</sup> / <sub>2</sub> in (6.5 mL)	Ultra-Clear polyallomer stainless steel	344088 326820 301099	346256* 346256* 305022	202212	03313 40 000 rpm 133 000 x g 1	120
¹/ <sub>2</sub> x 2¹/ <sub>2</sub> in	thickwall polycarbonate	355645	none	303313		120
(4 mL)	thickwall polyallomer	355644	none			
<sup>1</sup> / <sub>2</sub> x 1 <sup>5</sup> / <sub>8</sub> in (4 mL)	Ultra-Clear	344093	346256*	303402	43 000 rpm 134 000 x g	69
<sup>1</sup> / <sub>2</sub> x 1 <sup>1</sup> / <sub>4</sub> in (3 mL)	Ultra-Clear	344092	346256*	303401	43 000 rpm 126 000 x g	71
<sup>5</sup> / <sub>16</sub> x 1 <sup>15</sup> / <sub>16</sub> (2 mL)	Ultra-Clear	344091	303624	303376	40 000 rpm 136 000 x g	69

<sup>\*</sup>Includes the 344672 neoprene gasket and the 346246 stem.

<sup>†</sup> Floating spacers, part of the g-Max system of tube support, are made of Noryl, a registered trademark of General Electric.

<sup>&</sup>lt;sup>‡</sup> Minimum fill volume is 5 mL.

<sup>\*\*</sup> Read the section on Run Speeds before using stainless steel (SS) tubes.

Use only those caps listed in the tables. Do not use any substitutes. Cap O-rings and gaskets should be dry and free from lubricant during assembly. Aluminum or stainless steel caps should be tightened using a hex driver while the tube is held in the tube vise. (Refer to publication L5-TB-060.)

Polycarbonate bottles may be centrifuged completely filled, or if desired, partially filled with at least 5 mL of liquid. Use only caps listed in the table. Tighten caps by hand. Use the 335381 tool to remove bottles from the rotor.

Quick-Seal® tubes should be filled leaving a small air space at the neck of the tube. Do not leave a large air space. Too much air can cause the tube to deform, which may disrupt gradients or sample. See publications IN-163 and IN-18 for detailed information on the use and care of Quick-Seal tubes.

The g-Max<sup>™</sup> system uses a combination of short bell-top Quick-Seal tubes (see Table 1) and floating spacers (also referred to as g-Max spacers). This provides an advantage over using conventional sleeve-type adapters, in which the tubes are no longer at the maximum radius, resulting in a reduction of g force. In the g-Max system, floating spacers sit on top of the Quick-Seal tubes—thus, there is no reduction of radial distance and, therefore, no reduction of g force. Further, adapters cannot be run at maximum speed. Finally, the shorter pathlength of the tubes permits faster run times for gradient or pelleting separations. For more information on the g-Max system, see publication DS-709.

#### ROTOR PREPARATION

Before using the rotor, be sure it is equipped with the correct overspeed protection device. For instrument classes G, H, and R use the 36-sector overspeed disk. For all other classes, the mechanical cartridges (see the Supply List for part numbers) should be installed in the sides of the rotor base. Be sure the threads in the rotor handle and body are well lubricated with Spinkote I lubricant, and the O-rings are thinly coated with silicone vacuum grease. For runs at temperatures other than room temperature, always refrigerate or warm the rotor beforehand, since titanium is a poor conductor of heat. Tubes placed opposite each other in the rotor cavities should be filled to the same level with the same density liquid. Load the rotor symmetrically.

In Models L2-50/65 ultracentrifuges, use the stablizer level "50" for the Type 50 Ti rotor. In Models L2-65B/75B, use the fourth level (four dots).

#### **RUN TIMES**

The k factor of the rotor is a measure of the rotor's pelleting efficiency. (Beckman has calculated k factors for all of its preparative rotors at maximum rated speed and using full tubes.) The k factor is calculated from the formula:

$$k = \frac{\ln (r_{\text{max}}/r_{\text{min}})}{\omega^2} \times \frac{10^{-10}}{3.6}$$
 (1)

where  $\omega$  is the angular velocity of the rotor in radians per second ( $\omega = 0.105 \text{ x RPM}$ ),  $r_{\text{max}}$  is the maximum radius, and  $r_{\text{min}}$  is the minimum radius.

Use the k factor in the following equation to estimate the run time t (in hours) required to pellet particles of known sedimentation coefficient s (in Svedberg units).

$$t = \frac{k}{s} \tag{2}$$

Run times can be estimated for centrifugation at less than maximum speed by adjusting the k factor as follows:

$$k_{\text{adj}} = k \left( \frac{\text{rated speed of rotor}}{\text{actual run speed}} \right)^2$$
 (3)

Run times can also be estimated from data established in prior experiments if the k factor of the previous rotor is known. For any two rotors<sup>1</sup> a and b:

$$\frac{t_{\mathbf{a}}}{t_{\mathbf{b}}} = \frac{k_{\mathbf{a}}}{k_{\mathbf{b}}} \tag{4}$$

where the k factors have been adjusted for the actual run speed used.

Finally, use equation (1) above to calculate the k factor when the column of liquid is such that the operational  $r_{\max}$  and  $r_{\min}$  are significantly different from the  $r_{\max}$  and  $r_{\min}$  of the rotor (e.g., when using partially filled tubes or adapters). For more information on k factors see the Rotors and Tubes Manual and publication DS-719.

#### **RUN SPEEDS**

The centrifugal force at a given radius in a rotor is a function of the rotor speed. Comparisons of forces between different rotors are made by comparing the rotors' relative centrifugal fields (RCF). When rotational speed is selected so that identical samples are subjected to the same RCF in two different rotors, one may then describe the samples as having been subjected to the same centrifugal force field (refer to Table 3).

<sup>&</sup>lt;sup>1</sup> The method is most accurate when comparing two like rotors, i.e., two vertical tube rotors or two fixed angle rotors.

Table 3. Relative Centrifugal Fields. Entries in this table are calculated from the formula RCR = 1.12 r (RPM/1000)<sup>2</sup> and then rounded to three significant digits.

200000000000000000000000000000000000000	T	10 1111 00 018/01/100		1	
Rotor	Relative Centrifugal Field (x g)				
Speed (rpm)	At r max (80.8 mm)	At r <sub>av</sub> (59.1 mm)	At r <sub>min</sub> (37.4 mm)	k Factor*	
50 000	226 000	165 000	105 000	78	
46 000	191 000	140 000	88 600	92	
45 000	183 000	134 000	84 800	96	
43 000	167 000	122 000	77 500	106	
40 000	145 000	106 000	67 000	122	
37 500	127 000	93 100	58 900	139	
35 000	111 000	81 100	51 300	159	
31 500	89 800	65 700	41 600	197	
30 000	81 400	59 600	37 700	217	
25 000	56 600	41 400	26 200	312	
20 000	36 200	26 500	16 800	488	

Calculated for all Beckman preparative rotors as a measure of the rotor's relative pelleting efficiency in water at 20°C.

Rotor speeds may not be selected in excess of those provided in Tables 1 and 2 for particular tubes and caps. In addition, rotor speeds must be reduced in the following circumstances.

1. For centrifuging nonprecipitating solutions of density greater than 1.2 g/mL in *plastic tubes*, use the following square-root reduction formula to determine the allowable rotor speed.

reduced maximum speed = 
$$50\ 000\ \text{rpm}\sqrt{\frac{1.2\ \text{g/mL}}{\text{density of tube contents}}}$$
 (5)

This speed reduction will protect the rotor from excessive stresses due to the added load.

- 2. When CsCl or other self-forming gradient salt is centrifuged, the square-root reduction formula above usually will not guard against the precipitation of salt crystals. Solid CsCl has a density of 4 g/mL, and if precipitated during centrifugation will alter the density distribution, and therefore sample separation. Figures 1 and 2, together with the description and examples, show how to reduce rotor speed when using CsCl gradients.
- 3. When centrifuging *stainless steel tubes*, the maximum permissible rotor speed depends on the type of cap assembly used and the total load to be centrifuged.

- a. Maximum speeds based on the type of cap assembly are: for all-aluminum caps: 50 000 rpm for aluminum caps with stainless steel stems: 37 500 rpm
  - b. Maximum speeds based on the tube load should be determined from the equation:

reduced maximum speed = 
$$50\ 000\ \text{rpm}\ \sqrt{\frac{16.2\ \text{grams}}{M}}$$
 (6)

where M is the mass in grams of the capped stainless steel tube plus the liquid contents to be centrifuged. NOTE: In all cases, the lower speed determined in a and b is the limiting speed.

#### SELECTING CSCI GRADIENTS

Rotor speed is used to control the shape of a CsCl gradient, but precipitation of CsCl during centrifugation must be prevented. Speed and density combinations that intersect on or below the curves in Figure 1 ensure that CsCl will not precipitate during centrifugation of the Type 50 Ti rotor. Curves are provided at two temperatures: 20°C (black curves) and 4°C (gray curves). Note that for a given initial homogeneous CsCl solution density, the maximum allowable speed increases as the fill volume of solution in the tube decreases.

The curves in Figure 2 show gradient profiles at equilibrium. Each curve was generated for a single rotor speed using the maximum allowable homogeneous CsCl densities (one for each fill volume) that avoid precipitation at that speed and temperature.<sup>2</sup> Figure 2 can be used to approximate the banding positions of sample particles. In general, lower speeds provide better particle resolution, but will require longer run times to achieve particle equilibration.

**NOTE:** The curves in Figures 1 and 2 are for solutions of CsCl salt only. If other salts are present in significant concentrations, the overall CsCl concentration or the rotor speed must be reduced. This prevents precipitation of salts concentrated at the tube bottom.

Solutions can be centrifuged faster (therefore for shorter run times) when a tube is partially filled with gradient and sample. (For centrifuging partially filled thinwall tubes, fill the remainder of the tube with a low-density, immiscible liquid, such as mineral oil. (Do not use an oil overlay in Ultra-Clear tubes.) For example, a *full* tube of a 1.65 g/mL homogeneous CsCl solution can be centrifuged no faster than 40 000 rpm at 20°C (from Figure 1). Figure 2 presents the gradient profile, from 1.485 g/mL at the meniscus to 1.86 g/mL at the tube bottom. The same solution in a 3/4-filled tube can be centrifuged at 20°C at 44 000 rpm.

<sup>&</sup>lt;sup>2</sup>The gradients in Figure 2 can be generated from step or linear gradients, or from homogeneous solutions.





Figure 2. CsCl Gradients at Equilibrium. Centrifugation of homogeneous CsCl solutions at maximum allowable speeds (from Figure 1) results in the gradients presented. Density increases from the top to the bottom of the tube (from 37.4 to 80.8 mm from the axis of rotation).

Each square on the grid represents 0.010 g/mL by 1.27 mm.



Interpolation of Figure 2 between the 43 000 and 46 000 rpm curves to the 3/4-filled level gives the 44 000 rpm curve. The same solution in a 1/2-filled tube can be centrifuged at 46 500 rpm, and in a 1/4-filled tube at 50 000 rpm. Note that the 46 000 rpm and 50 000 rpm curves in Figure 2 are identical for 4°C and 20°C. The 1/4-filled volume at 50 000 rpm yields a gradient distribution of 1.28 to 1.44 g/mL at the tube bottom.

# TYPICAL EXAMPLES OF DETERMINING RUN PARAMETERS Example A: Starting with a homogeneous CsCl solution density (e.g., 1.50 g/mL) and approximate particle densities (e.g., 1.48 and 1.52 g/mL), where will particles band?

- 1. In Figure 1 find the curve that corresponds to the desired run temperature (20°C) and tube fill volume (3/4-filled). The maximum allowable rotor speed is determined from the point where this curve intersects the homogeneous CsCl density (46 000 rpm).
- 2. In Figure 2, sketch a horizontal line corresponding to each particle density.
- 3. Mark the point where each density intersects the curve corresponding to the maximum speed and selected temperature.
- 4. Particles will band at these points along the tube axis.

In this example, particles will band at about 65 and 68 mm from the axis of rotation (about 3 mm of interband separation at the 26-degree tube angle, or  $d_{\theta}$ ). When the tube is held upright, there will be about 7 mm of interband separation ( $d_{\text{up}} = d_{\theta} / \sin \theta$ ). About 16 hours of centrifugation is required for particle equilibration. If a full tube is used instead of a 3/4-filled tube, the maximum allowable speed for a 1.50 g/mL solution is 45 000 rpm and about 20 hours of centrifugation will be required.

# Example B: Knowing particle densities (e.g., 1.70 and 1.65 g/mL), how do you achieve good separation?

- 1. In Figure 2, sketch a horizontal line corresponding to each particle density.
- 2. Select the curve at the desired temperature (4°C) and tube volume (full) that gives good separation.
- 3. Note the speed indicated along the curve (30 000 rpm).
- 4. From Figure 1, determine the maximum allowable homogeneous CsCl density that corresponds to the selected run parameters (temperature, speed, and fill volume), 1.695 g/mL. These parameters will provide the particle separation selected in Step 2.

In this example, particles will band at about 50 and 61 mm from the axis of rotation (about 11 mm of interband separation at the tube angle). When the tube is held upright, there will be about 25 mm of interband separation.

#### MAINTENANCE

Routinely inspect the overspeed device. If the disk is scratched, damaged, or missing, replace it. If a mechanical cartridge has tripped, replace both cartridges, as the other has also fatigued. Instructions for installing a new device are given in the Rotors and Tubes Manual. Do not use sharp tools on the rotor. Store the rotor in a dry environment (not in the instrument) with the lid removed.

Silicone vacuum grease should be routinely applied to the O-rings in the handle and rotor body. Replace the O-rings every 6 months or whenever worn or damaged. Keep the rotor threads well lubricated with Spinkote lubricant. Refer to the Rotors and Tubes Manual for the chemical resistance of rotor and tube materials. Contact your Beckman Field Service Representative for information about the rotor repair program and the Field Rotor Inspection Program.

#### RETURNING A ROTOR

Before returning a rotor or accessory for any reason, prior permission (a Returned Goods Authorization form) must be obtained from Beckman Instruments. This RGA form may be obtained from your local Sales Office. It should contain the following information:

- serial number,
- a history of use (approximate frequency of use),
- the reason for the return,
- the original purchase order number, billing number, and shipping number, if possible,
- the name and phone number of the person to be notified upon receipt of the rotor or accessory at the factory, and
- the name and phone number of the person to be notified about repair costs, etc.

To protect our personnel, it is the customer's responsibility to ensure that the parts are free from pathogens and/or radioactivity. Sterilization and decontamination must be done before returning the parts. Smaller items (such as tubes, bottles, etc.) should be enclosed in a sealed plastic bag.


All parts must be accompanied by a note, plainly visible on the outside of the box or bag, stating that they are safe to handle and that they are not contaminated with pathogens or radioactivity. Failure to attach this notification will result in return or disposal of the items without review of the reported problem.

# SUPPLY LIST

# REPLACEMENT ROTOR SUPPLIES

See the Rotors, Tubes, and Accessories catalog (PL-174) for detailed information on reordering supplies. For your convenience, a partial list is given below.

nst is given below.
Type 50 Ti Rotor
and spacers see Tables 1 and 2
Rotor handle
Handle O-ring
Rotor O-ring
Overspeed disk (50 000 rpm)
Mechanical overspeed dvice (50 000 rpm)
OTHER SUPPLIES
Solution 555
Rotor Cleaning Kit
Silicone vacuum grease
Spinkote lubricant
Tube-sealing kit (with 60-Hz sealer)
Tube-sealing kit (with 50-Hz sealer)
Tube-sealer racks for Quick-Seal tubes:
tubes 344322 and 342413
tube 345830
tube 344622 344641
Tube Topper kit (60 Hz)
Tube Topper kit (50 Hz)
Tube rack for 5/8-in. diameter tubes (24 places) 348123
Tube removal tool for capped tubes
Tube-cap vise
Hex driver for 7/ <sub>16</sub> -in. (22 mm) cap nuts
Hex driver for 5/ <sub>16</sub> -in. (7.8 mm) cap nuts
(for tube cap 303624)
Tube removal tool for polycarbonate
bottles with Noryl <sup>4</sup> caps
Domes with though caps



36-sector Overspeed Disk (50 000 rpm)

<sup>&</sup>lt;sup>4</sup>Registered trademark of General Electric.