## OPTION 011 DATA OUTPUT


### INSTALLATION INSTRUCTIONS FOR OPTION 011 DATA OUTPUT FOR AE BALANCES

#### General information

By using an Option 011 Data Output, a METTLER TOLEDO AE balance can be connected to one or several of the following instruments.

- METTLER TOLEDO CL instruments, METTLER TOLEDO GA40 Printer (with transfer button at the 011 data output)
- Non-METTLER TOLEDO instruments with active 20 mA data interfaces
- Non-METTLER TOLEDO instruments with RS232C data interfaces (computers, terminals)

The METTLER TOLEDO brochure "The Option 011 Data Output for METTLER TOLEDO AE balances (CL/RS232C)", order No. 701091, provides information about the data traffic with non-METTLER TOLEDO instruments, as well as hardware definitions.


#### Opening the instrument

This is required only if it is necessary to select a configuration that is different from the way the instrument has been set at the factory.

(Factory setting corresponds to the standard METTLER TOLEDO CL configuration: 2400 baud/even parity/send on transfer. Cable shielding RS232C: not grounded).

- To open the instrument, loosen the four screws in the corners (arrows).

#### Configuring the data output

- Baud rate: according to the setting of the instrument(s) to be connected, set with coding plug (1) according to instructions on P/C board (factory setting 2400 baud).
- Operating modes: Before being shipped from the factory, all positions of the operating mode switch (2) are set to OFF.


Position	ON	OFF
1)	parity disable	parity enable
2)	odd parity [with 1) OFF] space [with 1) ON]	even parity [with 1) OFF] mark [with 1) ON]
3)	no function	no function
4)	send continuous	send on transfer

- Cable shielding (RS232C): Plug in coding plug (3) at center and right (4) → cable shielding is at housing and ground.

Plug in coding plug (3) at center and left (5) (or use no coding plug at all) → cable shielding is floating (not grounded).


Subject to technical changes and to the availability of the accessories supplied with the instruments.  
 Technische Änderungen und Änderungen im Lieferumfang des Zubehörs vorbehalten.  
 Sous réserve de modifications techniques et de disponibilité des accessoires.  
 Reservadas las modificaciones técnicas y la disponibilidad de los accesorios.

#### Connecting the data output

- Disconnect power cable from balance.
- Using clips (6), hook on data output to openings in rear wall of balance.
- Plug in connector (7) of flat ribbon cable. Flat ribbon cable must sag downwards.
- Push top of data output against the balance until power supply connector (8) is completely plugged in.
- Reconnect balance with data output to power line by plugging power line into socket (9).

#### Connecting data receivers

Data receivers equipped with a current loop interface and those with an RS232C interface can be connected at the same time to an Option 011 Data Output. (In so doing, please note that the configuration of the data output - baud rate, parity, transfer mode - is the same for both interfaces).

- Data receivers with a 20 mA current loop interface are connected to socket (11).
- Data receivers with an RS232C interface are connected to socket (12).
- If switch 4 of operating mode switch (2) is in the OFF position (send on transfer), a hand switch (42500) or a foot switch (46278) must be connected to sockets (10) in order to forward transfer instructions.

#### Removing the data output

- It is absolutely necessary to first remove the power cable from the data output.
- Carefully fold back Data Output.
- Pull out flat ribbon cable connector (7).
- Unhook clips (6) from rear wall of balance.
- Using the cover plate, close ribbon cable connector socket (7) at rear wall of balance.

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