Operating instructions
Bedienungsanleitung
Mode d'emploi

METTLER TOLEDO

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

TABLE OF CONTENTS

	<u>Page</u>
INTRODUCTION	2
A. CL DATA INTERFACE	3
A.1 Connection of the GA40 Printer	3
A.2 Connection of CL instruments	3
A.3 Connection of non-METTLER TOLEDO instruments	3
A.4 Limiting and operating data of the CL data interface	4
B. RS232C DATA INTERFACE	5
B.1 Connection of non-METTLER TOLEDO instruments with RS232C data interfacing	5
B.2 Limiting and operating data of the RS232C data interface	6
C. OPERATION	7
C.1 Data format	7
C.2 Transfer modes	8
C.3 What's wrong if. . . ?	9
C.4 Accessories	9


INTRODUCTION

General remarks

METTLER TOLEDO AE balances can be equipped with an 011 Data Output Option (with a CL and an RS232C data interface).

These data interfaces allow measuring results to be transferred to a data receiver. As a result, the balance can be used as a measuring value generator in a data system.

The alphanumeric printer, the METTLER TOLEDO CL System or a non-METTLER TOLEDO instrument with active 20 mA current loop interfacing (terminal, desk-top calculator) can be connected over the CL Data Interface.


The RS232C Data Interface allows a correspondingly equipped non-METTLER TOLEDO instrument (terminal, (mini-) computer) to be connected.

Further information / technical data

- Consult the particular balance Operating Instructions for information on balance operation.
- Connection and configuration of the data output are explained in the "Installation Instructions for Option 011 Data Output for AE balances". (These instructions are included in the data output kit).
- For a general understanding of the CL Data Interface (instrument and program technical information), we recommend you read the Technical Information Bulletin (TIB), "The METTLER TOLEDO CL Data Interface".


A. CL DATA INTERFACE

With the O11 Data Output Option, both data interfaces (RS232C and CL) are always in operation at the same time. Of course, care should be taken to assure that the configuration of the data output -- baud rate, parity, transfer mode -- is the same for both data interfaces.

A.1 Connection of the GA40 Printer


- Set switch (1) on the printer to "CL".
- Connect cable 42559 (0.5 m)
42561 (1.5 m)
42562 (5 m)
- To release the data transfer (weighing result), a transfer command must be given; this is done by using hand key 42500 or foot pedal (46278).
- The configuration of the data output must correspond to the factory setting (2400 baud/even parity/send on transfer).

A.2 Connection of CL instruments

The components of the METTLER TOLEDO CL Data Communications System can be connected to the CL output of the O11 Option: CL220 HP Adapter, CL241 CL-IEEE Data Interface, CL310 5-Channel Line Selector, CL410 Digital Process Interface, CL530 Alphanumeric Display, CL540 Data Terminal, CL620/CL630 SQC Controller.

- The configuration values must be the same on both the CL instrument to be connected and the O11 Option.
- Connect cable 42931 to the CL port of the data output and to the slave channel socket of the CL instrument.


A.3 Connection of non-METTLER TOLEDO instruments

The O11 Data Output Option has a passive transfer loop. This must be powered by a current source from the data receiver. There is galvanic separation between the transfer loop and the balance electronics.

- The configuration values must be the same in both the instrument to be connected and the O11 Data Output Option (see also below).

Some peripheral instruments (e.g. the GA24/GA44 Printers) already have a transfer button on them. When this button is pressed, the transfer command is sent over the 7-pole connection cable from the peripheral instrument to the O11 Data Output. In such cases, no special transfer mechanism (handkey, foot pedal) needs to be connected to socket (10). (See page 2 of Installation Instructions).

Connector and cable:


The data output has a 7-pole chassis coupling (female) with the following pin allocation (external view):

Pins 1 and 4:	transfer input (Pin 1: plus)
Pins 3 and 5:	transmit loop (Pin 5: plus)
Pins 6 and 7:	blank
Pin 2:	shielding

The connection cable must satisfy the following requirements:


- 4-wire, stranded in pairs
- shielded
- resistance approx. $125 \Omega/\text{km}$
- cable cross section (per conductor) 0.14 mm^2
- capacity approx. 130 nF/km
- admissible cable length up to 1200 Baud: 1000 m; 2400 Baud or more: 500 m

A .4 Limiting and operating data of the CL Data Interface

Transfer principle:

- unidirectional
- bit serial
- asynchronous
- one start bit and one stop bit
- 20 mA current loop (current in idle position)
- 7 bits per character and parity bit (ASCII)
- no handshake

Limiting data:


The transfer loop must be powered by an external current source. To prevent the CL Data Interface from being destroyed by the external current source, the limiting data given in the adjacent diagram absolutely must be observed: the U/I characteristic of the source must be within the shaded area. For LOW, the loop current must be $< 1 \text{ mA}$. The transmit loop in the 011 Data Output Option operates correctly with any polarity of the current source.

To assure distortion-free transfer, the following conditions must also be met:

- Source voltage: 15 V ($+10\%/-0\%$)
- Current (high) $\geq 18 \text{ mA}$ and $\leq 24 \text{ mA}$
- Curve slope $2 \dots 20 \text{ mA}/\mu\text{s}$
- Cable as defined above

B. RS232C DATA INTERFACE


With the 011 Data Output Option, both data interfaces (RS232C and CL) are always in operation at the same time. Of course, care should be taken to assure that the configuration of the data output -- baud rate, parity, transfer mode -- is the same for both data interfaces.

B.1 Connection of non-METTLER TOLEDO instruments with RS232C data interfacing

The 011 Data Output Option is set up as DCE (Data Communication Equipment).

The configuration values on the instrument to be connected and the 011 Data Output Option must be the same (see also below).

Connector and cable:

The data output has a 25-pole, D-subminiature socket (female) with the following pin allocation (external view).

- Pin 1: protective Ground, cable shield connected to housing and protective GND (it is possible to set it to "floating" by removing the jumper inside)
- Pin 3: R x D, transmit (+/-)
- Pin 5: CTS, clear to send (marking)
- Pin 6: DSR, data set ready (marking)
- Pin 7: GND, signal ground
- Pin 8: DCD, carrier detect (marking)
- Pin 20: DTR, data terminal ready (is interpreted)

Depending on whether the instrument to be connected is a DTE (video monitor, printer) or a DCE (personal-, desk-top- or minicomputer), the connection cable must be wired accordingly.

Distortion-free transfer is assured for the following maximum cable lengths: up to 1200 Baud: 5 m

2400 Baud or more: 2 m

In addition, the following conditions apply:

- Use of an external instrument that maintains the limiting and operating data as given in the RS232C EIA standard (see also below).
- The balance and the external instrument must be connected to the same power source.

B.2 Limiting and operating data of the RS232C data interface

- Transfer principle:
- bitserial
 - asynchronous
 - one start bit and one stop bit
 - voltage-controlled
 - 7 bits per character plus parity bit (ASCII)
 - negative voltage in the idle condition

The transmit lines: have the following characteristics:

MARKING / OFF (HIGH): at least -3 Volt

SPACING / ON (LOW): at least +3 Volt

These values are valid with a load resistance of at least $3\text{ k}\Omega$. The sources must be designed so that they are protected against external sources of up to $25\text{ V}/0.5\text{ A}$ (of any polarity).

The maximum permissible current in such cases is 0.5 A . With short-circuit to ground or to another line, maximum permissible current is also 0.5 A . With a load resistance of max. $7\text{ k}\Omega$, the amount of voltage at the balance may not exceed 15 volts.

Receiver lines: The input resistance of the circuit should be between $3\text{ k}\Omega$ and $7\text{ k}\Omega$. Input idle voltage (bias) should not be more than 2 V . The receive circuit must be set up for voltages of at least 25 volts. The maximum capacitance of the receiver circuit and cable may not exceed 2500 pF . Input voltages of -3 V or more negative are interpreted as MARKING / OFF (HIGH). Such voltages of $+3\text{ V}$ and more are interpreted as SPACING / ON (LOW). The absolute values of the voltages do not exceed 25 V . Voltages between -3 V and $+3\text{ V}$ are not defined.

Handshake: If the "Data Terminal Ready" line has a voltage of $\geq +3\text{ V}$, the 011 data output can transmit results. If the value is $\leq -3\text{ V}$, the data output cannot transmit.


If the level changes from positive to negative during a transfer, no more than 2 additional characters will be transferred. If the level changes back to positive, the transfer is continued.

If the line is not wired up, the data output behaves the same as if the external instrument would indicate readiness.

C. OPERATION

C.1 Data format

The results that are transmitted from the balance have a fixed format. The format can be broken down into three blocks. The string is terminated by CARRIAGE RETURN (CR) and LINE FEED (LF):

□: space

The transfer block begins with an identification block (I), which is composed of two characters. This I block is used for more exact identification of the transferred result.

The data block (D) contains actual measuring result. It is separated from the I block by a space. The measuring result is transferred flush right in the D block with sign in front, decimal place and the applicable number of digits. Length of the D block is 9 characters.

The unit block (U) is again separated from the D block by a space. This block contains the unit symbol "g" for grams.

Immediately after the last unit character, the entire transfer block is terminated by CARRIAGE RETURN (CR) and LINE FEED (LF).

Exact behavior of the balance is described in the following paragraphs.

I = Identification block: 2 characters

□ □ result released by the transfer key

S □ stable result

only with the

SD unstable result

Send Continuous

SI invalid result

transfer mode

D = Data block:

9 characters: (incl. decimal point and sign in front)

The result is transferred flush right in the data block.

Leading zeros are suppressed. Plus sign is also suppressed.

The minus sign is always immediately in front of the first digit.

When the DeltaDisplay is switched on, the last two digits are represented as blanks when DeltaDisplay has responded.

U = unit block: 0 to 5 characters:

This block indicates the unit that belongs to the weight value, namely g.

Immediately after the last unit character, the entire transfer block is terminated by CR and LF.

As special cases, the [[SI]] and [[TA]] transfers are immediately terminated by CR and LF.

Special operating mode		
Standby		BREAK
Calibrate		HI
OFF (Display)		HI
Configuration mode		HI
Error		SI SI SI
Lamp test		BREAK
Overload		SI SI SI
Underload		SI SI SI
Tare		SI SI TA

C.2 Transfer modes

The 011 Data Output Option can be set to two different transfer modes by means of a P/C board switch (see Installation Instructions). These modes are:

Send on Transfer: When the TRF button is pressed, the next stable result is transferred. The moment the transfer is made can be seen from the brief blanking out of the display. These manually released results are always identified by [[U]].

Send Continuous: The balance transmits a measuring result after each display cycle, at a maximum of every 0.125 seconds (with baud rates of less than 2400 baud, correspondingly less frequently). There are at least 25 bits "HIGH" between two strings.

Stable measuring results are identified by [[S U]]; unstable results by [[SD]].

If the balance cannot provide a reasonable result (e.g. when there is overload or underload), [[SI]] is transferred.

If the TRF key is pressed while in this transfer mode, the balance transmits the next stable result once with the identification [[U]].

The transfer can be recognized by the brief blanking out of the display.

If tare is pressed while in this transfer mode, the message [[TA]] is transferred automatically after conclusion of the tare operation, followed immediately by CR and LF. [[TA]] is also transferred when the balance is switched on.

C .3 What's wrong if . . . ?

(See also balance Operating Instructions, Chapter, "What's wrong if. . . ?").


... the balance provides no results?

- Incorrect cabling
- The connected instrument is not powering the current loop.
(With CL instruments: see the section on "Connection of non-METTLER TOLEDO instrument").

... the result received are garbled?

- Transfer rate of the external instrument is not the same as the one set for the data output.
- Parity setting of the external instrument is not the same as the one set for the data output.

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