
## OPTION 012 DATA INTERFACE

### INSTALLATION INSTRUCTIONS FOR THE 012 DATA INTERFACE OPTION FOR AE BALANCES

#### General information

A METTLER TOLEDO AE BALANCE can communicate with other instruments (full duplex) by using the bidirectional 012 Data Interface. With this option, two standard data interfaces are available: 20 mA current loop and RS232C data interfacing.


You can obtain information about data transfer with non-METTLER TOLEDO instruments and hardware and software definitions from the data interface description entitled, "The 012 Option bidirectional data interface for AE balances (CL/RS232C)"; Order No. 701323 from METTLER TOLEDO. This description is very strongly recommended if you intend to develop your own programs for balance control.

#### Opening the housing

This is necessary only if one of the configurations that deviates from the factory configuration must be set. (Factory setting corresponds to the METTLER TOLEDO CL standard configuration: 2400 baud/even parity/send on transfer/RS232C cable shielding: not grounded).

- To open the housing, remove the four corner screws (arrows). Turn the housing over and lift off cover.
- Remove the four P/C board fastener screws, carefully lift up P/C board; pull the ribbon cable connector through the housing rear wall.

#### Configuring the data interface

- Baud rate: set according to the instrument to be connected by using jumper (1).
- Operating modes: at the factory, all switches on the operating mode switch (2) are set to the OFF position.


| Switch | ON                                                   | OFF                                                  |
|--------|------------------------------------------------------|------------------------------------------------------|
| 1      | parity disable                                       | parity enable                                        |
| 2      | odd parity<br>[if 1 is OFF]<br>space<br>[if 1 is ON] | even parity<br>[if 1 is OFF]<br>mark<br>[if 1 is ON] |
| 3      | out only *)                                          | bidirectional                                        |
| 4      | send continuous                                      | send on transfer                                     |

Subject to technical changes and to the availability of the accessories supplied with the instruments.  
 Technische Änderungen und Änderungen im Lieferumfang des Zubehörs vorbehalten.  
 Sous réserve de modifications techniques et de disponibilité des accessoires.  
 Reservadas las modificaciones técnicas y la disponibilidad de los accesorios.

\*) Data output: with connection of a printer, the data interface can be set to "out only" using switch 3 (setting: ON → only data output is active).

- Cable shielding (RS232C): the factory setting of jumper (4) is on Pins 2 and 3 (floating). If cable shielding must be placed to the housing/power supply ground, the jumper must be placed on Pins 1 and 2.


#### Adusting the power supply voltage

Note: The power supply voltage must be set on the housing of the data interface.


The voltage selector switch can be moved to the proper position by using a screwdriver:

| Switch position               | Voltage range |
|-------------------------------|---------------|
| 115 V                         | 92...132 V    |
| 220 V                         | 184...265 V   |
| 100 mA slow-blowing microfuse |               |

#### Installing the data interface housing

- Disconnect the power cable from the balance.
- Insert housing with the tabs (6) into the openings in the rear wall of the balance.
- Plug in the ribbon cable connector (7). The ribbon cable must lead in from the bottom.
- Move the top of the housing towards the balance until the power connector (8) is completely plugged in.
- Reconnect the balance with the data interface to the power supply: plug power cable into socket (9).

#### Connecting the instruments

The data interfaces of the 012 Option can be used simultaneously. If this is done, please be aware that the configuration (baud rate, parity, operating mode) is the same for both data interfaces. The input currently unused must be in the BREAK mode to prevent the other interface from being blocked (RS232C: space, CL: no current).

- Instruments with 20 mA current loop data interfacing are connected to socket (11).
- Instruments with RS232C data interfacing are connected to socket (12). (Non-METTLER TOLEDO instrument is DTE).
- If the data transfer is to be triggered by hand ("send on transfer" mode), handkey (42500) or foot pedal (46278) is connected to sockets (10).

#### Removing the data interface

- It is absolutely necessary to first disconnect the power cable from the data interface housing.
- Carefully swing housing straight back and down.
- Disconnect the ribbon cable connector (7).
- Unhook the tabs (6) from the balance rear wall.
- Cover the ribbon cable socket on the balance rear wall with the cover.

## Brief information about bidirectional operation

With the 012 Option, AE balances are equipped with full duplex data interfaces, i.e., the balance can not only transmit weighing results to a data receiver, but also receive, evaluate and carry out certain control instructions.

These instructions are listed below. A detailed description of these instructions can be found in the data interface description, "The 012 Option bidirectional data interface for AE balances (CL/RS232C)", Order No. 701323 from METTLER TOLEDO.

□: space

| INSTRUCTION SET for control of the AE balance |                                                                               |                                            |
|-----------------------------------------------|-------------------------------------------------------------------------------|--------------------------------------------|
|                                               | <u>Send value</u><br>(Send next stable weighing result)                       |                                            |
|                                               | <u>Send Immediate value</u><br>(Immediately send a weighing result)           | <u>Request of measuring value</u>          |
|                                               | <u>Send Immediate value and Repeat</u><br>(Send all weighing results)         |                                            |
|                                               | <u>Tare</u>                                                                   | Tare                                       |
|                                               | <u>Display text</u><br>(Bring "text" to balance display)                      | <u>Control of balance display</u>          |
|                                               | <u>Display text; symbol</u><br>(Bring "text" and "symbol" to balance display) |                                            |
|                                               | <u>Display reset</u><br>(Return to display of weighing result)                |                                            |
|                                               | <u>Remote enable (1)</u><br>(Switch on remote control)                        | <u>Influence of control bar on balance</u> |
|                                               | <u>Remote disable (0)</u><br>(Switch off remote control)                      |                                            |
|                                               | <u>Clear</u><br>(Reset balance and data interface)                            | <u>Reset of all instructions</u>           |

| DATA FORMAT at the data output of the AE balance |                                                                                                                                                                       |
|--------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                  | Transfer triggered by pressing button<br>S Stable result<br>SD Unstable result<br>SI Invalid result<br>Unit, g<br>Data block, weighing result<br>Identification block |

| Balance ERROR MESSAGES | Possible cause                                                                         |
|------------------------|----------------------------------------------------------------------------------------|
| ET Transfer error      | Incorrect baud rate, parity, stop bits, electr. disturbances                           |
| ES Syntax error        | Incorrect instruction set, see above                                                   |
| EL Logistics error     | Instruction sequence to balance too fast or more than seven characters sent to display |

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