**Operating instructions**  
**Bedienungsanleitung**  
**Mode d'emploi**

**METTLER TOLEDO**

**Bidirectional data interface of AE balances**

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
**Option 012 (CL/RS232C)**


**TABLE OF CONTENTS**

	<u>Page</u>
<b>INTRODUCTION</b>	<b>2</b>
<b>A. THE CL DATA INTERFACE</b>	<b>3</b>
A.1 Connection of METTLER TOLEDO CL instruments	3
A.2 Connection of non-METTLER TOLEDO instruments with current loop data interfacing	3
A.3 Limiting and operating data of the CL Data Interface	4
<b>B. THE RS232C DATA INTERFACE</b>	<b>5</b>
B.1 Connection of non-METTLER TOLEDO RS232C instruments	5
B.2 Limiting and operating data of the RS232C Data Interface	6
<b>C. OPERATION</b>	<b>7</b>
C.1 General information	7
C.2 Data format at the data output of AE balances	7
C.3 Instruction set for control of AE balances	8
C.4 Error messages from the balance	11
C.5 Transfer modes	12
<b>D. WHAT'S WRONG IF...?</b>	<b>13</b>
<b>E. ACCESSORIES</b>	<b>14</b>
<b>APPENDIX</b>	<b>43</b>


## INTRODUCTION

### General remarks


METTLER TOLEDO AE balances can be equipped with a bidirectional 012 Data Output Option.

AE balances can transfer weighing results to a data receiver (computer, terminal, etc.) via the two data interfaces (20 mA current loop and RS232C). At the same time, however, they can also receive and carry out instructions (full duplex operation). This makes it possible for AE balances to be integrated into controlled weighing systems.

The following can be connected to the current loop data interface:

- a METTLER TOLEDO CL instrument
- a METTLER TOLEDO printer (**only** with option 011)
- a non-METTLER TOLEDO instrument (computer, terminal, printer).

A correspondingly equipped non-METTLER TOLEDO instrument (terminal, computer, printer, etc.) can be connected to the RS232C data interface.

### Additional information / technical data


Balance operation: consult the particular balance Operating Instructions

Connection and configuration of the data output: all necessary information can be found in the Installation Instructions for the 012 Option (included in the option).

METTLER TOLEDO CL data interface: for a general understanding of the CL Data Interface (technical instrument and programming data), we call your attention to the Technical Information Bulletin (TIB) entitled, "The METTLER TOLEDO CL Data Interface".

## A. THE CL DATA INTERFACE


The data outputs of both interfaces (CL and RS232C) are always in operation at the same time. Of course, it is necessary to note that the communication parameters (baud rate, parity, transfer mode) for both interfaces are the same.

### A.1 Connection of METTLER TOLEDO CL instruments

The components of the METTLER TOLEDO CL Data Communication System can be connected to the CL socket (5-pole) of the 012 Option: The CL220 HP Adapter, the CL241 CL-IEEE Data Interface Adapter, the CL310 5-Channel Line Selector, the CL410 Digital Process Interface, the CL530 Alphanumeric Display, the CL540 Data Terminal, the CL620/CL630 SQC Controllers.

- The communication parameters on the CL instrument to be connected and the 012 Data Interface Option must agree.
- Connect the 5-pole connection cable to the CL socket and to the slave channel socket of the CL instrument.


### A.2 Connection of non-METTLER TOLEDO instruments with current loop data interfacing

The 012 Option data interface has a passive transmit and a passive receive loop. These loops must be powered by the connected instrument.

The current loops and the balance electronics are galvanically separated.

- The communication parameters on the instrument to be connected and on the 012 Option must agree.

Connectors and cables:


The 012 Option has a 5-pole chassis connector (male) with the following pin allocations (external view):

- |              |                        |
|--------------|------------------------|
| Pin 1 and 2: | Receive loop (Pin 2+)  |
| Pin 3:       | Cable shield           |
| Pin 4 and 5: | Transmit loop (Pin 4+) |


The connection cable must meet the following requirements:

- 4-wire, two twisted pairs
- shielded
- resistance approx.  $125 \Omega/\text{km}$
- cable cross section (per conductor)  $0.14 \text{ mm}^2$
- capacitance approx.  $130 \text{ nF/km}$
- maximum admissible cable length up to 1200 baud: 1000 meters  
2400 baud or more: 500 meters

### A .3 Limiting and operating data of the CL Data Interface

- Transfer principle:
- full duplex with no automatic echo
  - bit serial
  - asynchronous
  - one start and one stop bit
  - $20 \text{ mA}$  loop current (current in idle condition)
  - 7 bits per character plus parity bit (ASCII)
  - handshake accomplished

#### Limiting data:

Both transfer loops must be powered by an external current source. To prevent the balance CL data interface from being destroyed by an external current source, it is absolutely necessary that the limiting data given in the adjacent diagram be observed:

The U/I characteristics of the sources must be within the crosshatched area.

- Transfer loop (current level)
  - HIGH: see limiting data
  - LOW: less than  $1 \text{ mA}$
- Receive loop (response threshold)
  - HIGH: more than  $15 \text{ mA}$
  - LOW: less than  $5 \text{ mA}$



Both loops are protected against damage caused by incorrect polarity of the current sources.

To assure that transfer will be free of disturbance, the following conditions must also be met:

- Voltage level of the source:  $15 \text{ V} (+10 \%/-0 \%)$
- Current (HIGH) between  $18\dots24 \text{ mA}$
- Slope  $2\dots20 \text{ mA}/\mu\text{s}$
- Cables defined above

## B. THE RS232C DATA INTERFACE

The data outputs of both data interfaces (CL and RS232C) are always in operation at the same time. Of course, it is necessary to note that the communication parameters for both data interfaces are the same (baud rate, parity, transfer mode).

### B.1 Connection of non-METTLER TOLEDO RS232C instruments

The RS232C data interface of the 012 Option is set up as DCE (data communication equipment). The pin allocations that result from this setup are listed below.

The communication parameters in the instrument to be connected and in the RS232C data interface of the 012 Option must agree.

Connectors and cables:

The standard 25-pole, D-subminiature socket (female) has the following pin allocations (external view).

- Pin 1: Protective ground, cable shielding to housing/power-line ground (can be configured internally to "floating" by using a short-circuit plug).
- Pin 2: TxD, receive line (+/-)
- Pin 3: RxD, transmit line (+/-)
- Pin 4: RTS, Request to send (is evaluated)
- Pin 5: CTS, Clear to send (marking)
- Pin 6: DSR, Data set ready (marking)
- Pin 7: GND, Signal ground
- Pin 8: DCD, Carrier detect (marking)
- Pin 20: DTR, Data terminal ready (is evaluated)

The following maximum cable lengths must be observed to assure disturbance-free transfer:

up to 1200 baud:	5 meters
2400 baud and more:	2 meters

In addition, the following conditions apply:

- Use of an external instrument that operates according to the limiting and operating data as stated in the RS232C EIA Standard.
- The balance and the external instrument must be connected to the same power source.

Note: The RS232C connection cable from METTLER TOLEDO (89164 or 89165) only satisfies a minimum configuration. Only the following lines are complete: Pins 1...7 and Pin 20. The "carrier detect" line, Pin 8, for example is not connected.

## B.2 Limiting and operating data of the RS232C data interface

- Transfer principle:
- bit serial
  - asynchronous
  - one start and one stop bit
  - voltage-controlled data interface
  - 7 bits per character plus parity bit (ASCII)
  - negative voltage in idle condition

The transmit lines have the following characteristics:

MARKING / OFF (HIGH): at least -3 Volt

SPACING / ON (LOW): at least +3 Volt

These values are valid at a load resistance of at least  $3\text{ k}\Omega$ . The sources are set up so that they are protected against external current sources up to 25 V / 0.5 A (of any polarity).

The maximum permissible current in such cases is 0.5 A. With short-circuit to ground or to another line, the maximum permissible current is also 0.5 A. With a load resistance of a maximum of  $7\text{ k}\Omega$ , the voltage level at the output of the balance does not exceed 15 volts.

Receiver lines: The input resistance of the circuit lies between  $3\text{ k}\Omega$  and  $7\text{ k}\Omega$ . The idle input voltage (bias) is not more than 2 V. The receive circuit is set up for voltages of at least 25 volts. The maximum capacitance of the receiver circuit and cable does not exceed 2500 pF. Input voltages of -3 V and less are interpreted as MARKING / OFF (HIGH). Those of +3 V and more are interpreted as SPACING / ON (LOW). The absolute values of the voltages must not exceed 25 V. Voltages between -3 V and +3 V are not defined.

Handshake: If the “Data Terminal Ready” line indicates a voltage of  $\geq +3\text{ V}$ , the 012 data output transmits results. If the value is  $\leq -3\text{ V}$ , the data output will not transmit. If the level changes from positive to negative during a transmission, a maximum of 2 characters will be sent. If the level changes back to positive, the transmission is continued.

If the line is not used, the data output operates in the same way as when the external instrument registers readiness.

## **C. OPERATION**

### C.1 General information


The 012 Option allows AE balances to be equipped with maximum capability data interfaces (current loop and RS232C), i.e., measuring results are not only transferred to a data receiver, but certain control instructions can also be sent to the balance in the command direction. For the most part, these instructions are composed of control information which can be used for upper and lower case letters. Each instruction must always be terminated by the character sequence CARRIAGE RETURN (CR) and LINE FEED (LF). In addition, all information transferred from the balance to the data receiver is also terminated with the same character sequence, CR LF.

Every interruption of the data interface (BREAK) puts the balance back into the normal weighing mode. Instructions that are not carried out are cleared.

If the balance should be addressed with a quick series of instructions, it can be briefly overloaded. To prevent this, the data interfaces of the 012 Option are provided with controlled operation (handshake).

### C.2 Data format at the data output of AE balances

Every measuring results is simultaneously passed on to both data outputs (CL and RS232C) in a uniform format. The transfer data string can be divided into three blocks. Each data string is terminated by CARRIAGE RETURN and LINE FEED.

□: space

The transfer data string begins with an identification block (I), which is composed of two characters. This I block is used for exact identification of the transferred result.

The data block (D) contains the actual measuring result. It is separated from the I block by a space. The measuring result is transferred in the D block with sign in front, decimal point and the particular applicable number of places. The length of the D block is nine characters.

The unit block (U) is also separated from the D block by a space. This block contains the unit symbol «g» for grams.

The entire transfer data string is terminated by CARRIAGE RETURN (CR) and LINE FEED (LF) immediately after the last unit characters.

A precise description of the behavior of the balance is given in the following sections.

I = Identification block: 2 characters

□ □ result triggered by pressing the transfer key

S □ stable result

results triggered

SD unstable result

by

SI invalid result

instructions

D = Data block:

9 characters: (incl. decimal point and sign in front)

The result is transferred right justified. Leading zeros are suppressed. A plus sign (for positive results) is also suppressed. The minus sign always appears immediately in front of the first digit.

When DeltaDisplay is switched on, the last two places are blanked out when DeltaDisplay responds.

U = unit block:

0 to 5 characters:


This block indicates the unit that corresponds to the weight value, namely g.

The entire transfer block is terminated by CR and LF immediately after the last unit character.

Transfer of [[SI]] and [[TA]] occur as special cases, terminated immediately by CR and LF.

Special operating conditions:	Standby	BREAK
	Calibrate	HI
	OFF (display)	HI
	Configuration mode	HI
	Error	SI SI SI
	Lamp test	BREAK
	Overload	SI SI SI
	Underload	SI SI SI
	Tare	SI SI TA

### C.3 Instructions set for control of AE balances

Send value

By using this instruction, the controlling instrument can request single measuring values via the data interface. After receipt of «S», the balance transmits the next stable weight result. The moment of transfer to the data receiver can be recognized by the brief blinking of the balance display.

A data transfer triggered in this manner is identified in the data string by «S\_», and differs from transfers triggered by pressing a handkey or foot pedal, whereby two blanks «\_ \_» appear in the identification block.

If the balance cannot make a meaningful weighing result available, it transmits «SI» (I for invalid) at this point. This identification occurs, for example, when the weighing pan is at overload.

Several results can be requested by transmitting the instruction «S» several times. In this situation, it is to be noted that no more than one result is transmitted per display cycle and that instructions are not stored. For example, if an instruction has not been carried out, it is written over by the next one to arrive. It is also cleared if a BREAK occurs in the data interface.

S	I	C <sub>R</sub>	L <sub>F</sub>
---	---	----------------	----------------

Send Immediate value

For dynamic measuring, it is possible to make a request to the balance for weighing results that have not been released by the stability. A measuring result is transferred at the conclusion of the current display cycle by the instruction «SI». To differentiate dynamic from stable measuring results, the identification block can also contain the character sequence «SD» (D for dynamic) instead of «S<sub>U</sub>». The display does not blink briefly when measuring value transfers are triggered by this instruction.

Because the «SI» instruction request current results, it is possible that an «SI» will be transmitted instead of the weighing result under special operating conditions (see also Chapter C.2, Special operating conditions).

While weighing-in, the last two decimal places are replaced by blanks (DeltaDisplay) for dynamic results.

S	I	R	C <sub>R</sub>	L <sub>F</sub>
---	---	---	----------------	----------------

Send Immediate value and Repeat

The instruction «SIR» is an order for the balance to transmit the particular measuring result after each display cycle. In its effect, the instruction has the same meaning as a quick succession of «SI» instructions (see above).

The data rate can get very high when working with this instruction. The rate depends on the length of the display cycle and whether or not DeltaDisplay is switched on. If transfer of all results is wished via the data interfaces, the baud rate must be set to at least 2400 baud.

The «SIR» instruction can be written over by transmitting an «S» or «SI» instruction, or cleared by transmitting «C» (Clear). Interruption of the data interface (BREAK) also cancels out «SIR».

T	C <sub>R</sub>	L <sub>F</sub>
---	----------------	----------------

Tare

A tare procedure can be triggered by the instruction «T» over the data interface instead of pressing the balance single control bar.

If the balance display shows "OFF", the instruction is also accepted. In this special case, the tare instruction «T» causes the balance to switch over to the normal weighing mode.

If an attempt is made to trigger tare while the balance is in overload/underload range, the error message «EL» is sent (Logistic Error, see Chapter on error messages).

D	█	text	C <sub>R</sub>	L <sub>F</sub>
---	---	------	----------------	----------------


Display text

Generally speaking, the balance display is used to display weighing results. However, by using the instruction «D<sub>U</sub>text» messages can also be displayed, e.g. to guide the person operating the balance or to display values that have been requested and converted. By using the «D<sub>U</sub>text» instruction, a "text" character sequence can be written in the balance display, right justified. As text characters, all printable characters in the ISO 646 Code Table are permitted except semicolon. The only point to be considered is the limited illustration capabilities of the 7-segment display (Table in Appendix).

A maximum of seven characters can be illustrated, corresponding to the number of decimal places the balance can display. An exception to this is the decimal point which can also appear after each character. If these rules are not maintained, the balance reacts with the error message, «EL».

Missing text, i.e., the «D» instruction, is interpreted as a blank text and leads to a balance display that is blanked out. Contrary to this, «D» leads to return to the normal display mode (see further below).

The balance operates again in a completely normal fashion after receipt of a display instruction, i.e., measuring results can be requested by using request instructions, tare can be triggered, etc. without these instructions being affected by the missing display of the weighing result.

In addition to the displayed text, or converted weighing result, a symbol can also be added, separated from the text by a semicolon. This character controls the display location on the left of the display.

The following symbols can be used for sym:

	(space)	→	(blank)
+	(plus)	→	(blank)
-	(minus)	→	-
o	(letter)	→	o often used for converted weighing results

Furthermore, a unit character can be added, again separated by a semicolon. However, this character is ignored by AE balances.

With this instruction, the display is again released for a weighing result. All displays triggered by display instructions are cleared.

An interruption of the data interface (BREAK) has the same effect as the above.

The balance single control bar can be disabled with the «R1» instruction. The balance can then accept tare instructions only by way of the data interface. This remote control is switched off again by the instruction «R0».

This instruction has the same effect as switching the balance on/off. As with the «T» instruction, a switch out of the "OFF" mode can be accomplished by using the «C» instruction.

#### C.4 Error messages from the balance

In certain situations, the balance is not in a condition to carry out the instructions received (e.g. an instruction that cannot be encoded, violation of certain limitations). In such cases, the instruction received is rejected and not carried out. To let the instruction transmitter know about this, the balance transmits an error message via the data interface.

An error message is composed of two letters:

«ET» Transmission Error

«ES» Syntax Error

«EL» Logistic Error

#### Transmission errors, «ET»

The «ET» error message is transmitted when the character sequence received from the balance is not correct. Character sequences are declared not correct when the parity of one or more characters is not correct or the formatting of the individual characters is not correct (missing stop bit or similar problems). Such errors occur in operation only very sporadically if the baud rate and parity settings are correct. Their cause is usually the result of electrical disturbances in the transfer lines; malfunctions of the electronics are also possible and should be taken into consideration.

#### Syntax errors, «ES»

Correctly received instructions are declared syntax errors when they do not hold exactly to the instruction form defined in the preceding chapters. For example, the request instruction «S1R» is acknowledged by the balance with a syntax error because the correct form is «SIR».

Syntax errors typically occur only in the design phase of the system during test runs. In this phase, the programs of the controlling computer are developed and tested. In a system that has been tested, program routines should be worked out so that no syntax errors occur during operation.

#### Logistic errors, «EL»

Even though the balance has received an instruction that is correct in its transfer and syntax, it still does not mean that it can be carried out. This would be the case, for example, if more than the permitted 7 characters (plus any decimal points) are transmitted to the balance display. In such a case, the instruction is ignored and an error message is transferred to the instruction transmitter.

If the operating speeds of the various balance elements are not regarded, there is a special problem. With instruction sequences that are too fast, the balance is not capable of working out the sequence. For this reason, it will suppress a number of instructions and acknowledge them with the error message «EL».

## C.5 Transfer modes

The O12 Option data output can be set to two different transfer modes by means of a switch on the P/C board (see Installation Instructions). These are:

Send on Transfer: After pressing the TRF button, the next stable result is transferred. The moment of transfer can be recognized by the brief blinking of the balance display. These manually triggered results are always identified by [[UU]].

Send Continuous: The balance transmits a measuring result after each display cycle, however a maximum of every 0.125 seconds (with baud rates of less than 2400 baud, correspondingly slower).

Stable results are identified by [[S U]]; unstable results by [[SD]].

If the balance cannot provide a reasonable result (e.g. when there is overload or underload), [[SI]] is transferred.

If the TRF button is pressed while in this transfer mode, the balance transmits the next stable result once with the identification [[UU]].

Transfer can be recognized by the brief blinking of the display.

If tare is pressed while in this transfer mode, the message [[TA]] followed immediately by CR and LF is automatically transferred after conclusion of the tare routine. [[TA]] is also transmitted when the balance is switched on.

## **D. WHAT'S WRONG IF . . . ?**

(see also the balance Operating Instructions, the chapter entitled, "What's wrong if. . . ?").

... the balance transmits no results?

- The connected instrument is not powering the current loop or current is not sufficient.
- The data input was put out of order by setting "OUT ONLY" (see also the Option 012 Installation Instructions).
- The handshake lines do not have the proper level (e.g. line 8 of the METTLER TOLEDO RS232C cable is not connected).
- The baud rates of the instruments do not agree.
- Transmit and receive loops have been crossed (see TIP, "The METTLER TOLEDO CL Data Interface").

... the balance transmits only garbled results?

- Incorrect baud rate has been set. (see the 012 Option Installation Instructions)
- Incorrect parity setting.

... the balance is not accepting instructions?

- The 012 Option has been configured to "OUT ONLY" (ST5).
- Baud rate, parity incorrectly set (see Option 012 Installation Instructions).
- Ribbon cable between the balance and the 012 Option not plugged in. (see Option 012 Installation Instructions).

... the balance responds to each instruction  
«EL»

- The balance is in the "OFF" mode. First transmit «T» or «O» or with trigger tare by hand. (see balance Operating Instructions).

... the balance always transmits «SI»?

- The balance is outside the permissible weighing range (overload or underload).

... the balance does not transmit a measuring value when the transfer button is pressed?

- The balance pan has not stabilized (decimal point all the way to the left of the displays is still on).
- Transfer button defective or incorrectly connected.
- Transmit and receive loops have been crossed. (see TIB, "The METTLER TOLEDO CL Data Interface").
- The connected instrument is not powering the current loop.

## E. ACCESSORIES

		<u>Order No.</u>
	CL connection cable 5-pole/5-pole	0.5 m      42554 1 m      42560 2 m      42555 5 m      42556 15 m      42557 30 m      42558
	RS232C connection cable, 2 meters coupling - plug (female- male)	89164
	RS232C connection cable, 2 meters plug- plug (male - male)	89165
	Handkey	42500
	Foot pedal	46278

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