# **OPTION 013 DATA INTERFACE**

# INSTALLATION INSTRUCTIONS FOR THE 013 DATA INTERFACE OPTION FOR AE BALANCES

## General information

A METTLER TOLEDO AE balance can be integrated into an IEEE488 bus system by using the bidirectional 013 Data Interface Option.

The data interface description entitled, "The 013 Bidirectional Data Interface for AE balances (IEEE/HB-IB)" Order No. 701326 available from METTLER TOLEDO can provide you with information about data transfer with non-METTLER TOLEDO instruments as well as about hardware and software definitions. The description is very strongly recommended if you intend to develop programs for balance control on your own.

# Opening the housing

This is necessary only if a configuration that is different from the factory setting must be set.

Factory setting: data exchange with Service Request (SRQ)/send on transfer/ hardware reset if malfunction occurs/cable shielding grounded.

- To open the housing, remove the screws from the four corners (arrows). Turn the housing over and lift off cover.

#### Configuring the data interface

- Operating mode switches (1): do not change the position of switch 1.
- Switches 2 and 4 are set to OFF at the factory: switch 3 is set to ON.

Switch	ON	OFF
1	If program is on single chip micro- processor	If program is on external EPROM
2	Data exchange with no SRQ (Service Request)*	Data exchange with SRQ (Service Request)
3	Hardware reset if malfunction occurs	No hardware reset if malfunction occurs
4	Send Continuous	Send on Transfer

- \*) Data transfer controlled by evaluation of the balance status byte (see Table below).
- Cable shielding (Pin 12): coding plug (2) is placed on Pins 2 and 3 at the factory, i.e., the cable shielding is grounded. It can also be placed on Pins 1 and 2 (floating).

© Mettler-Toledo GmbH 2000

Subject to technical changes and to the availability of the accessories supplied with the instruments.

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English

9

# Changing the power-line voltage setting

Note: The power-line voltage must be set on the housing of the data interface.

The voltage selector switch can be set to the proper range by using a screwdriver:

Switch position	Voltage range
115 V 220 V	92132 V 184265 V
Microfuse: 100 mA slow-blowing	

# Installing the data interface housing

- Disconnect the power cable from the balance.
- Insert the tabs (6) on the housing into the openings on the rear of the balance.
- Plug in the ribbon cable connector (7). The loop of the ribbon cable must lead in <u>from the bottom</u>.
- Swing housing up against the balance until the power connector (8) is plugged in completely.
- Reconnect source by plugging the power cable back into connector (9).

# Connecting to the IEEE488 bus

The 013 Data Interface Option conforms to the IEEE488 standard. The balance can be addressed from a control instrument (Controller) as a "Talker" and a "Listener".

 The instrument address or the "Talker Only" function is defined with bit switches (11) (binary 0...30 and TON).

**Device Address 15** 

- TON (Talker Only): the balance can only transmit (Talker). Addressing cannot be carried out.
- Connect the IEEE488 cable to socket (12).
- If data transfer is to be triggered by hand ("Send on Transfer" mode), connect the handkey (42500) or foot pedal (46278) to sockets (10).

# Removing the data interface

- It is <u>absolutely</u> necessary to first disconnect the power cable from the data interface housing.
- Carefully swing housing straight back and down.
- Disconnect the ribbon cable connector (7).
- Unhook the tabs (6) from the balance rear wall.
- Cover the ribbon cable socket on the balance rear wall with the cover.

## Brief information about bidirectional operation

AE balances are equipped with a full duplex data interfaces with the 013 Option, i.e., the balance can not only transmit weighing results to a data receiver, but also receive, evaluate and carry out certain control instructions.

These instructions are listed below. A detailed explanation of these instructions can be found in the data interface description entitled, "The 013 Bidirectional Data Output Option for AE balances (IEEE488/HP-IB)", Order Number 701326 available from METTLER TOLEDO.

ப: space

INSTRUCTION SET for control of the AE balance		
S C <sub>R</sub> L <sub>F</sub>	<u>S</u> end value (Transmit next stable weighing result)	
S I C <sub>R</sub> L <sub>F</sub>	<u>S</u> end <u>I</u> mmediate value (Immediately send a weighing result)	Request of measuring
S I R C <sub>R</sub> L <sub>F</sub>	Send Immediate value and Repeat (Send all weighing results)	value
T C <sub>R</sub> L <sub>F</sub>	<u>T</u> are	Tare
D _ text $C_R L_F$	<u>D</u> isplay text (Bring "text" to balance display)	
D = text ; sym, $C_R L_F$	<u>D</u> isplay text; symbol (Bring "text" and "symbol" to balance display)	Conrol of balance display
D C <sub>R</sub> L <sub>F</sub>	<u>D</u> isplay reset (Return to display of weighing result)	
R 1 C <sub>R</sub> L <sub>F</sub>	<u>R</u> emot enable (1) (Switch on remote control)	Influence of control bar on balance
R 0 C <sub>R</sub> L <sub>F</sub>	<u>R</u> emote disable (0) (Switch off remote control)	
C C <sub>R</sub> L <sub>F</sub>	<u>C</u> lear (Reset balance and data interface)	Reset of all instructions

STATUS BYTE of AE balances			
X + + X X X X	Ø – 0 or 1 · X –	don't care (not defined)	
Ready for command (RFC) Message available (MAV) SERVICE REQUEST (SR)		Balance is ready for instruction Balance holding result ready "Service Request" bit (IEEE488)	

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