**Operating instructions
Bedienungsanleitung
Mode d'emploi**

METTLER TOLEDO

**Bidirectional data interface of AE balances
013 option (IEEE488/HP-IB)**

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**Interface bidirectionnel des balances AE
Option 013 (IEEE488/HP-IB)**


TABLE OF CONTENTSPage**INTRODUCTION****2****A. IEEE488 DATA INTERFACE****3**

- A.1 General description of the IEEE488 Bus System
- A.2 The IEEE488 Data Interface of AE balances
- A.3 Connectors and cables

3

4

6

B. OPERATION**7**

- B.1 General information
- B.2 Data format at the data output of AE balances
- B.3 Commands used for the control of AE balances
- B.4 Error messages from the balance
- B.5 Transfer modes

7

7


8

11

12

C. WHAT IF...?**13****D. ACCESSORIES****13****E. DATA TRANSFER FROM 013 OPTION TO PERSONAL COMPUTER (IEEE488/HP-IB)****14****APPENDIX****43**


INTRODUCTION

General review

METTLER TOLEDO AE balances can be equipped with the bidirectional Option 013 Data Interface (IEEE488).

This interface can be used by AE balances to transfer weighing results to one or more data receiver (computers, terminals, printers, etc.). It can also receive commands and carry them out. This makes it possible to integrate AE balances in an IEEE488-controlled bus system.

The IEEE488 Bus is used for the exchange of data between two or more instruments which are all connected to the same data bus. In so doing, not all the instruments are forced to participate in the data exchange.

Additional information / technical specifications

Balance operation:

Consult the appropriate Operating Instructions

Connection and configuration
of the data interface:

All the necessary information is packed with the Installation Instructions
for Option 013.

IEEE488 Data Interface:

Various authors
IEEE488 Standard Digital Interface for Programmable Instrumentation

A. IEEE488 DATA INTERFACE


Standard: IEEE Standard Digital Interface for Programmable Instrumentation IEEE488

The Option 013 Interface for AE balances corresponds to the above mentioned standard. Other designation for this standard could be GPIB (General Purpose Instrument Bus) or HP-IB (Hewlett-Packard Interface Bus). Only the instrument socket does not conform to the standard for the connection to an IEC bus system. However, as far as the electronics and the functions are concerned, these interfaces are identical.

A.1 General description of the IEEE488 Bus System

The control instrument (CONTROLLER) defines a data source (TALKER) as well as one or several data receivers (LISTENERS) by means of addressing them.

TALK ONLY is a special function which requires no addressing of instruments.

For the above mentioned circuits, negative logic applies. This is to say:

LOW \cong voltage $\leq 0,8 \text{ V} \cong$ logic 1

HIGH \cong voltage $\geq 2,0 \text{ V} \cong$ logic 0

A.2 IEEE488 Data Interface of AE balances

With Option 013, any AE balance is equipped with a interface that conforms to the IEEE488 standard. The balance supports the following interface functions:

| | |
|----------------------|---|
| TALKER (T) | <p>As a matter of principle, a data transfer from the sending storage onto the bus takes place after every talker address call.</p> <p>If the sending storage is empty, or if no concluding character exists, no data string can be transmitted. The bus remains blocked. Before every talker address call, it is recommended to use a SERIAL POLL to determine whether a concluded data string exists (→ STATUS BYTE). The end of the string is marked by the IEEE488 bus circuit EOI.</p> |
| LISTENER (L) | <p>The instructions are accepted by the balance after the listener address has be recalled. The end of the instruction block (concluding character) is recognized from the program. The EOI circuit is not evaluated.</p> |
| SERVICE REQUEST (SR) | <p>If the complete data string and the concluding character are present in the IEEE488 sending storage, the SRQ circuit of the bus becomes active.</p> |
| SERIAL POLL (SP) | <p>The control instrument (CONTROLLER) can use the SERIAL POLL to request information concerning the status of the instruments connected to the bus. The instruments thus addressed send a STATUS BYTE.</p> <p>With an active SRQ circuit, a SERIAL POLL request is generally initiated.</p> |

| STATUS BYTE of AE balances | | | | | | | |
|----------------------------|---|---|---|--|----------------------------------|---|---|
| T | R | S | 4 | 5 | 6 | 7 | 0 |
| X | φ | φ | φ | X | X | X | X |
| | | | | φ = ∅ or 1; X = don't care (not defined) | | | |
| | | | | If bit 4 = 1 | Balance is ready for instruction | | |
| | | | | If bit 5 = 1 | Balance holding result ready | | |
| | | | | If bit 6 = 1 | "Service Request" bit (IEEE488) | | |

| | |
|--------------------------------------|--|
| SOURCE & ACCEPTOR HANDSHAKE (SH, AH) | <p>The mutual interchange of data or instructions only takes place if the receiving instrument is ready for it. As a result of this controlled interface operation, the loss of data is avoided.</p> |
|--------------------------------------|--|

If more than two instruments are operated by the bus, it is necessary to allocate a "Device Address" to each instrument.

In AE balances equipped with an Option 013, the device address is selected with a bit switch (see Installation Instructions for the Option 013).

The table in the Appendix makes it easier to convert the decimal address (0...30) to the corresponding switch position and to program the TALKER and LISTENER addresses (ASCII characters).


Special operating modes

Various operating modes can be selected on the housing of the Option 013 Interface. The switches used for this purpose are indicated in the Installation Instructions.

| | |
|---|---|
| TALKER ONLY (TON) | The balance is the only data source. It can only send. There is <u>no</u> addressing of the instruments. |
| Data operation <u>without</u> SRQ | If the SERVICE REQUEST circuit is not being evaluated by the controller, the data traffic can also be built up without this interface function. For this purpose, the STATUS BYTE sent out by the balance must be evaluated (see SERVICE REQUEST and SERIAL POLL). |
| Hardware Reset in case of disruption (Watchdog) | Periodically, the microprocessor program runs through a loop in which a counter is reset. If this program gets hung up somewhere - for example, because of an electrical disruption in the power supply line or in the data circuit - this counter is no longer reset. Consequently, the microprocessor indicates a "Hardware Reset" when the counter has reached its end value (Watchdog). |
| Internal ROM | The position of slide switch 1 at the operating mode switch should <u>not</u> be changed by the user (correct setting: "OFF"). |
| Grounding concept | The shield (SHIELD, pin 12) of the IEEE488 cable is connected at the balance interface with the protective ground. If this shield is already grounded at some other point, it can be set to "floating" (see Installation Instructions for Option 013). This makes it possible to avoid ground loops. |


A.3 Connector and cable

The connector at the housing of the Option 013 Interface conforms to the IEEE488 standard. Therefore, it is not possible to connect a cable which is conform to DIN IEC 66.22 without using an adapter. Such an adapter can be prepared with the help of the connector layout indicated below.

Pin

| | | |
|----------|-----------|-----------------------------------|
| 1 ... 4 | DIO 1...4 | Data circuits (Data Input/Output) |
| 5 | EOI | End or Identify |
| 6 | DAV | Data Valid |
| 7 | NRFD | Not Ready for Data |
| 8 | NDAC | Not Data Accepted |
| 9 | IFC | Interface clear |
| 10 | SRQ | Service Request |
| 11 | ATN | Addresses or data on DIO circuits |
| 12 | SHIELD | Cable shielding |
| 13 ...16 | DIO5...8 | Data circuits (Data Input/Output) |
| 17 | REN | Remote Enable |
| 18 ...24 | GND | Signal Grounds |

Up to 15 instruments can be connected to the IEEE488 bus. At the same time, cable length is to be held to a minimum (starshaped wiring). Make sure that the total length of the entire cable system is no larger than the number of instruments connected to the bus x 2 m, but at the most 20 m.

- Carry out the configuration of the balance interface according to the Installation Instructions.
- Using the bus cable, connect the instruments with each other.


B. OPERATION

B.1 General information

The data interface of Option 013 makes possible the bidirectional data exchange between the AE balances and the IEEE488 Bus System. This means that not only measuring results are transferred from the balance in the command direction. These commands or instructions generally consist of a control information in the form of an ASCII character string for which both upper and lower case letters can be used. Every command must always be concluded with the characters CARRIAGE RETURN (CR) and LINE FEED (LF). Every information transmitted from the balance to the IEEE488 Bus is concluded in the same manner with the character sequence, CR and LF.

B.2 Data format at the data output of AE balances

Every measuring results is transferred in uniform formatting as an ASCII string to the data bus. The transfer string can be divided into three blocks. It is always concluded with CARRIAGE RETURN and LINE FEED.

□: space

The transfer string begins with an identification block (I) which consists of two characters. This block is used to more specifically identify the transferred result.

Data block (D) contains the actual measuring result. It is separated from the I-block by a space. In the D block, the measuring result, including the sign, the decimal point and the applicable number of decades is transferred right justified. The D-block is 9 characters long.

Unit block (U) is again separated from the D-block by a space. This block contains the unit symbol «g» for gram.

Immediately after the last unit symbol, the entire transfer string is concluded with CARRIAGE RETURN (CR) and LINE FEED (LF)

The exact behavior of the balance is described in the following sections.

I = Identification block: 2 characters

□ □ Result triggered by the transfer button

S □ Stable result

SD Unstable result

SI Invalid result

Results

triggered by

instructions


D = Data block: 9 characters: (including decimal point and sign).
 The result is transferred in the data block right justified. Zeros that come before the result are suppressed. The same holds true for the plus sign. The minus sign is always placed immediately before the first number.
 If the DeltaDisplay is switched on, the last two numbers are indicated as blanks, provided the DeltaDisplay has reacted.

U = unit block: 0...5 characters:
 This block indicates the unit of weight that belongs to the weight value, i.e., grams.
 Immediately after the last unit character has been sent, the entire transfer block is concluded with CR and LF.

The [[SI]] and [[TA]] transfers, which are concluded immediately with the CR and LF, are special cases.

| | | |
|-------------------------------|--------------------|----------|
| Special operating conditions: | Standby | BREAK |
| | Calibrate | HI |
| | OFF (display) | HI |
| | Configuration mode | HI |
| | Error | SI SI SI |
| | Lamptest | BREAK |
| | Overload | SI SI SI |
| | Underload | SI SI SI |
| | Taring | SI SI TA |

B.3 Commands used for the control of AE balances


Send value

With this command, the control instrument - acting through the data interface - can request information about individual measurement values. After receiving an «S», the balance sends the next weighing result after the weighing pan has become stable. The moment at which the transmission takes place can be recognized by the brief blinking-out of the balance display.

A data transfer triggered in this manner is identified in the data string with «S_L», in contrast to a data transfer triggered with a hand or foot control. When either of these two is actuated, two blanks «_{LL}» appear in the identification block.

If the balance cannot provide a meaningful weighing result, it transmits at this place an «SI» (I for invalid). This identification occurs if, for example, the weighing pan stabilizes with an overload on the pan.

If the command «S» is sent several times, several results can be requested. In this situation, it should be noted that at most one result can be sent for each display cycle, and that commands cannot be stored. Thus, if a command has not been carried out, it will be overwritten by the next arriving command.


Send ImmEDIATE value

In dynamic measuring processes, the balance can also be asked for weighing results which have not been released by the stability detector. Upon the command «SI», a measurement value is transferred immediately at the end of the current display cycle. To differentiate between dynamic and stable measuring results, the identification block may contain the character sequence «SD» (D for dynamic) instead of the character sequence «SL». When the transfer of measurement values is triggered by this command, the display is not blanked out.

Since the command «SI» asks for current results, it is possible that under certain special operating conditions an «SI» will be sent instead of the weighing result (see also Section B.2, Special operating conditions).

During the weighing-in process, the last two digits of the dynamic results are blanked out (DeltaDisplay).


Send ImmEDIATE value and Repeat

The command «SIR» instructs the balance to transfer the pertinent measuring result after every display cycle. The effect of this command is similar to a rapid sequence of «SI» commands (see above).

When working with this command, the rate of data can become very large. It depends on the duration of the display cycle and on whether the DeltaDisplay is switched on or not.

The command «SIR» can be written over by sending an «S» or an «SI», or it can be cancelled with the character «C» (Clear).


Tare

Instead of taring the balance with the single control bar, this can also be done via the data interface with the command «T».

If the balance is in the "OFF" position (display), this command is accepted. In this special case, the tare command «T» causes the balance to switch over to the standard weighing mode.

If an attempt is made to tare the balance in the overload/underload range, the balance reacts with the error message «EL» (Logistic Error, see section on Error Messages).

Display text

The display field of the balance is used primarily to display weighing results. However, by using the command «D└text» it is also possible to make messages visible, e.g., to instruct the balance operator, or to display values that have been determined by request and then converted.

With the command «D└text», a "text" sequence can be written right justified, on the balance display. Admissible for these texts are any desired printable characters of the ISO 646 code table, except for the semicolon (;). The only thing that needs to be taken into account is the limited ability of the 7-segment display to represent these characters (see Table in the Appendix).

A maximum of seven characters can be represented, depending on the number of decades in the balance display. The only exception is the decimal point which can appear as an addition to every single character. If these rules are not observed, the balance reacts with the error message, «EL».

Missing texts, i.e., the command «D□», are interpreted as empty texts and, therefore, cause the balance display to go blank. Contrary to this, the command «D» causes the balance to return to its normal display mode (see further information below).

Even after receiving a display command, the balance continues to operate normally. Thus, it is possible to request measuring results, to tare, etc., without these commands being affected by the missing display.

| | | | | | | |
|---|---|------|---|-----|----------------|----------------|
| D | □ | text | : | sym | C _R | L _F |
|---|---|------|---|-----|----------------|----------------|

Display text; symbol [: unit]

In addition to the displayed text, or to the converted weighing result, it is also possible to add a symbol which is separated by a semicolon. This character controls the location all the way to the left of the display.

The following symbols can be used for sym:

- (Space) → (dark)
- + (Plus) → (dark)
- (Minus) → -
- o (letter) → □ is used frequently for converted weighing results

In addition, a unit symbol - again separated by a semicolon - may be added. However, this symbol is ignored by the AE balance.

| | | |
|---|----------------|----------------|
| D | C _R | L _F |
|---|----------------|----------------|

Display Reset

This command serves to release the display for weighing results. All readings caused by display commands are erased.

| | | | |
|---|---|----------------|----------------|
| R | 1 | C _R | L _F |
|---|---|----------------|----------------|

Remote

The single control bar of the balance can be disabled with the command «R1». In that case, the balance can accept tare commands only via the data interface. This remote control operation is again switched off with the command «R0».

| | | |
|---|----------------|----------------|
| C | C _R | L _F |
|---|----------------|----------------|

Clear

This command has the same effect as the switching-off and then on again of the balance. Just as with the command «T», the command «C» can also be used to switch the balance on by starting from the "OFF" mode.

B.4 Error messages

In certain situations, the balance can not carry out a command it has received (e.g. a command that cannot be deciphered, an infraction against existing limitations). In these cases, the command thus received is rejected and is not carried out. To enable the sender of the command to recognize this fact, the balance sends an error message via the data interface.

Such an error message consists of two characters:

«ES» Syntax Error

«EL» Logistic Error

Syntax error «ES»

Syntactic errors are considered those commands which have been correctly received, but do not adhere exactly to the command forms defined in the preceding sections. For example, the request command «SIR» will be considered a syntactic error by the balance, because the correct command form is «SIR».

Typically, syntactic errors occur only in the build-up phase of the system during test runs. During this phase, the programs of the controlling calculator are being developed and tested. The program sequences of the fully tested system should no longer contain any syntactic errors.

Logistic error «EL»

Even though a command may have been correctly received from a transfer and syntactic point of view, the balance may still not be able to carry it out. For example, more than the seven permitted characters (plus decimal points, if any) may have been sent to the balance display. In this case, the command is ignored and the error message is forwarded to the sender of the command.

A special problem exists if the working speed of the various balance facilities has not been respected. If the commands follow each other too fast, the balance can not work off the backlog. It will therefore suppress some command and come up with the error message «EL».

B.5 Transfer modes

The Option 013 Data Output can be set by means of the bit switch (see Installation Instructions) to two different transfer modes. These are:

Send on transfer: After the TRF button is pressed, the next stable result is transferred. The moment at which this transfer takes place can be recognized by the brief blanking-out of the display. These manually triggered results are always identified with `[[L L]]`.

Send continuous: The balance sends a measuring result after every display cycle, but in no less than every 0.125 s.

Measurement values that are stable are identified with `[[S L]]`, those that are not with `[[SD]]`.

If the balance cannot supply a meaningful result (e.g. in overload or underload), an `[[SI]]` is sent.

If the TRF button is activated during this transfer mode, the balance sends the next stable result once with the identification `[[L L]]`.

The transfer can be recognized by the brief blanking-out of the display.

If the balance is tared during this transfer mode, the message `[[TA]]` will be sent automatically after the tare process is completed, immediately before CR and LF. After the balance is switched on, `[[TA]]` is also transmitted.

C. WHAT IF . . . ?

(see also section "WHAT IF. . . ?") in the Operating Instructions of the balance.

- | | |
|--|---|
| ... the balance supplies no results? | <ul style="list-style-type: none">– Incorrect voltage is selected on the housing of the data interface.– Data input was put out of order with the selection of the "TALKER ONLY" mode. (see Installation Instructions of Option 013).– An incorrect Device Address was selected.– An incorrect Operating Mode was selected.– No concluding symbol, or an incorrect concluding symbol, was used.– Position of slide (1) at mode switch was changed (correct position: "OFF"). |
| ... the data transfer is disrupted? | <ul style="list-style-type: none">– IEEE488 bus cables are too long.– Grounding concept should be checked. |
| ... the balance accepts no command? | <ul style="list-style-type: none">– Option 013 is configured to "TALKER ONLY". |
| ... the balance answers every command with «EL» | <ul style="list-style-type: none">– The balance is in the "OFF" mode. First send a «T» or a «C», or tare by hand. |
| ... the balance always sends «SI»? | <ul style="list-style-type: none">– The balance is outside of the admissible weighing range (overload or underload). |
| ... after the transfer button is activated, the balance does not transfer a measurement value? | <ul style="list-style-type: none">– The weighing pan did not stand still.– The transfer button is defective or incorrectly connected. |
| ... the bus cable does not fit the housing of the data interface? | <ul style="list-style-type: none">– Use a standard IEEE488 connector (<u>not</u> the DIN IEC 66.22 standard). |

D. ACCESSORIES

IEEE488 Bus cable

89181

Hand control

42500

Foot control

46278

E. Data transfer from 013 Option to personal computer (IEEE488/HP-IB)

Program examples with HP85 (Hewlett-Packard desktop computer)

To operate the IEEE488 interface one needs to observe a few basic facts. If attention is paid to these, error-free data transfer is quite straightforward. As a quick guide for the user of Option 013, two basic demonstration programs are listed below with the aim of illustrating these facts.

A. Option 013 with Service Request Interrupt (factory-set)

State of mode switches (1) inside the interface housing: (see Installation Instructions, page 1)

Switch 1: OFF (do not change, factory-set)

2: OFF Data exchange with Service Request

3: ON Hardware reset if malfunction occurs

4: OFF Send on Transfer mode

HP85 program sequence:

```
10 Z=715
20 CLEAR 715
30 OUTPUT Z USING "K" ; "SIR"
40 ON INTR 7 GOSUB 100
50 ENABLE INTR 7:8
60 !
65 GOTO 65
70 !
100 STATUS 7,1 ; S
110 A=SPOLL (Z)
120 ENTER Z ; S$
125 DISP S$
130 ENABLE INTR 7:8
140 RETURN
```

Key points in this program:

- Line 10 Definition of Select Code, incl. device address of Option 013 (factory setting: 15)
- Line 30 Sends a command to the balance (e.g. SIR) (not applicable when using Transfer button)
- Line 40 Opens the interrupt, states starting address of interrupt routine
- Line 50 Enables the interrupt, this is triggered by an activated SRQ line
- Line 65 This loop corresponds to the main program used in practice.
- Line 100 Interrupt routine: Interrogates the status byte of the computer interface, this resets the SRQ bit
- Line 110 Interrogates the status byte of Option 013, so enabling data transmission
- Line 120 Reads in the data string
- Line 140 Repeats interrupt enable (as line 50)

B. Option 013 without Service Request

State of mode switches (1) inside the interface housing: (see Installation Instructions, page 1)

Switch 1: OFF (do not change, factory-set)

2: ON Data exchange without Service Request

3: ON Hardware reset if malfunction occurs

4: OFF Send on Transfer mode

HP85 program sequence:

```
10 Z=715
20 CLEAR 715
30 OUTPUT Z USING "K" ; "SIR"
40 !
100 A=SPOLL (Z)
110 IF BIT (A,5)=0 THEN 100
120 ENTER Z ; S$
130 DISP S$
140 GOTO 100
```

Key points in this program:

- Line 10 Definition of Select Code, incl. device address of Option 013 (factory setting: 15)
- Line 30 Sends a command to the balance (e.g. SIR) (not applicable when using Transfer button)
- Line 100 Interrogates status byte of Option 013; if data are present (see next step) transmission is enabled
- Line 110 Tests bit 5 of status byte; this bit indicates whether data are ready, otherwise return to status interrogation at line 100.

The two programs given above cover the two most common practical cases. With other HP computers, use the equivalent commands with the same meaning.

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