
**Operating instructions**  
**Bedienungsanleitung**  
**Mode d'emploi**  
**Instrucciones de manejo**  
**Istruzioni per l'uso**

**METTLER** **TOLEDO**

**METTLER TOLEDO**  
**33360**  
**210260**

**Density determination kit**

Page 3

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


$$\rho_2 = \frac{A}{P} \cdot \rho_0 \text{ (g/cm}^3\text{)}$$

# Density determination kit

to determine the density of solids and liquids with top-loading electronic METTLER TOLEDO balances with pan diameters of 80 and 130 mm.

- for the density determination of solids, density determination kit 33360.
- for the density determination of liquids, density determination kit 33360  
+ displacement body 210260.

## 1. Definition of density

- to determine the density of solids and liquids with top-loading electronic METTLER TOLEDO balances with pan diameters of 80 and 130 mm.

The density  $\rho$  of an homogeneous body is the relationship of its mass (m) to its volume (V) or its mass per unit of volume.

$$\rho = \frac{m}{V} \quad \begin{matrix} [g] \\ [cm^3] \end{matrix}$$

$\rho$  = Density  
 $m$  = Mass  
 $V$  = Volume

According to DIN (German Engineering Standards) 1305, the term "Weight" can be used instead of "Mass" (m).

## 2. Principle of density determination

According to the Archimedean Principle, a solid body immersed in a liquid apparently loses as much of its own weight as the weight of the liquid it has displaced. This makes it possible to determine the unknown value. Depending on whether the body in question is liquid or solid, a slightly different procedure is involved.

The **density of a liquid** is determined by using a sinker of known **volume**. To this end, the sinker is first weighed while in air, and then while immersed. From these two weighings (made in grams), density  $\rho_1$  of the liquid is calculated as follows:

$$\rho_1 = \frac{A_s - B_s}{V_s} + L_s$$

$\rho_1$  = Density of liquid to be tested at a given temperature T

$A_s$  = Weight of sinker in air

$B_s$  = Weight of sinker when immersed in liquid

$V_s$  = Volume of sinker

$L_s$  = Air buoyancy per ml of sinker (correction of  $A_s$  or  $P_s$ ): about + 0,001 g/ ml; see also Section 3.2, Air Buoyancy

$P_s$  = Buoyancy of sinker in liquid (instead of  $A_s - B_s$ ): can be read directly on METTLER TOLEDO electronic balances; see Section 4.1

or simply

$$\rho_1 = \frac{P_s}{V_s} + L_s$$

The **density of a solid body** is determined by using a liquid of known **density**. To this end, the solid body is first weighed in air and then immersed. From these two weighings (made in grams or in carats), density  $\rho_2$  is calculated as follows:

$$\rho_2 = \frac{A}{A - B} \cdot \rho_o$$

$\rho_2$  = Density of solid body

$A$  = Weight of solid body in air

$B$  = Weight of solid body when immersed in test liquid

$\rho_o$  = Density of test liquid at a given temperature T

$P$  = Buoyancy of solid body in test liquid (instead of  $A - B$ ): can be read directly on METTLER TOLEDO electronic balances; see Section 4.2

or simply

$$\rho_2 = \frac{A}{P} \cdot \rho_o$$

The air buoyancy of a solid body is not always taken into account. If necessary, result  $\rho_2$  can be corrected by about + 0.001 g/ cm<sup>3</sup> similar to what was done in the first case; see also Section 3.2, Air Buoyancy.

### **3. Accuracy of density determination**

The required accuracy of the result determines whether any factors will be taken into account. The following listing should make it possible to quantitatively evaluate the influence of these factors.

### 3.1 Temperature

In the case of solid bodies, the change in density caused by a change in temperature is generally so small that the temperature of a solid body can be disregarded as far as density determination is concerned (this also applies to sinkers).

In the case of liquids, however, the density may change on the order of 0.1 ... 1 % pro °C, and this may already appear at the third decimal place.

Example:	Distilled water: Hydrocarbons + alcohols:	Density change about 1 ‰ per 5 °C Density change about 1 ‰ per 1 °C
----------	--	--

It should be noted that such a change in the density of a liquid is directly taken into account when the density of **solid bodies** is being determined because of the test liquid in which the solid body has to be immersed. For this reason, the temperature of the liquid must always be taken into consideration when the density is to be determined with better than 1 % accuracy.

### **Density table for distilled water:**

(according to "Handbook of Chemistry and Physics" 66th Ed. 1985-1986, F4-F5)

Temperature (°C)	Density (g/ ml)
15.0	0.9991
15.5	0.9990
16.0	0.9990
16.5	0.9989
17.0	0.9988
17.5	0.9987
18.0	0.9986
18.5	0.9985
19.0	0.9984
19.5	0.9983
20.0	0.9982
20.5	0.9981
21.0	0.9980
21.5	0.9979
22.0	0.9978
22.5	0.9977
23.0	0.9975
23.5	0.9974

Temperature (°C)	Density (g/ ml)
24.0	0.9973
24.5	0.9972
25.0	0.9970
25.5	0.9969
26.0	0.9968
26.5	0.9966
27.0	0.9965
27.5	0.9964
28.0	0.9962
28.5	0.9961
29.0	0.9959
29.5	0.9958
30.0	0.9956
30.5	0.9955
31.0	0.9953
31.5	0.9952
32.0	0.9950

For all other liquids, the density at temperature T must be taken from a book of tables.

### 3.2 Air Buoyancy

Depending on its physical conditions, each cubic centimeter ( $\text{cm}^3$ ) of air weighs 1 ... 1.2 mg. Thus, any object that is being weighed in air is subject to this kind of buoyancy for each  $\text{cm}^3$  of its volume. This means that – with a density of 1 g/  $\text{cm}^3$  – an error of approximately 0.1 % would occur if the air buoyancy is not taken into consideration. If a result with 3 or 4 decimal places is required, the result will have to be corrected for air buoyancy.

The true density is about 0.001 g/  $\text{cm}^3$  more than the calculated density.

When the density of liquids is determined, the operating procedure (division by 10 ml, corresponding to the volume of the sinker), provides a fourth-place result, even if the balance display shows only a three-place indication. It is therefore advisable – in this case – to make a general buoyancy correction.

### 3.3 Volume tolerance of the float

German Weights and Measures Regulation EO 13-4 paragraph 9.21

The volume of the float, together with the lower portion of the suspension wire, must be adjusted so that a 30 g float arrangement does not produce a measuring error exceeding  $\pm 0.00059 \text{ g/cm}^3$  when determining the water density at a temperature of 20 °C.

### 3.4 Immersion depth of the sinker or of the gem holder

The **sinker** (immersion body) is suspended from a platinum wire which has a diameter of 0.2 mm. In the water, the wire thus experiences a buoyancy of 0.3 mg when 10 mm of the wire are immersed.

If the fluid level stands 10 mm above the eylet of the sinker, about 20 mm of the wire are immersed. With a density (of liquid) of about 1, this would result in a buoyancy of 0.6 mg. But since this would again have to be divided by 10 ml, this influence can be disregarded.

The submersible part of the **gem holder** is made of a wire that has a diameter of 0.8 mm. With a liquid density of approx. 1, this results in a buoyancy of about 5 mg for each 10 mm that is submerged. But since the gem holder remains submerged while the solid body is weighed in air, and since – with electronic balances – the immersion depth does not change from one weighing to the next (in spite of the difference in weight), the buoyancy of the gem holder remains constant and can thus be disregarded. Condition: Do not change the liquid level. (The change of the liquid level caused by submerging the solid body can – in most cases – be disregarded. A solid body with a volume of 1  $\text{cm}^3$  causes the liquid to rise by about 0.5 mm. This is equivalent to a buoyancy of about 0.15 mg, i.e. a density error of 0.15 mg/  $\text{cm}^3$ ).

### 3.5 Surface tension of liquid

Since the liquid adheres to the suspension wire, an apparent weight increase occurs. In the case of the sinker (wire diameter 0.2 mm), and when water is used as the liquid, this force amounts to about 1 mg. By using wetting agents or organic liquids, this force can be reduced to 0.3 mg. However, this value will also be divided by 10 ml so that the resulting density error amounts to no more than 0.0001 g/ ml.

Because of the larger diameter of the wire submerged in the water, a force of up to 3 mg acts upon the gem holder. Here, too, similar to what is described in Section 3.4, any influence on the result will be virtually eliminated by having the gem basket submerged during both weighings (A and B). For very high accuracy requirements, reducing the surface tension would constitute an additional precautionary measure (see also Section 3.6).

### **3.6 Air bubbles**

With poorly wetting liquids such as water without wetting agents, for example, it is possible that air bubbles will adhere to the submerged solid body, sinker or wire basket of gem holder. Because of their buoyancy, these bubbles could affect the result. Thus, a bubble 1 mm in diameter would cause a buoyancy of 0.5 mg, while a bubble with a diameter of 2 mm would result in a buoyancy of as much as 4 mg.


#### **Precautionary measures:**

- Degrease solids that are resistant to solvents.
- Periodically clean gem holder and sinker; do not touch immersible parts with your hands.
- Gently shake gem holder when first immersing it in liquid, before suspending it from hook, so as to loosen air bubbles that might stick to it.
- Use wetting agents or organic solvents as auxiliary liquid (e.g. Kodak Photo-Flo, Pervitro 75% 72409).  
(The change of density caused by the addition of a wetting agent to distilled water is negligible. If, for example, 0.1 ml wetting agent with a density of 1.2 is added to 250 ml of water, the overall density changes by 0.001 g/ml.)

### **3.7 Porous bodies**

Since the submersion of porous bodies does not generally cause a one hundred percent displacement of air within the pores by the liquid, errors will occur as a consequence. Thus the density of the body can only be approximated.

## 4. How to perform density determinations

### Arrangement

- 1 Weighing pan of METTLER TOLEDO balance
- 2 Bracket attached to weighing pan ...
- 3 ... clamping screws hold bracket in place
- 4 Bridge placed on balance housing across weighing pan.  
Warning: Weighing pan 1 and bracket 2 must not brush against bridge:
- 5 250 ml beaker placed on bridge
- 6 Sinker or gem holder ...
- 7 ... suspension; the sinker or gem holder should hang freely in center of beaker without touching.  
Please note: The first weighing is dry!
- 8 Thermometer

The listing of the various parts is repeated on page 33 (Parts List)

### Reminders:

- It is recommended to check the calibration.
- Along with the text below, please check the illustrated instructions which are enclosed with this description.

### 4.1 Determining the density of liquids

- Use displacement body 210260 + density determination kit for liquids 33360.
- Place empty beaker on bridge 4, centered below suspension 7; insert thermometer 8.
- Attach sinker 6. It must not touch the beaker or the thermometer.
- Tare balance to show exactly zero.
- Pour liquid to be tested into beaker up to about 10 mm above the eyelet of the sinker. No air bubbles should adhere to the sinker. If necessary, brush off bubbles with a fine brush or the like.
- The balance displays (with minus sign in front) buoyancy  $P_s$  of the sinker. Value  $P_s$  displayed by the balance must now be divided by 10. After that, add 0.001 g/ ml.

The result is density  $\rho_1$  (g/ ml) at a given temperature T (to be read).

$$\rho_1 = \frac{P_s}{V_s} + L_s = \frac{P_s \text{ [g]}}{10 \text{ ml}} + 0.001 \text{ g/ ml at } T \text{ }^\circ\text{C} \quad (P_s = A_s - B_s)$$

If the direct buoyancy indication is to be dispensed with,  $P_s$  can be obtained from two separate weighings:  $A_s$  = sinker in air;  $B_s$  = sinker in liquid. Procedure is the same as above, but the balance is tared before the sinker is attached. With the sinker attached, the balance then indicates value  $A_s$ . After the liquid is filled in, the balance indicates (without additional taring!) value  $B_s$ .

**Example with buoyancy indication:** Determining the density of acetone.

Measured values:	$P_s = 7.898 \text{ g}$	$T = 20 \text{ }^\circ\text{C}$	$\rho_1 = \frac{7.898 \text{ g}}{10 \text{ ml}}$	$+ 0.001 \text{ g/ ml} = 0.7908 \text{ g/ ml at } 20 \text{ }^\circ\text{C}$
------------------	-------------------------	---------------------------------	--	--

The density of the acetone tested amounts to 0.7908 g/ ml at 20 °C.

## 4.2 Determining the density of solids

- Use density determination kit for solids (= gem holder **6**) 33360.
  - Pour enough auxiliary liquid into beaker so that the solid body after being placed in wire basket of gem holder will be covered with at least 10 mm of liquid. Insert thermometer. Place beaker on bridge **4** centrally located below the suspension **7**.
  - Attach gem holder (watch for air bubbles, especially at the wire basket).
  - Tare balance to read exactly zero.
  - Place dry solid in upper cup of gem holder. Write down weight displayed by the balance (in grams or carats). This is weight A of the solid when weighed in air.
  - Again tare balance (with solid body in cup).  
Balance display indicates exactly zero.
  - Remove solid body from cup and place it in lower wire basket (in liquid).  
Balance now displays buoyancy P of solid body (with minus sign in front).

Divide first indication A by second indication P, then multiply this intermediate result by the density  $\rho_0$  of the auxiliary liquid (at the given temperature T).

The final result is density  $\rho_2$  of the solid body. The unit of measurement is the same as that used for density  $\rho_0$  of the auxiliary liquid, e.g. g/cm<sup>3</sup> (identical to g/ml).

$$\rho_2 = \frac{A}{P} \cdot \rho_o = \frac{A \text{ [g]}}{P \text{ [g]}} \cdot \rho_o \text{ [g/cm}^3\text{]}$$

Since the units of A and P cancel each other, the weighing can be performed just as easily in carats as in grams. Result  $\rho_2$  always has the same unit of measurement as result  $\rho_o$ .

If the direct buoyancy reading is to be dispensed with, P can also be obtained from A - B. The balance will display B instead of P if it has not been tared after result A was displayed or if it was tared with the solid body **removed**.

**Example with buoyancy indication:** Determining the density of a coin with the help of distilled water.

$$\text{Measurement values: } A = 3.011 \text{ g} \\ P = 0.336 \text{ g} \\ T = 25.5^\circ \text{C} \rightarrow \rho_0 = 0.997 \text{ g/cm}^3$$

$$\rho_2 = \frac{3.011}{0.336} \cdot 0.997 \text{ g/cm}^3 = 8.93 \text{ g/cm}^3$$

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




For fast, reliable and efficient determination of density of solid, porous or fluid substances we also recommend the METTLER TOLEDO LabWare "Density".

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To protect your METTLER TOLEDO product's future:

METTLER TOLEDO Service assures the quality, measuring accuracy and preservation of value of all METTLER TOLEDO products for years to come.

Please send for full details about our attractive terms of service.

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