# Agilent 2100 Bioanalyzer

Installation and Safety Guide

**Agilent Technologies**

# Notices

© Agilent Technologies, Inc. 2017

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

## Manual Part Number

G2938-90007 Rev. D

## Edition

02/2017

Printed in Germany

Agilent Technologies  
Hewlett-Packard-Strasse 8  
76337 Waldbronn

## For Research Use Only

Not for use in diagnostic procedures

## Warranty

**The material contained in this document is provided “as is,” and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.**

## Technology Licenses

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with the terms of such license.

## Restricted Rights Legend

If software is for use in the performance of a U.S. Government prime contract or sub-contract, Software is delivered and licensed as “Commercial computer software” as defined in DFAR 252.227-7014 (June 1995), or as a “commercial item” as defined in FAR 2.101(a) or as “Restricted computer software” as defined in FAR 52.227-19 (June 1987) or any equivalent agency regulation or contract clause. Use, duplication or disclosure of Software is subject to Agilent Technologies’ standard commercial license terms, and non-DOD Departments and Agencies of the U.S. Government will

receive no greater than Restricted Rights as defined in FAR 52.227-19(c)(1-2) (June 1987). U.S. Government users will receive no greater than Limited Rights as defined in FAR 52.227-14 (June 1987) or DFAR 252.227-7015 (b)(2) (November 1995), as applicable in any technical data.

## Safety Notices

### CAUTION

A **CAUTION** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.


### WARNING

A **WARNING** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a **WARNING** notice until the indicated conditions are fully understood and met.

# Contents

<b>1</b>	<b>Safety Information</b>	<b>5</b>
	Power Cords	7
	Operation	8
	Safety Symbols	9
	Laser Safety	10
	Chemical and Biological Safety	11
	Legal Notice	13
<b>2</b>	<b>Site Requirements for the Agilent 2100 Bioanalyzer</b>	<b>75</b>
	Power Considerations	76
	Bench space	77
	Environment	78
	Sound Emission	79
	Physical Specifications of the Agilent 2100 Bioanalyzer	80
<b>3</b>	<b>Unpacking the Agilent 2100 Bioanalyzer System</b>	<b>81</b>
	Damaged Packaging	82
<b>4</b>	<b>Installing the Agilent 2100 Bioanalyzer System</b>	<b>85</b>
	Setting up your PC	86
	Setting up the Agilent 2100 Bioanalyzer	87
<b>5</b>	<b>Starting the Agilent 2100 Bioanalyzer</b>	<b>91</b>
	Turning on the Agilent 2100 Bioanalyzer	92
	Connecting the Bioanalyzer via the USB/Serial Cable	93
	Starting the Agilent 2100 Bioanalyzer Software	94
	Preparing the Assay	95
<b>6</b>	<b>Maintenance of the Agilent 2100 Bioanalyzer</b>	<b>97</b>
<b>7</b>	<b>Spare Parts and Accessories</b>	<b>99</b>

## Contents

# 1 Safety Information

Power Cords	7
Operation	8
Safety Symbols	9
Laser Safety	10
Chemical and Biological Safety	11
Legal Notice	13

The following general safety precautions must be observed during all phases of operation, service, and repair of the Agilent 2100 Bioanalyzer.

All safety instructions should be read and understood before installation, operation and maintenance of the instrument. Failure to comply with these precautions or with specific warnings elsewhere in this manual violates safety standards of design, manufacture and intended use of the instrument.

Agilent Technologies assumes no liability for the customer's failure to comply with these requirements.

## WARNING

**A WARNING notice denotes a hazard.**

**It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death.**

→ Do not proceed beyond a WARNING notice until the indicated conditions are fully understood and met.

## 1 Safety Information

### CAUTION

A CAUTION notice denotes a hazard.

It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data.

→ Do not proceed beyond a CAUTION notice until the indicated conditions are fully understood and met.

---

### NOTE

Instrument Safety: This is a Safety Class I instrument (provided with a terminal for protective grounding) and has been manufactured and tested according to international safety standards.

---

## Power Cords

Different power cords are offered as options with the Agilent 2100 Bioanalyzer. The female end of all power cords is identical. It plugs into the power-input socket at the rear of the instrument. The male end of each power cord is different and designed to match the wall socket of a particular country or region.

### WARNING

#### Electric Shock

**The absence of ground connection and the use of an unspecified power cord can lead to electric shock or short circuit.**

- Never operate your instrument from a power outlet that has no ground connection.
  - Never use a power cord other than the Agilent Technologies power cord designed for your region.
  - Never use cables other than those supplied by Agilent Technologies to ensure proper functionality and compliance with safety or EMC regulations.
-

## Operation

### WARNING

**Ensure the proper usage of the equipment**

**The protection provided by the equipment may be impaired.**

→ The operator of this instrument is advised to use the equipment in a manner as specified in this manual.

---

Before the instrument is switched on, all protective ground terminals, extension cords, auto-transformers, and devices connected to it must be connected to a protective ground socket. Any interruption of the protective grounding will cause a potential shock hazard that could result in serious personal injury. Whenever it is likely that the protection has been impaired, the instrument must be made inoperative and be secured against any intended operation.

To operate the instrument safely:

- Do not remove any cover of the Agilent 2100 Bioanalyzer
- Avoid any maintenance of the instrument under voltage
- Do not let liquid drip into the Agilent 2100 Bioanalyzer. It could cause a shock or it could damage the Agilent 2100 Bioanalyzer.
- Do not replace components of the instruments (e.g. electrode cartridge) with power turned on.
- Do not operate the instrument in the presence of flammable gases or fumes. Operation of any electrical instrument in such an environment constitutes a definite safety hazard.
- Do not install substitute parts or make any unauthorized modification to the instrument.

# Safety Symbols

Symbol	Description
	Electrosensitive device.
	The apparatus is marked with this symbol when the user should refer to the <i>Installation and Safety Manual</i> in order to prevent risk of harm to the operator and to protect the apparatus against damage.
	Indicates dangerous voltages.
	Indicates a protected ground terminal.
	Laser Radiation. When open avoid exposure to beam.

## Laser Safety

The Agilent 2100 Bioanalyzer is classified as a "Laser Class 1" product (IEC825-1, CFR1040.10). During normal operation of the Agilent 2100 Bioanalyzer no laser light is accessible to the user.

When the laser is in use, the laser light source emits light at a power that may be harmful to the eyes. To prevent operator exposure to potentially harmful laser light, an interlock mechanism turns off the laser whenever the lid of the Agilent 2100 Bioanalyzer is not closed. Furthermore, the lid can't be completely opened and has a black light absorbing surface in order to prevent any potential laser light reflection.

### **WARNING**

#### **Harmful laser light**

**The laser light source emits light at a power that may be harmful to the eyes.**

- Never look into the beam or direct it towards someone else.
  - System interlocks should never be disabled.
-

## Chemical and Biological Safety

### WARNING

#### Toxic and hazardous solvents

**The handling of solvents and reagents can hold health risks.**

- When working with solvents observe appropriate safety procedures (for example, goggles, safety gloves and protective clothing) as described in the material handling and safety data sheet supplied by the solvent vendor, especially when toxic or hazardous solvents are used.
- 

### WARNING

#### Pathogenic, toxic, or radioactive samples

**Handling and use of pathogenic, toxic, or radioactive samples and of genetically modified organisms holds risks for health and environment.**

- Ensure that all necessary safety regulations, guidelines, precautions and practices are adhered to accordingly.
  - Consult the laboratory safety officer for advise on the level of containment required for the application, and proper decontamination or sterilization procedures to follow if fluids escape from containers.
- 

To operate the instrument safely:

- Observe all cautionary information printed on the original solution containers prior to their use.
- Because leaks, spills, or loss of sample may generate aerosols, observe proper safety precautions.
- Agilent 2100 Bioanalyzer covers have not been designed as bioseals for aerosol or liquid containment.
- Handle body fluids with care because they can transmit disease. No known test offers complete assurance that they are free of micro-organisms. Some of the most virulent – Hepatitis (B and C) and HIV (I-V) viruses, atypical mycobacteria, and certain systemic fungi – further emphasize the need for aerosol protection.

## 1 Safety Information

### Chemical and Biological Safety

- Always follow local state and federal biohazard handling regulation when disposing of biohazardous waste material.
- Handle all infectious samples according to good laboratory procedures and methods to prevent spread of disease.
- Dispose of all waste solutions and products according to appropriate environmental health and safety guidelines.

#### **WARNING**

#### **Harmful chemical and biological substances**

**Residues of clinical samples may contain chemicals and biological substances that are dangerous and harmful to persons working on the 2100 Bioanalyzer.**

- The Agilent 2100 Bioanalyzer and any accessories must be decontaminated before requesting service by an Agilent Technologies Field Service Representative and before returning the 2100 Bioanalyzer to Agilent Technologies for repair or replacement.
-

## Legal Notice

The Agilent 2100 Bioanalyzer is sold for research use only and not intended for diagnostic procedures.

By purchasing this instrument, the purchaser is granted the limited right to use only this instrument. Purchase of this instrument does not include any right, express or implied, to use any other patented product, method or process, or to use any other portion or component of any patented system or systems, software, microfluidic devices or reagents, either alone or in conjunction with this product, unless use of such method, process or other portion or component is separately authorized.

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**1 Safety Information**  
Legal Notice

## 2

# Site Requirements for the Agilent 2100 Bioanalyzer

Power Considerations 76

Bench space 77

Environment 78

Sound Emission 79

Physical Specifications of the Agilent 2100 Bioanalyzer 80

## Power Considerations

The Agilent 2100 Bioanalyzer power supply has wide ranging capabilities and accepts any line voltage in the range 100 - 240 V  $\pm$ 10% with a line frequency of 50 - 60 Hz  $\pm$ 5% (see also [Table 1](#) on page 80 for more details). Consequently, there is no voltage selector in the rear of the instrument.

### CAUTION

Damage through high line voltage

Connecting the Agilent 2100 Bioanalyzer to a line voltage higher than specified may result in damage of the instrument.

→ Connect the Agilent 2100 Bioanalyzer to the specified line voltage only.

---

## Bench space

The Agilent 2100 Bioanalyzer dimensions and weight (see [Table 1](#) on page 80) allow you to place it on almost any desk or laboratory bench. It needs approximately 8 cm (3.1 inches) of clear space at the rear for air circulation and electric connections.

## Environment

Your Agilent 2100 Bioanalyzer will work within the specifications of ambient temperatures and relative humidity described in [Table 1](#) on page 80.

### CAUTION

#### Condensation

Condensation will damage the electronics.

- Do not store, ship or use the instrument under conditions where temperature fluctuations could cause condensation within the instrument.
  - If your instrument was shipped in cold weather, leave it in its box and allow it to warm up slowly to room temperature to avoid condensation.
- 

### NOTE

The Agilent 2100 Bioanalyzer is designed to operate in a controlled electromagnetic environment (EN61326/A1) where RF transmitters such as mobile telephones should not be used in close proximity.

---

### NOTE

In order to function properly, do not place the Agilent 2100 Bioanalyzer on a vibrating surface or near vibrating objects.

---

## Sound Emission

### Manufacturers Declaration

This statement is provided to comply with the requirements of the German Sound Emission Directive of 18 January 1991.

This product has a sound pressure emission (at the operator position) < 70 dB.

- Sound Pressure  $L_p < 70$  dB (A)
- At Operator Position
- Normal Operation
- According to ISO 7779:1988/EN 27779/1991 (Type Test)

## Physical Specifications of the Agilent 2100 Bioanalyzer

Table 1 on page 80 lists the physical specifications of the Agilent 2100 Bioanalyzer.

**Table 1** Physical specifications

Type	Specification	Comment
Weight	10 kg (22 lbs)	
Dimensions	162 — 412 — 290 mm (6.4 — 16.2 — 11.4 in)	Width — depth — height
Line voltage	100-240 VAC ±10%	Wide-ranging capability
Line frequency	50-60 Hz ±5%	
Power consumption	30 W / 60 VA	Maximum
Ambient operating temperature	5-40 °C (41-104 °F)	The assays set-up to be analyzed on a chip temperature of 30 °C. If the ambient temperature is above 30 °C the assays are analyzed using that temperature.
Ambient non-operating temperature	-40-70 °C (-40-158 °F)	
Humidity	<70%, at 15-27 °C (59-81 °F)	Non-condensing
Operating altitude	Up to 2000 m (6500 ft)	
Non-operating altitude	Up to 4600 m (14950 ft)	For storage of the Agilent 2100 Bioanalyzer
Safety standards: IEC, CSA, UL	Installation Category II, Pollution Degree 2	

### 3

## Unpacking the Agilent 2100 Bioanalyzer System

Damaged Packaging	82
Delivery Checklist	82

## Damaged Packaging

Upon receipt of your Agilent 2100 Bioanalyzer and a potential computer inspect the shipping containers for any signs of damage. If containers or cushioning material are damaged, save them until the contents have been checked for completeness and the 2100 Bioanalyzer and computer have been mechanically and electrically checked. If the shipping container or cushion material is damaged, notify the carrier as well as Agilent Technologies. Save the shipping material for the carrier's inspection.

### CAUTION

If there are signs of damage to the Agilent 2100 Bioanalyzer:

- Do not attempt to install the instrument.
- Contact your local Agilent Technologies sales and service office.

## Delivery Checklist

Ensure that all parts and materials have been delivered with your 2100 Bioanalyzer.

Please report any missing or damaged parts to your local Agilent Technologies sales and service office.

**Table 2** Agilent 2100 Bioanalyzer (G2939B) delivery checklist

Description	Quantity
Agilent 2100 Bioanalyzer instrument	1
Electrode cartridge	1
Test chip bundle	1
Declaration of conformity	1
Installation and Safety Manual	1
RS-232 Connector cable	1
USB/serial adapter cable	1

**Table 2** Agilent 2100 Bioanalyzer (G2939B) delivery checklist

<b>Description</b>	<b>Quantity</b>
Fuses	2
Chip priming station	1
IKA vortex mixer	1
Agilent 2100 Bioanalyzer Expert Software, including declaration of system validation.	1
Instrument control license and electrophoresis license	1

### **3 Unpacking the Agilent 2100 Bioanalyzer System**

#### **Damaged Packaging**

## 4

# Installing the Agilent 2100 Bioanalyzer System

Setting up your PC 86

Setting up the Agilent 2100 Bioanalyzer 87

After unpacking and checking the completeness of the shipment, the Agilent 2100 Bioanalyzer system is ready to be installed. Depending on system configuration, the installation requires up to three steps:

- 1 Setting up the PC.
- 2 Setting up the Agilent 2100 Bioanalyzer.

### NOTE

If you have problems installing your Agilent 2100 Bioanalyzer, refer to the Troubleshooting section of the *Maintenance and Troubleshooting Guide*.

## **4** Installing the Agilent 2100 Bioanalyzer System

### Setting up your PC

# Setting up your PC

For setting up the PC, please refer to the instructions supplied with the PC.

## Setting up the Agilent 2100 Bioanalyzer

To set-up the Agilent 2100 Bioanalyzer instrument:

- 1 Locate bench space that meets the requirements described in “[Bench space](#)” on page 77.
- 2 Place the instrument on the bench in an upright position.

### CAUTION

Damage by electrostatic discharge

The heater plate and electrodes are sensitive to electrostatic discharge.

- Make sure that you are grounded before inserting a chip to prevent the Agilent 2100 Bioanalyzer from being damaged by electrostatic discharge. The lid is grounded to leak off electrostatic energy of the operator.
-


## 4 Installing the Agilent 2100 Bioanalyzer System

### Setting up the Agilent 2100 Bioanalyzer

- 3 Open the lid of the 2100 Bioanalyzer and slide the electrode cartridge into the lid and move the lever into a flat (closed) position.

- 4 Connect one end of the serial cable to the RS-232 port located at the rear of the instrument.

- 5 Connect the other end of the serial cable to the serial port (COM-port) of your PC.

- 6 If your PC is not equipped with a serial port, use the USB/serial adapter cable that is provided with the instrument.
- 7 Ensure that the power switch located at the rear of the instrument is in off position.


**NOTE**

The power input socket accepts a line voltage of 100-240 VAC  $\pm$ 10% with a line frequency of 50-60 Hz. Maximum power consumption is 60 VA. There is no voltage selector at the rear of the Agilent 2100 Bioanalyzer because the power supply has a wide-ranging capability.

- 8 Plug the female end of the power cord into the power input socket of the Agilent 2100 Bioanalyzer instrument. Plug the male end of the power cord into the electrical outlet.

## **4 Installing the Agilent 2100 Bioanalyzer System**

### **Setting up the Agilent 2100 Bioanalyzer**

## 5 Starting the Agilent 2100 Bioanalyzer

Turning on the Agilent 2100 Bioanalyzer 92

Connecting the Bioanalyzer via the USB/Serial Cable 93

Starting the Agilent 2100 Bioanalyzer Software 94

Preparing the Assay 95

When starting the 2100 Expert Software, wait a few seconds until the self-test procedures for the Agilent 2100 Bioanalyzer are finished. Otherwise the software may not recognize the Agilent 2100 Bioanalyzer.

### CAUTION

Condensation within the instrument

Temperature fluctuations could cause condensation within the instrument.

- Do not use the Agilent 2100 Bioanalyzer under conditions with temperature fluctuations.
- Let the instrument equilibrate to room temperature for one day to avoid damage of the electronics.

## Turning on the Agilent 2100 Bioanalyzer

- 1 Turn on your Agilent 2100 Bioanalyzer. The power switch is located at the rear of the instrument, where the power cable plugs in (see step 4 on page 88).
- 2 The status indicator lamp will light green if power is present and all instrument self-tests have passed successfully.

Status Indicator

- Green - Ready
- Green, flashing - Busy
- Orange - Self-test
- Red - Error

If the status indicator lamp lights red, this indicates a hardware error. In this case, please refer to the Troubleshooting section of the *Agilent 2100 Bioanalyzer Maintenance and Troubleshooting Guide* or the *Online Help*.

- 3 Once your Agilent 2100 Bioanalyzer has passed the self- tests, it is ready to run samples.
- 4 Start your PC by pushing the appropriate power buttons as outlined in the documentation.

## Connecting the Bioanalyzer via the USB/Serial Cable

If you connect the 2100 Bioanalyzer instrument via the USB/serial adapter cable, the driver needs to be installed on the PC. This is preinstalled on bundle PCs but needs to be installed manually on 3rd party PCs.

The driver can be found on the 2100 Expert Software CD. For installation instructions please refer to the readme file of the 2100 Expert Software.

## Starting the Agilent 2100 Bioanalyzer Software

The Agilent 2100 Bioanalyzer software is pre-installed on bundle PCs. However, the full functionality of the software becomes available only if you add the required license key(s). License keys are supplied with the 2100 Bioanalyzer.

- 1 Open the 2100 Expert Software by double-clicking the icon on the desktop.
- 2 Select **Help > Registration** and highlight the **Add License** tab.
- 3 Add the appropriate license key.

### NOTE

For more details of adding license keys, refer to the *Online Help* .

---

- 4 To establish instrument-PC communication, assign a COM-port in the instrument context. When the communication is established, instrument specific information like serial number, vendor information and firmware version is displayed.

### NOTE

For detailed information and in case of problems, please refer to the Troubleshooting section of the *Agilent 2100 Bioanalyzer Maintenance and Troubleshooting Guide* or the *Online Help*.

---

## Preparing the Assay

For information about the assay preparation, please refer to the corresponding reagent kit guide. All kit guides are included in the online help of the 2100 Expert software. The quick start guide is provided with every analysis kit. Use the reagent kit guides as a reference when preparing and running an assay with the Agilent 2100 Bioanalyzer.

The Agilent Technologies Genomics Web site offers support and useful information about current developments, products and technology:

<http://www.agilent.com/genomics/bioanalyzer>

## **5 Starting the Agilent 2100 Bioanalyzer**

### **Preparing the Assay**

## 6 Maintenance of the Agilent 2100 Bioanalyzer


For detailed maintenance procedures of the Agilent 2100 Bioanalyzer instrument and its accessories please refer to:

- the *Online Help* of the Agilent 2100 Bioanalyzer software or
- the Maintenance part of the *Troubleshooting and Maintenance Guide*

## **6 Maintenance of the Agilent 2100 Bioanalyzer**

### **Preparing the Assay**

## 7 Spare Parts and Accessories

The following table provides a list of spare parts and accessories that are available for the Agilent 2100 Bioanalyzer. To buy parts, please refer to the Agilent Online Store: <http://www.agilent.com/home/buyonline.html>

**Table 3** Spare Parts and Accessories for the Agilent 2100 Bioanalyzer

<i>Reorder number</i>	<i>Part</i>	<i>Description</i>
RS232-61601	RS-232 Cable	communication cable between PC and instrument
8121-1013	USB/Serial Adapter Cable	connects RS-232 cables to USB PC ports (for PCs without serial port)
2110-0007	Fuse	fuse for power supply
G2938-68716	Gasket Kit	contains spare parts for chip priming station: 1 plastic adapter, 1 ring and 10 gaskets
5042-1398	Adjustable Clip	used in combination with a syringe to apply defined pressure for chip priming
5065-9966	Vortex Mixer Adapter	replacement part for IKA MS2 Vortexer (3 mounting screws)
5065-4413	Electrode Cartridge	no extra electrode pin set; pin-set not separately re-orderable
5065-4401	Chip Priming Station	includes priming station, timer and 1 syringe clip
G2938-68300	Test Chip Kit for Electrophoresis Assays	comprises autofocus chip, electrode/diode test chip and documentation
5065-4428	IKA Vortexer MS3 type	Bundled to all Bioanalyzer systems but NOT available from Agilent standalone. Replacement vortexers must be ordered directly from IKA. Description: "MS 3 basic S36 Agilent" IKA Part Number: "3617036"

# Index

## A

accessories 99  
assay preparation 95

## B

bench space 77  
bioanalyzer software 83  
bioanalyzer  
    accessories 99  
    delivery checklist 82  
    maintenance 97  
    spare parts 99  
biological safety 11  
body fluids 11

## C

cable 82  
chemical safety 11  
chip priming station 83  
communication 94  
com-port 94

## D

declaration of conformity 82  
decontamination 12  
delivery checklist 82  
dimensions 80

## E

environment 78

## F

flammable gases 8  
frequency 76, 80  
fuse 83

## G

Genomics Web site 95

## H

humidity 80

## I

installation 85  
interlock mechanism 10

## K

kit guide 95

## L

laser safety 10  
laser 10  
legal notice 97  
license key 94  
license 94

## M

maintenance 97

## O

online store 99  
operating temperature 80

## P

packaging 82  
physical specifications 80  
power considerations 76  
power consumption 80  
power cord 7  
power supply 76

## R

reagent kit guide 95  
RS232 cable 82

## S

safety information 5  
safety standards 80  
safety symbols 9  
setup  
    bioanalyzer 87  
    PC 86  
software license 94  
software 83  
sound emission 79  
spare parts 99  
start  
    bioanalyzer software 94  
    bioanalyzer 91  
status indicator 92

**T**

turning on the Agilent 2100 Bioanalyzer 92

**V**

voltage 80

vortex mixer 83

**W**

weight 80

[www.agilent.com](http://www.agilent.com)

© Agilent Technologies 2017

Printed in Germany  
02/2017

G2938-90007  
Rev. D