HEAVY DUTY VORTEX MIXER SPECIFICATIONS

Overall dimensions (L x W x H): 24.1 x 16.8 x 16.5cm

Electrical (50/60 Hz): 230V | 0,13 amps | 30 watts Fuses: 5mm x 20mm, 5 amp 250 volt

quick acting

Speed range On Mode: 300-2500rpm

Touch Mode: 300-3500rpm

Orbit: 4.9mm

Duty rating: Continuous duty

Controls: Standard: Touch/off/on - rocker switch

Speed knob: variable 1 to 10 dial marks

Advanced: See page 5

Capacity: 1,1kg Ship weight: 6.8kg

MICROPLATE VORTEX MIXER SPECIFICATIONS

Overall dimensions (L x W x H): 26.6 x 13.7 x 11.4cm

Electrical (50/60 Hz): 230V | 0,13 amps | 30 watts 5mm x 20mm, 5 amp 250 volt Fuses:

quick acting

Speed range On Mode: 300-2500rpm

Touch Mode: 300-3500rpm

Orbit: 3.5mm

Continuous duty **Duty rating:**

Controls: Standard: Touch/off/on - rocker switch


Speed knob: variable 1 to 10 dial marks

Advanced: See page 5

Capacity: 2 Microplates


5,4kg Ship weight:

STANDARD VORTEX MIXER CONTROL PANEL



The front panel of the Standard Vortex Mixer contains all the switches and control knob needed to operate the unit.

- A. 3-way rocker switch: touch/off/on rocker switch starts/stops the vortexing function.
- B. Speed control: The speed control knob controls the vortexing speed. User should turn knob clockwise until the vortexer reaches the desired speed.



The front panel of the Advanced Vortex Mixer contains all the switches, controls and displays needed to operate the unit.

- **A. 3-way rocker switch:** touch/standby/on rocker switch starts/stops the vortexing function.
- **B. Time display:** Displays accumulated time (continuous mode) or how much time is remaining (timed mode). The display range is from 0 to 9,999 minutes in one (1) second increments. The display will indicate minutes and seconds until the timer reaches 99 minutes and 59 seconds (99:59), then the display will automatically display minutes up to 9,999.
- C. Up/down arrows for set-point control.
- **D. Speed display:** Displays the speed of the Vortex Mixer.
- **E.** Up/down arrows for set-point control.