2000005495

MS 3 basic/digital_082017

IKA®

MS 3 digital

MS 3 basic


BETRIEBSANLEITUNG	DE	3
OPERATING INSTRUCTIONS	EN	10
MODE D'EMPLOI	FR	17
INDICACIONES DE SEGURIDAD	ES	24
VEILIGHEIDINSTRUCTIES	NL	25
NORME DI SICUREZZA	IT	26
SÄKERHETSANVISNINGAR	SV	27
SIKKERHEDSHENVISNINGER	DA	28
SIKKERHETSANVISNINGER	NO	29
TURVALLISUUSOHJEET	FI	30
INSTRUÇÕES DE SERVIÇO	PT	31
WSKAZÓWKI BEZPIECZEŃSTWA	PL	32
BEZPEČNOSTNÍ POKYNY	CS	33
BIZTONSÁGI UTASÍTÁSOK	HU	34
VARNOSTNA OPOZORILA	SL	35
BEZPEČNOSTNÉ POKYNY	SK	36
OHUTUSUJUHISED	ET	37
DROŠĪBAS NORĀDES	LV	38
NURODYMAI DĒL SAUGUMO	LT	39

IKA-Werke, Germany

MS 3 basic

Pos.Bezeichnung	Designation	Désignation
MS 3 basic/digital		
A Taste "Power"	Button "Power"	Bouton "Power"
B Standardaufsatz	Standard attachment	Support standard
C Drehknopf Drehzahleinstellung	Rotary knob of speed adjustment	Bouton rotatif pour la réglage de la vitesse de rotation
D Gerätefuss	Base	Pied
E Taste "Start/Stop"	Button "Start/Stop"	Bouton "Start/Stop"
MS 3 digital		
F Tasten "Timer"	Buttons "Timer"	Boutons "Minuteur"
G Display	Display	Indication

MS 3 digital

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DE

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Contents

EN

	Page
Declaration of conformity	2
Safety instructions	10
Unpack	11
Correct use	11
Commissioning	12
Accessories	14
Maintenance	15
Technical data	16
Warranty	16

Safety instructions

- **Read the operating instructions in full before starting up and follow the safety instructions.**
- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- The voltage stated on the nameplate must correspond to the mains voltage.
- The device must only be operated with the original plug-in power supply unit.
- Wear your personal protective equipment in accordance with the hazard category of the medium to be processed. Otherwise there is a risk of:
 - splashing liquids
 - projectile parts
 - body parts, hair, clothing and jewellery getting caught.
- Set up the appliance in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- The feet of the appliance must be clean and undamaged.
- Position the knob (C) at the left stop before starting up. Gradually increase the speed.

- Keep hand/arm vibration to a minimum (Directive 2002/44/EC) when operated manually:
 - Keep the sample container as vertical as possible.
 - Only exert the minimum pressure necessary for the mixing process.
 - Use attachments for multiple samples when there is a large number of samples.
- If the shaking movement of the appliance increases (resonance), reduce the speed or pass through the critical phase as quickly as possible.
- Reduce the speed if
 - the medium splashes out of the vessel because the speed is too high
 - the appliance is not running smoothly
 - the appliance begins to move around because of dynamic forces.
- Firmly secure the accessories and vessels in place, otherwise shaking vessels could be damaged or projected out.
- Check the appliance and accessories beforehand for damage each time you use them. Do not use damaged components.
- Position one single shaking vessel in the centre and several shaking vessels so that they are evenly spread out.
- Sharp edged vessels cause wear debris on the attachments.
- Beware of the risk of
 - flammable materials
 - glass breakage as a result of mechanical shaking power
- Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.
- **Do not** operate the appliance in explosive atmospheres, with hazardous substances or under water.
- Safe operation is only guaranteed with the accessories described in the "Accessories" chapter.
- Always disconnect the plug before fitting accessories.
- The appliance does not start up again automatically following a cut in the power supply.
- Protect the appliance and accessories from bumps and impacts.
- The appliance may heat up when in use.
- The appliance may only be opened by experts.

Unpack

- **Unpack**

- Please unpack the device carefully
- In the case of any damage a fact report must be set immediately (post, rail or forwarder)

- **Delivery scope**

MS 3 basic

- Shaker MS 3 basic
- Power supply
- Standard attachment MS 3.1
- Universal attachment MS 3.3
- One-hand insert MS 1.21
- Operating instructions

MS 3 digital

- Shaker MS 3 digital
- Power supply
- Standard attachment MS 3.1
- Universal attachment MS 3.3
- One-hand insert MS 1.21
- Microtiter attachment MS 3.4
- Test tube insert MS 1.32
- Operating instructions

Correct use

- **Use**

- ▶ For mixing liquids
- Touch mode for single test tubes
- Continuous mode for single or several vessels




- **Range of use**



- Laboratories
- Pharmacies
- Schools
- Universities

This device is suitable for use in all areas except industrial areas.

The protection for the user is not ensured no more, if the equipment with accessories is operated, which are supplied or are not recommended not by the manufacturer or if the equipment in not correct use is operated against the manufacturer default.

Commissioning

Setting	Action	Display
	 <p>Position the rotating knob at the left stop before starting up</p>  <p>If these conditions are met, the unit is ready for service when the mains plug has been plugged in</p>	
<p>1 Switching on</p> <p>Change to continuous operation</p>	 <p>Mode A Operating mode "Touch operation" Speed upto 3000 rpm (= Factory setting)</p>  <p>Mode A Operating mode "Continuous operation" Speed upto 1300 rpm</p>	  
<p>2 Setting timer / Change timer setting</p> <p>Not available on MS 3 BASIC</p> <p>Change to minutes mode</p>	 <p>Seconds mode The value previously set will be displayed (Factory setting=000)</p> <p>Countdown running: An audible signal will be given when the countdown is complete</p>  <p>Switch off the device</p>   <p>Switching on the device: Press and hold button "1" and also press the "Power" button</p>  <p>Minutes mode</p>	  





Setting	Action	Display
		<p>Countdown running: An audible signal will be given when the countdown is complete</p>
3 Speed adjustment		
4 Change mode A → B		<p>Switch off the device</p> <p>Switching on the device: Press and hold button "Start/Stop" and also press the "Power" button</p> <p>Mode B Operation modes "Touch" and "Continuous operation" Speed upto 3000 rpm</p>   
B → A		<p>Switch off the device and switch on again, see 1</p>

Operating modes

	Modus A (with speed limiter)	Modus B (without speed limiter)
Touch mode <i>(with standard attachment MS 3.1)</i>	max. 3000 rpm	max. 3000 rpm
Continuous mode <i>(with all attachments)</i>	max. 1300 rpm	max. 3000 rpm


Accessories

Possible uses and permitted speed ranges of attachments:


Designation	Description	Ident-Nr.	Touch mode	Continuous mode	Speed range (rpm)	Mode
 MS 3.1 Standard attachment	For test tubes and small vessels upto ø 50 mm	3426300	x - -	- x x	0/200-3000 0/200-1300 0/200-3000	A and B A B
 MS 3.3 Universal attachment	For different rubber foam inserts	3426600	-	x	0/200-1300 0/200-3000	A B
 MS 3.4 Microtiter attachment	For one microtiter plate	3426400	-	x	0/200-1300	A
 MS 3.5 PCR-Plate attachment	For one 96-well PCR-plate	3428000	-	x	0/200-1300	A
 MS 3.51 PCR-Plate insert	- For inserting into the PCR-Plate attachment	3428700	-	x	0/200-1300	A
 MS 1.21 One-hand insert	For inserting into the universal attachment	L001540	-	x	0/200-1300 0/200-3000	A B
 MS 1.30 Test tube insert	- For inserting into the universal attachment - For 24 tubes ø 7,8 mm	25005776	-	x	0/200-1300	A
 MS 1.31 Test tube insert	- For inserting into the universal attachment - For 14 test tubes ø 10 mm	L001840	-	x	0/200-1300	A
 MS 1.32 Test tube insert	- For inserting into the universal attachment - For 6 test tubes ø 12 mm	L001850	-	x	0/200-1300	A
 MS 1.33 Test tube insert	- For inserting into the universal attachment - For 4 test tubes ø 16 mm	L001860	-	x	0/200-1300	A
 MS 1.34 Test tube insert	- For inserting into the universal attachment - You can make holes as you like	L001830	-	x	0/200-1300	A

Changing attachments


MS 3.1

MS 3.3
MS 3.4

Using inserts

Maintenance

Cleaning

The appliance is maintenance-free.

For cleaning, disconnect the main plug!

Only use cleansing agents which have been recommended by **IKA**. Use to remove:

Dyes	isopropyl alcohol
Construction materials	water containing tenside/ isopropyl alcohol
Cosmetics	water containing tenside/ isopropyl alcohol
Foodstuffs	water containing tenside
Fuels	water containing tenside

- Do not allow moisture to get into the appliance when cleaning
- Wear protective gloves during cleaning the devices.
- Before using another than the recommended method for cleaning or decontamination, the user must ascertain with **IKA** that this method does not destroy the instrument.

Spare parts order

When ordering spare parts, please give:

- Machine type
- Manufacturing number, see type plate
- Item and designation of the spare part, see spare parts diagram and spare parts list, see www.ika.com.

Repair

Please only send devices in for repair that have been cleaned and are free of materials which might present health hazards.

For this, use the “**certificate of compliance**” form which you can obtain from **IKA** or can download a version for printing from the **IKA** website at www.ika.com.

If your appliance requires repair, return it in its original packaging. Storage packaging is not sufficient when sending the device - also use appropriate transport packaging.

Technical data

Power supply

Input	V	100 - 240
	A	0,8
	Hz	50-60

Output	Vdc	24
		24W.LPS (Limited power source)
Protection class		2 (double insulated) <input type="checkbox"/>

Shaker

Operating voltage	Vdc	24
	mA	800
Power consumption, normal operation	W	20
Power consumption, standby operation	W	2
Motor-output power	W	8
Drive		EC - motor
Speed range	rpm	0/200 – 3000 infinitely adjustable
Speed adjustment		rotating knob in front of the appliance
Speed display		scale 0 – 3000
Agitation stroke	mm	4,5
Shaking motion		horizontal, circular
Perm duration of operation	%	100
Timer mode seconds		1 - 999 sec
mode minutes		1 - 999 min
Timer display		digital
Perm ambient temperature	°C	+5 bis +40
Perm. relative humidity	%	80
Protection type acc. to DIN EN 60529		IP 21
Contamination level		2
Overvoltage category		II
Operation at a terrestrial altitude	m	max. 2000

Dimensions total W x D x H	mm	148 x 205 x 63
without attachment		
Weight without supported load	kg	2,9
Max. supported load incl. attachment	kg	0,5

Subject to technical changes!

Warranty

In accordance with **IKA** warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

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IKA®-Werke GmbH & Co.KG