


JA-20.1 FIXED ANGLE ROTOR

SPECIFICATIONS



Maximum speed 20 000 rpm*
 Critical speed range† 600 to 800 rpm
 Density rating at maximum speed. 1.2 g/mL
 Maximum allowable imbalance of opposing loads. 2 grams
 Relative Centrifugal Field‡ at maximum speed

Outer row of tubes

At r_{max} (115 mm) $51\,500 \times g$
 At r_{max} with adapter 342327 (107 mm) $47\,900 \times g$
 At r_{min} (64 mm) $28\,700 \times g$

Inner row of tubes

At r_{max} (98 mm) $43\,900 \times g$
 At r_{max} with adapter 342327 (90 mm) $40\,300 \times g$
 At r_{min} (47 mm) $21\,100 \times g$

k factors at maximum speed

Outer row of tubes 371
 Outer row of tubes with adapter 342327 325
 Inner row of tubes 465
 Inner row of tubes with adapter 342327 411

Conditions requiring speed reductions see RUN SPEEDS

Number of tube cavities 32

Available tubes see Table 4

Nominal tube dimensions (largest tube) 18×99 mm

Nominal tube capacity (largest tube) 15 mL

Nominal rotor capacity 480 mL

Approximate acceleration time to maximum speed

(rotor fully loaded) $1\frac{1}{2}$ min

Approximate deceleration time from maximum speed

(rotor fully loaded) $2\frac{1}{2}$ min

Weight of fully loaded rotor 7.5 kg (16.5 lb)

Rotor material aluminum

Rotor entry code for Avanti J-E, J2 and J6 series

microprocessor-controlled centrifuges. 20.1

* Maximum speed at 2°C in the Avanti J-E 17 000 rpm

† The critical speed range is the range of speeds over which the rotor shifts so as to rotate about its center of mass. Passing through the critical speed range is characterized by some vibration.

‡ Relative Centrifugal Field (RCF) is the ratio of the centrifugal acceleration at a specified radius and speed ($r\omega^2$) to the standard acceleration of gravity (g) according to the following formula:

$$RCF = \frac{r\omega^2}{g}$$

where r is the radius in millimeters, ω is the angular velocity in radians per second ($2\pi \text{ RPM} / 60$), and g is the standard acceleration of gravity (9807 mm/s^2). After substitution:

$$RCF = 1.12 r \left(\frac{\text{RPM}}{1000} \right)^2$$