

5. Specifications

- **The Boekel Scientific Tube Rotator provides:**

- Digital Microprocessor Control.
- Adjustable speed provides gentle or vigorous mixing action.
- Adjustable angled head offers mixing angles from 90 to 180 degrees
- Five Carousel options provide greater versatility for most tube mixing needs
- Timer feature offers time-down or time-up counting feature
- 2 Year warranty

- **Speed range**4 to 80 rpm

- **Carousel angle**90-180 degrees

- **Timer**Up to 99 Hrs and 59 minutes

- **Carousel Capacity**

Re-order Number

14 X 27 – 34 mm tubes (50 ml Tubes)	260653
46 X 9 -11 mm tubes (1.5/2.0 ml Tubes)	260651
28 X 14 – 18 mm tubes (15 ml Tubes)	260652
63 X 14 mm glass culture tubes (drum rotor)	260655
63 X 18.5 mm glass culture tubes (drum rotor)	260654

- **Power**..... 24VDC, 2.70A via 100-240 volt AC adaptor