



**GS-3**  
 Aluminum  
 6-Place, 500 ml  
 Bottle Dimensions: 70 × 190 mm  
 Fixed-Angle: 20°  
 Maximum Speed: 9,000 rpm  
 K Factor: 4203  
 Rotor Mass: 21.2 kg  
 Autoclavable: 121°C

	RCF (g)	R (cm)
Max.	13,689	15.13
Min.	3,565	3.94

		RC-28S	RC-26 Plus	RC-24 <sup>1</sup>	RC-5C Plus RC-5C <sup>1</sup>	RC-5B Plus RC-5B <sup>1</sup> RC-5 <sup>1</sup>	RC-2B <sup>1</sup> RC-2 <sup>1</sup>	SS-3 <sup>1</sup> SS-4 <sup>1</sup>	RC-3B <sup>1</sup> RC-3B Plus RC-3C <sup>1</sup> RC-3C Plus
<b>GS-3</b>	<b>Max Speed (rpm)</b>	9,000	9,000	9,000	9,000	9,000	9,000	n/a	n/a
	<b>Max RCF (g)</b>	13,689	13,689	13,689	13,689	13,689	13,689	n/a	n/a

**Catalog No. 07017**

Includes:  
 GS-3 Rotor  
 Rotor Cover **07025**  
 Rotor Extractor Tool **01047**  
 Instruction Manual **28505**

**B-2**

**TUBES AND BOTTLES**

CAPACITY (ml)		DESCRIPTION	MAX SPEED (rpm)	DIMENSIONS (mm)	TUBE CAT. NO. (Qty/Pkg)	SEALING ASSEMBLY*		ADAPTER REQ'D		NOTES
NOMINAL	FILL					CAT. NO. (Qty/Pkg)	DESCRIPTION	CAT. NO. (Qty/Pkg)	NO. OF PLACES	
500	480	Polypropylene DRY-SPIN	9,000	70 × 190	<b>06829</b> (6)	Included with Bottle	Polypropylene, Sealing	—	—	Replacement Cap <b>06841</b> (12/pkg). Replacement O-ring sealing assembly <b>03545</b> (24/pkg).
500	425	Polypropylene Oak Ridge Bottle	8,000	70 × 160	<b>03943</b> (6)		Polypropylene Sealing	—	—	Nonsealing Cap <b>03928</b> (12/pkg) available.
500	425	Polycarbonate Oak Ridge Bottle	9,000	70 × 160	<b>03944</b> (6)		Polypropylene Sealing	—	—	Nonsealing Cap <b>03928</b> (12/pkg) available.
250	280	Polypropylene DRY-SPIN	9,000	61 × 46	<b>03069</b> (12)	Included with Bottle	Polypropylene, Sealing	<b>00516</b> (1) <b>00614</b> (1)	1 1	Replacement Cap <b>06841</b> (12/pkg). Replacement O-ring sealing assembly <b>03545</b> (24/pkg). Adapters <b>00516</b> and <b>00614</b> must be used together.
250	270	Polycarbonate DRY-SPIN	9,000	61 × 146	<b>03098</b> (12)	Included with Bottle	Polypropylene, Sealing	<b>00516</b> (1) <b>00614</b> (1)	1 1	Replacement Cap <b>06841</b> (12/pkg). Replacement O-ring sealing assembly <b>03545</b> (24/pkg). Adapters <b>00516</b> and <b>00614</b> must be used together.
35	—	Polyallomer Nalgene® Conical Oak Ridge Tube	—	29 × 105	—	—	—	<b>03072</b> (1) <b>00614</b> (1)	1 1	Adapters <b>03072</b> and <b>00614</b> must be used together.
35	—	Polycarbonate Nalgene® Conical Oak Ridge Tube	—	29 × 105	—	—	—	<b>03072</b> (1) <b>00614</b> (1)	1 1	Adapters <b>03072</b> and <b>00614</b> must be used together.
315-3	**	**	**	**	**	**	**	—	—	Additional Adapter <b>00614</b> required.

<sup>(1)</sup>Instrument no longer available.

\*Sealing Assembly replacement parts may be available. See section E for availability.

\*\*All tubes and adapters listed under GSA rotor may be used in GS-3 with additional adapter **00614**. See pages B-39 and B-40 of this catalog.

All trademarks are the property of their respective owners.