



## VWR® General Purpose Water Baths with Digital Controllers

The latest generation of VWR International, LLC offers superior temperature control, range, and uniformity. Bath fluids can be controlled at temperatures as high as 100°C (60°C without cover) with 0.1°C precision and  $\pm 0.2^\circ\text{C}$  uniformity. Bath temperature is displayed continuously on a bright, easy-to-read LED panel in your choice of  $^\circ\text{C}$  or  $^\circ\text{F}$ . Set point temperature is recalled with just the touch of a button. Three user-defined temperature preset buttons allow for quick selection of frequently used temperature set points.

Dual thermostats provide optimum protection for your work and water bath. The high limit alarm alerts you if bath temperature exceeds your pre-set limit. A secondary Safety Set thermostat guards against thermal runaway, automatically disconnecting heater power should bath temperature get too high or the liquid level drop too low.

VWR Water Baths are also designed for operating convenience. The steeply gabled, polycarbonate cover accommodates glassware of varying heights and tilts out of your way when loading or removing samples, allowing condensate to drain neatly back into the bath.



- Temperature Uniformity:  $\pm 0.2^\circ\text{C}/\pm 1.8^\circ\text{F}$
- Temperature Stability:  $\pm 0.25^\circ\text{C}$
- Display Resolution: 0.1°



### Digital Controller

- Accurate Microprocessor Control
- Three User-defined Temperature Preset Buttons
- Redundant Safety Backup
- Front Panel Lockout
- Electronic Calibration
- 0.1° Resolution

Capacity	2 liter	2 liter shallow	5 liter	10 liter
<b>Digital Water Baths Specifications</b>				
Temperature Range	Ambient +5°C to 100°C			
Working Access (L x W x D)	3.9 x 4.3 x 5.5 in 9.9 x 10.9 x 14 cm	5 x 10.8 x 1.9 in 12.7 x 27.4 x 4.8 cm	5 x 10.8 x 5.5 in 12.7 x 27.4 x 14 cm	10.6 x 11.6 x 5.5 in 26.9 x 29.5 x 14 cm
Overall Dimensions (L x W x H)	10.2 x 8 x 11.8 in 25.9 x 20.3 x 30 cm	10.6 x 13.6 x 12 in 26.9 x 34.5 x 30.4 cm	10.6 x 13.6 x 12 in 26.9 x 34.5 x 30.4 cm	16.7 x 14.8 x 11.8 in 42.4 x 37.6 x 30 cm
Heater (60Hz & 50Hz)	300W	300W	500W	500W
Cat. No. (120VAC/60Hz)	97025-126	97025-138	97025-130	97025-134



Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit [www.vwr.com](http://www.vwr.com) to view our privacy policy, trademark owners and additional disclaimers. © 2012 VWR International, LLC. All rights reserved.