

VWR™ Hotplate Stirrers - VMS



VWR™ Magnetic Stirrers with Heating

- Seamless ceramic top plate
- 4 x 4, 7 x 7, 10 x 10 models
- Temp range 50 - 500°C
- Wide speed range (100 - 1500 rpm)



82026-776

Digital Display - Allows for accurate and repeatable temperature settings.

Microprocessor Control - Sends feedback to the motor to maintain constant speed.

Temperature Controller - Connection for VT-4 contact thermometer enables precise temperature control of the medium.

Ceramic Top - Seamless one-piece ceramic top provides excellent resistance to chemicals.

Elevated Control Panel - Minimizes contact with spills and allows for easy reading and adjustment of controls.

Powerful Motor - Provides excellent slow speed control.

3-Year Warranty - Limited parts and labor warranty from time of purchase.

SPECIFICATIONS	VMS - C4	VMS - C7	VMS - C10
VWR Cat #	82026-772	82026-776	82026-780
Voltage	115 v	115 v	115 v
Stirring Quantity (H ₂ O)	5 l	10 l	15 l
Speed Display	scale	scale	scale
Speed Range	100 - 1500 rpm	100 - 1500 rpm	100 - 1500 rpm
Heat Power	250 W	1000 W	1500 W
Temperature Range	50 - 500°C	50 - 500°C	50 - 500°C
Connection for Temp Controller	no	yes	yes
Fixed Safety Circuit	550°C	550°C	550°C
Material Heating Plate	glass-ceramic	glass-ceramic	glass-ceramic
Dimensions Heating Plate	4 x 4 in	7 x 7 in	10 x 10 in
Dimensions (W x D x H)	150 x 260 x 105 mm	220 x 335 x 105 mm	300 x 415 x 105 mm
Weight	approx. 6.5 lb	approx. 11 lb	approx. 13 lb
Certifications	UL / CUL / CE	UL / CUL / CE	UL / CUL / CE
LIST PRICE	\$338.00	\$358.00	\$478.00