

---

# Performance Characteristics and Specifications

<b>Model</b>	<b>M48725</b>	<b>M48720-33 &amp; M48720-33CN</b>
Volts	120	240
Amps	.05	.02
Hz	50/60	50
Watts	3	3
Design	Designed for continuous service.	Designed for continuous service.
Tube Diameter	6 mm to 20 mm	6 mm to 20 mm
Tube Length	Up to 150 mm	Up to 150 mm
Tube Capacity	16	16
Pad Dimensions	4-1/2" x 12"	4-1/2" x 12"
Ambient Operating Temperature	-4°C - 60°C	-4°C - 60°C
Ambient Relative Humidity	80%	80%
Maximum Weight Capacity of Load	1-1/4 lbs.	1-1/4 lbs.