






Multichannel Models

| VistaLab Catalog # | | Color | Volume Range | Dispensing Increments | Accuracy | Precision |
|--------------------------|--------------------------|---|--------------|-----------------------|--|--|
| 8 Channel | 12 Channel | | | | | |
| 1160-0020 1160-0020L* | 1060-0020 1060-0020L* |  | 0.5 – 20µL | 0.05µL | ±1.0% at 20µL ±1.5% at 10µL ±3.0% at 2µL | 0.40% at 20µL 0.80% at 10µL 1.80% at 2µL |
| 1160-0125 1160-0125L* | 1060-0125 1060-0125L* |  | 2 – 125µL | 0.5µL | ±0.8% at 125µL ±1.0% at 62.5µL ±3.0% at 12.5µL | 0.25% at 125µL 0.30% at 62.5µL 1.25% at 12.5µL |
| 1160-0250 1160-0250L* | 1060-0250 1060-0250L* |  | 5 – 250µL | 1µL | ±0.8% at 250µL ±1.0% at 125µL ±2.5% at 25µL | 0.25% at 250µL 0.40% at 125µL 1.00% at 25µL |
| | 1060-0850 1060-0850L* |  | 25 – 850µL | 5µL | ±1.0% at 850µL ±1.25% at 425µL ±3.5% at 85µL | 0.30% at 850µL 0.40% at 425µL 1.00% at 85µL |
| 1160-1250 1160-1250L* | |  | 25 – 1250µL | 5µL | ±0.8% at 1250µL ±1.0% at 625µL ±3.5% at 125µL | 0.20% at 1250µL 0.30% at 625µL 0.80% at 125µL |

