







VistaLab Product No.	Color	Volume Range	Dispensing Increments	Accuracy	Precision
1057-0002		0.2 – 2 μ L	0.002 μ L	\pm 1.5% at 2 μ L \pm 2.7% at 1 μ L \pm 9% at 0.2 μ L	0.9% at 2 μ L 1.8% at 1 μ L 7.6% at 0.2 μ L
1057-0010		1 – 10 μ L	0.01 μ L	\pm 1% at 10 μ L \pm 1.5% at 5 μ L \pm 2.5% at 1 μ L	0.4% at 10 μ L 0.6% at 5 μ L 1.2% at 1 μ L
1057-0020		2 – 20 μ L	0.02 μ L	\pm 1% at 20 μ L \pm 1.5% at 10 μ L \pm 7% at 2 μ L	0.3% at 20 μ L 0.5% at 10 μ L 2% at 2 μ L
1057-0100		10 – 100 μ L	0.1 μ L	\pm 0.8% at 100 μ L \pm 0.8% at 50 μ L \pm 3.5% at 10 μ L	0.15% at 100 μ L 0.24% at 50 μ L 1% at 10 μ L
1057-0200		20 – 200 μ L	0.2 μ L	\pm 0.8% at 200 μ L \pm 0.8% at 100 μ L \pm 2.5% at 20 μ L	0.15% at 200 μ L 0.25% at 100 μ L 1% at 20 μ L
1057-1000		100 – 1000 μ L	1.0 μ L	\pm 0.8% at 1000 μ L \pm 0.8% at 500 μ L \pm 3% at 100 μ L	0.15% at 1000 μ L 0.2% at 500 μ L 0.6% at 100 μ L
1057-0999	Ovation Quick-Set 3-piece Set (Includes 2-20 μ L, 20-200 μ L, 100-1000 μ L, 3 shelf hangers and 3 racks of tips.)				

