

PARTS LIST 1324 & 1325F

Description	115V	220V
Main Temperature Control	100026	100026
Over Temperature Control	1750615	1750615
Cord Set	1800516	101990
Door Handle	300142	300142
Heating Element 1324	9570771	9570801
Heating Element 1325F	9570777	9570804
On/Off Switch	103351	103351
Knob Main Temperature	X1000771	X1000771
Knob Over Temperature	4450506	4450506
Pilot Light Green	200021	200021
Pilot Light Red	200020	200020
Shelf Clips	200116	200116
Shelf	5500629	5500629
Thermometer Clip	5080865	5080865
Thermometer	200063	200063
Fan Motor 1325F	4880527	4880528
Fan Blade 3" 1325F	104073	104073
Fan Blade 4.75" 1325F	2600502	2600502

UNIT SPECIFICATIONS

Weight	Shipping	Net
1324	lbs.	lbs.
1325F	125 lbs.	85 lbs.

Dimensions	Exterior WxDxH	Interior WxDxH
1324	20 X 24 X 31"	16 X 19 X 22"
1325F	20.5 X 25.25 X 31.5"	16 X 17.5 X 22"

Capacity	Cubic Feet
1324	3.9
1325F	3.6

Temperature	Range	Uniformity
1324	Amb+5° to 200°C	+4°C
1325F	Amb+5° to 200°C	+4°C