Physical Specifications

Table 2 Physical Specifications (1100 Series modules)

Туре	Specification	Comments
Weight	11.0 kg (25.0 lbs) 11.0 kg (25.0 lbs) 10.2 kg (22.5 lbs) 17.0 kg (38.0 lbs)	G1310A ISO Pump G1314A VWD G1316A TCC G1362A RID
Dimensions (width × depth × height)	$345 \times 435 \times 140$ mm (13.5 × 17 × 5.5 inches) $345 \times 435 \times 180$ mm	G1310A, G1314A, G1321A G1362A
	(13.5 × 17 × 7.0 inches) 410 × 435 × 140 mm (16.1 × 17 × 5.5 inches)	G1316A
Line Voltage	100 – 120 or 220 – 240 VAC, ± 10 %	Wide-ranging capability
Line frequency	50 or 60 Hz, ± 5 %	
Power consumption	320 VA maximum 160 VA maximum	G1362A
Ambient operating temperature	0 – 55 °C (32 – 131 °F)	
Ambient non-operating temperature	-40 - 70 °C (-4 - 158 °F)	
Humidity	<95 %, at 25 – 40 °C (77 – 104 °F)	Non-condensing
Operating altitude	Up to 2000 m (6500 ft.)	
Non-operating altitude	Up to 4600 m (14950 ft.)	For storing the instrument
Safety standards: IEC, CSA, UL, EN	Installation Category II, Pollution Degree 2	

For your PC specifications refer to its documentation.

Table 1 Phy	vsical Specifications	
Туре	Specification	Comments
Weight	11 kg (25 lbs)	
Dimensions (height × weight × depth)	140 × 345 × 435 mm (5.5 × 13.5 × 17 inches)	
Line voltage	100–120 or 220–240 VAC, \pm 10%	Wide-ranging capability
Line frequency	50 or 60 Hz, \pm 5%	
Power consumption	220 VA	Maximum
Ambient operating temperatur	re 4–55 °C (41–131 °F)	
Ambient non-operating temper	rature -40–70 °C (-4–158 °F)	
Humidity	< 95%, at 25–40 °C (77–104 °F)	Non-condensing
Operating Altitude	Up to 2000 m (6500 ft)	
Non-operating altitude	Up to 4600 m (14950 ft)	For storing the quaternary pump
Safety standards: IEC, CSA, UL	Installation Category II, Pollution Degree 2	

Introduction to the Quaternary Pump Setting the 8-bit Configuration Switch

In the non-volatile memory the parameters are kept, regardless of whether you turn the instrument off and on again. They will be kept until the same set of parameters is subsequently changed and power is reset. All other previously stored configuration settings will still remain in the non-volatile memory.

In this way you can store more than one set of parameters using the same 8-bit configuration switch twice,, for example, for both GPIB and RS-232C.

GPIB Default Addresses

If you just want to change the GPIB address and need a detailed procedure, refer to the *Installing Your ChemStation System* handbook. Default GPIB address is set to the following addresses::

Default Addresses for Agilent Series 1100 Modules

Module	Address	Binary Address
Pump	22	00010110
FLD	23	00010111
VWD	24	0 0 0 1 1 0 0 0
Agilent 8453A	25	00011101
DAD/MWD	26	00011010
Column compartment	27	00011011
Autosampler	28	00011100
RID	29	00011101

where 0 means that the switch is down and 1 means that the switch is up.

Communication Settings for RS-232C Communication

The communication protocol used in the instrument supports only hardware handshake (CTS/RTR).

Table 36

Performance Specifications

Table 44

Performance Specification Agilent 1100 Series Quaternary Pump

Туре	Specification
Hydraulic system	Dual plunger in series pump with proprietary servo-controlled variable stroke drive, floating plungers and active inlet valve
Setable flow range	0.001 – 10 ml/min, in 0.001 ml/min increments
Flow range	0.2 – 10.0 ml/min
Flow precision	< 0.3 % RSD (typically 0.15 %), based on retention time, at 1 ml/min
Pressure	Operating range 0 – 40 MPa (0 – 400 bar, 0 – 5880 psi) up to 5 ml/min
	Operating range 0 – 20 MPa (0 – 200 bar, 0 – 2950 psi) up to 10 ml/min
Pressure pulsation	< 2 %amplitude (typically < 1 %), at 1 ml/min isopropanol, at all pressures > 1 MPa (10bar)
Compressibility compensation	User-selectable, based on mobile phase compressibility
Recommended pH range	1.0-12.5,solvents with pH >2.3 should not contain acids which attack stainless steel
Gradient formation	Low pressure quaternary mixing/gradient capability using proprietary high-speed proportioning valve Delay volume 800 — 1100 µI, dependent on back pressure
Composition Range	0 – 95 % or 5 – 100 %, user selectable
Composition Precision	< 0.2 % SD, at 0.2 and 1 ml/min
Control and data evaluation	Agilent ChemStation for LC

Specifications Performance Specifications

Table 44

Performance Specification Agilent 1100 Series Quaternary Pump

Analog output	For pressure monitoring, 2 mV/bar, one output	
Communications	Controller-area network (CAN), GPIB, RS-232C, APG Remote: ready, start, stop and shut-down signals, LAN optional	
Safety and maintenance	Extensive diagnostics, error detection and display (through control module and Agilent ChemStation), leak detection, safe leak handling, leak output signal for shutdown of pumping system. Low voltages in major maintenance areas.	
GLP features	Early maintenance feedback (EMF) for continuous tracking of instrument usage in terms of seal wear and volume of pumped mobile phase with user-settable limits and feedback messages. Electronic records of maintenance and errors.	
Housing	All materials recyclable.	

Installing the Autosampler

Site Requirements

Environment

Your autosampler will work within specifications at ambient temperatures and relative humidity as described in Table 1.

CAUTIONDo not store, ship or use your autosampler under conditions where
temperature fluctuations may cause condensation within the autosampler.
Condensation will damage the system electronics. If your autosampler was
shipped in cold weather, leave it in its box, and allow it to warm up slowly
to room temperature to avoid condensation.

Table 1 Physical Specifications - Autosampler (G1313A / G1329A / G1389A / G2260A)

Туре	Specification	Comments
Weight	14.2 kg (31.3 lbs)	
Dimensions (height \times width \times depth)	200 × 345 × 435 mm (8 × 13.5 × 17 inches)	
Line voltage	100 – 120 or 220 – 240 VAC, ± 10 %	Wide-ranging capability
Line frequency	50 or 60 Hz, ± 5 %	
Power consumption (apparent power) Power consumption (active power)	300 VA 200 W	Maximum Maximum
Ambient operating temperature	4-55 °C (41-131 °F)	see WARNING on page 16
Ambient non-operating temperature	-40 to 70 °C (-4 to 158 °F)	
Humidity	<95 %, at 25 $-$ 40 °C (77 $-$ 104 °F)	Non-condensing
Operating Altitude	Up to 2000 m (6500 ft)	
Non-operating altitude	Up to 4600 m (14950 ft)	For storing the autosampler
Safety standards: IEC, CSA, UL	Installation Category II, Pollution Degree 2	

WARNING

Using the autosampler at environmental temperatures higher than 50 $^\circ C$ (122 $^\circ F)$ may cause the rear panel to become hot.

Performance Specifications

Table 62

Performance Specifications Agilent 1100 Series Autosampler (G1313A) and Thermostatted Autosampler (G1329A). Valid when standard 100 μI metering head installed.

Туре	Specification
GLP features	Early maintenance feedback (EMF), electronic records of maintenance and errors
Communications	Controller-area network (CAN). GPIB (IEEE-448), RS232C, APG-remote standard, optional four external contact closures and BCD vial number output
Safety features	Leak detection and safe leak handling, low voltages in maintenance areas, error detection and display
Injection range	0.1 – 100 μl in 0.1 μl increments Up to 1500 μl with multiple draw (hardware modification required)
Replicate injections	1 – 99 from one vial
Precision	Typically < 0.5 % RSD of peak areas from 5 – 100 μl, Typically < 1 % RSD of peak areas from 1 – 5 μl
Minimum sample volume	1 μl from 5 μl sample in 100 μl microvial, or 1 μl from 10 μl sample in 300 μl microvial
Carryover	Typically < 0.1 %, < 0.05 % with external needle cleaning
Sample viscosity range	0.2 – 50 cp
Replicate injections per vial	1 - 99
Sample capacity	100×2 -ml vials in 1 tray
	40×2 -ml vials in $\frac{1}{2}$ tray
	15×6 -ml vials in ½ tray (Agilent vials only)
Injection cycle time	Typically 50 s depending on draw speed and injection volume

Specifications Performance Specifications

Table 63

Performance Specifications Agilent 1100 Series Autosampler (G1313A) and Thermostatted Autosampler (G1329A). Valid when standard 900 µl metering head installed.

Туре	Specification
Pressure	Operating range 0-20 MPa (0-200 bar, 0-2950 psi)
GLP features	Early maintenance feedback (EMF), electronic records of maintenance and errors
Communications	Controller-area network (CAN). GPIB (IEEE-448), RS232C, APG-remote standard, optional four external contact closures and BCD vial number output
Safety features	Leak detection and safe leak handling, low voltages in maintenance areas, error detection and display
Injection range	0.1 – 900 μl in 0.1 μl increments (recommended 1 μl increments) Up to 1800 μl with multiple draw (hardware modification required)
Replicate injections	1 – 99 from one vial
Precision	Typically < 0.5 % RSD of peak areas from 5 – 2000 μl, Typically < 1 % RSD of peak areas from 2000 – 5000 μl, Typically < 3 % RSD of peak areas from 1 – 5 μl
Minimum sample volume	1 μl from 5 μl sample in 100 μl microvial, or 1 μl from 10 μl sample in 300 μl microvial
Carryover	Typically < 0.1 %, < 0.05 % with external needle cleaning
Sample viscosity range	0.2 – 50 cp
Sample capacity	100 × 2-ml vials in 1 tray
	40 × 2-ml vials in $\frac{1}{2}$ tray
	15 × 6-ml vials in $\frac{1}{2}$ tray (Agilent vials only)
Injection cycle time	Typically 50 s, depending on draw speed and injection volume

Table 64

Performance Specifications Agilent 1100 Series Thermostatted Micro Autosampler (G1389A)

Туре	Specification
Sample capacity	$100 \ x \ 2$ ml vials in 1 tray. Microvials (100 or 300 μ l) with sleeves (reduced cooling performance with microvials)
Settable injection	0.01 to 8 µl with small loop capillary
volume	0.01 to 40 μl with extended loop capillary
Precision	Typically < 0.5 % RSD from 5 to 40 μl
	Typically < 1 % RSD from 1 to 5 μ l
	Typically < 3 % RSD from 0.2 to 1 μ l
Minimum sample volume	1 μl from 5 μl sample in 100 μl microvial, or 1 μl from 10 μl sample in 300 μl microvials
Carryover	Typically < 0.1 % without automated needle wash
	Typically < 0.05 % with external needle cleaning and 1 μl injection volume
Sample viscosity range	0.2 - 5 ср
Recommended pH-range	1.0 - 8.5 solvent with pH < 2.3 should not contain acids which attack stainless steel. Upper pH range is limited by fused silica capillaries
Material in contact with solvent	Stainless steel, sapphire, PTFE, PEEK, fused silica, Vespel
GLP features	Early maintenance feedback (EMF), electronic records of maintenance and errors
Communications	Controller-area network (CAN). GPIB (IEEE-448), RS232C, APG-remote standard, optional four external contact closures and BCD vial number output
Safety features	Leak detection and safe leak handling, low voltages in maintenance areas, error detection and display
Housing	All material recyclable

Specifications Performance Specifications

Table 65

Performance Specifications Agilent 1100 Series Preparative Autosampler (G2260A)

Туре	Specification
Pressure	Operating range 0-40 MPa (0-400 bar, 0-5800psi)
GLP features	Early maintenance feedback (EMF), electronic records of maintenance and errors
Communications	Controller-area network (CAN). GPIB (IEEE-448), RS232C, APG-remote standard, optional four external contact closures and BCD vial number output
Safety features	Leak detection and safe leak handling, low voltages in maintenance areas, error detection and display
Injection range	0.1 – 900 μ l in 0.1 μ l increments (recommended 1 μ l increments)
	Up to 1800 µl with multiple draw (hardware modification required)
	Up to 5000 µl with multiple draw (hardware modification required)
Replicate injections	1 – 99 from one vial
Precision	Typically < 0.5 % RSD of peak areas from 5 – 2000 μl, Typically < 1 % RSD of peak areas from 2000 – 5000 μl, Typically < 3 % RSD of peak areas from 1 – 5 μl
Minimum sample volume	1 μl from 5 μl sample in 100 μl microvial, or 1 μl from 10 μl sample in 300 μl microvial
Sample viscosity range	$0.2 - 50 { m ~cp}$
Sample capacity	100 × 2-ml vials in 1 tray
	15 × 6-ml vials in $\frac{1}{2}$ tray (Agilent vials only)
Injection cycle time	Typically 50 s, depending on draw speed and injection volume

Performance Specifications Agilent 1100 Series DAD

Reference conditions for data of Table 92:

- cell path length 10 mm, response time 2 s,
- flow 1 ml/min LC-grade Methanol,
- slit width 4 nm.

Linearity measured with caffeine at 265 nm.

Туре	Specification	Comments
Detection type	1024-element photodiode array	
Light source	Deuterium and tungsten lamps	
Wavelength range	190 – 950 nm	
Short term noise (ASTM) [*] Single and Multi-Wavelength	Typically ± 3 × 10 ⁻⁵ AU at 254 nm at flow rates <100 μl/min	For the 500 nl flow cell the noise is 2-3 times higher than with standard flow cell
Drift	2 × 10 ⁻³ AU/hr at 254 nm	
Linear absorbance range	> 2 AU (upper limit)	
Wavelength accuracy	± 1 nm	Self-calibration with deuterium lines, verification with holmium oxide filter
Wavelength bunching	1 – 400 nm	Programmable in steps of 1 nm
Slit width	1, 2, 4 , 8, 16 nm	Programmable slit
Diode width	< 1 nm	
Flow cell	500 nanoliter: 0.5 μl volume, 10 mm cell path length and 50 bar (725 psi) pressure maximum	

Туре	Specification	Comments
Maximum pressure	50 bar	
Control and data evaluation	Agilent ChemStation for LC	
Analog outputs	Recorder/integrator: 100 mV or 1 V, output range 0.001 – 2 AU, two outputs	
Communications	Controller-area network (CAN), GPIB, RS-232C, APG Remote: ready, start, stop and shut-down signals, LAN optional	
Safety and maintenance	Extensive diagnostics, error detection and display (through control module and ChemStation), leak detection, safe leak handling, leak output signal for shutdown of pumping system. Low voltages in major maintenance areas.	
GLP features	Early maintenance feedback (EMF) for continuous tracking of instrument usage in terms of lamp burn time with user-settable limits and feedback messages. Electronic records of maintenance and errors. Verification of wavelength accuracy with built-in holmium oxide filter.	
Housing	All materials recyclable.	

 Table 92
 Performance Specifications Agilent 1100 Series Diode Array Detector

* ASTM: "Standard Practice for Variable Wavelength Photometric Detectors Used in Liquid Chromatography".

For specification on the 500 nl flow cell refer to Table 85 on page 153.

Performance Specifications Agilent 1100 Thermostatted Column Compartment.

All specifications in Table 91 are valid for distilled water at ambient temperature (25 $^{\circ}$ C), set point at 40 $^{\circ}$ C and a flow range from 0.2–5 ml/min.

For flow rates below 100 $\mu l/min$ the column bracket must be installed

Туре	Specification	
Temperature range	10 degrees below ambient to 80 °C	
Temperature stability	± 0.15 °C	
Column capacity	Three 25 cm - NOTE: With fused silica capillaries connected, length limited by bend radii of capillary	
Warm-up/cool-down time	5 minutes from ambient to 40 °C 10 minutes from 40 – 20 °C	
Internal volume	3 μl left heat exchanger 6 μl right heat exchanger	
Communications	Controller-area network (CAN), GPIB, RS-232C, APG Remote: ready, start, stop and shut-down signals, LAN optional	
Safety and maintenance	Extensive diagnostics, error detection and display (through control module and Agilent ChemStation), leak detection, safe leak handling, leak output signal for shutdown of pumping system. Low voltages in major maintenance areas.	
GLP features	Column-identification module for GLP documentation of column type, see "Column-Identification System"	
Housing	All materials recyclable.	

 Table 91
 Performance Specifications Agilent 1100 Series Thermostatted Column Compartment

7 Performance Specifications

Performance Specifications Agilent 1100 Series Micro Vacuum Degasser

Туре	Specification	
Flow rate	0–5 ml/min per channel (5–10 ml/min at reduced degassing performance)	
Number of channels	4	
Internal volume per channel	Typically 1 ml per channel	
Materials in contact with solvent	PTFE – FEP – PEEK	
oH range	1 – 14	
Analog output (AUX)	For pressure monitoring, range $0-3$ V	
Evaporation of solvents into the atmosphere	< 200 µg/m ³ Acetonitrile and Methanol certification by IAS.	

 Table 88
 Performance Specifications Agilent 1100 Micro Vacuum Degasser