

Technical data

Power supply unit

Input

VAC 100 - 240


A 0,3

Hz 50 - 60

VDC 12

7,2W.LPS

(limited power source)

2 (double insulated) 

Base plate

Diameter base plate

mm 120

Material

glass

Output

Protection class

Color Squid

Operating voltage

VDC 12

mA 300

°C +5 to +40

% 80

% 100

IP54

2

II

Operation at a terrestrial altitude m max. 2000 above sea level

Dimensions mm 160 x 145 x 45

Weight kg 0,55

Motor

Speed range

rpm 0 - 2500

0 - 2000rpm

at 50 rpm steps

2000 - 2500rpm

at 100 rpm steps

Power consumption

W 3

power output

W 2

Stirring quantity max. (H₂O)

l 1