
SPECIFICATIONS

ELECTRICAL

3518: 120 VAC, 50/60 Hz, 100 Watts
3519: 240 VAC, 50/60 Hz, 100 Watts
3520: 120 VAC, 50/60 Hz, 100 Watts
3521: 240 VAC, 50/60 Hz, 100 Watts

SHAKER MOTION

All Models: 3/4-inch circular motion.

SHAKER SPEED

All Models: 40 - 400 rpm.

TIMER

All Models: 1 - 60 minutes or continuous.

DIMENSIONS

11"W x 13-1/2"D x 7"H (28 x 34 x 18 cm)

NET WEIGHT

3518, 3520: 34 lb (15.4 kg)
3519, 3521: 38 lb (17.3 kg)