

Appendix A Specifications

Output specifications:	500 V, 2500 mA, 500 W
Output range (programmable):	10–500 V, fully adjustable in 1 V increments 10–2500 mA, fully adjustable in 1 mA increments 1–500 W, fully adjustable in 1 W increments
Type of output:	Constant voltage, current, or power with automatic crossover
Output terminals	Four pair of recessed banana jacks in parallel
Timing mode	Continuous, timed to 99:59Hrs, or 99999 volt-hours.
Pause/resume function:	Yes
Display:	128 x 64 Monochrome Backlit LCD
Operating Conditions:	0–40°C; 0–95% humidity in the absence of condensation
Safety compliance:	EN61010
EMI:	Conforms to CE Standards for emissions and immunity class A.
Safety features:	No load detection; sudden load change detection; ground leak detection; overload/short circuit detection; overvoltage protection; overheating protection
Input protection:	Fuse on hot and neutral
Input power(nomial)	100–120/220–240 VAC, 50/60 Hz auto-switching
Dimensions:	27.5cm X 9.8cm X 34.5 cm
Weight:	2.5 kg