

UNIT SPECIFICATIONS

These units are 115 volt or 220 volt. Please refer to the unit data plate for its individual specification.

Weight	Shipping	Net
2325	280 lbs.	221 lbs.
2375	600 lbs.	442 lbs.
2425	280 lbs.	221 lbs.
2475	600 lbs.	442 lbs.

Dimensions	Exterior WxDxH (IN.)	Interior WxDxH (IN.)
2325	26.5x26x40	20.5x20x28
2375	26.5x26x80	20.5x20x28(x2)
2425	26.5x26x40	20.5x20x28
2475	26.5x26x80	20.5x20x28 (x2)

* each chamber

Capacity	Cubic Feet
2325	6.7
2375	6.7 (x2)
2425	6.7
2475	6.7 (x2)

Temperature	Range	Uniformity	CO2
2325	Amb. 8°C to 60°C	±0.25°C @ 37°C	0 – 20%
2375	Amb. 8°C to 60°C	±0.25°C @ 37°C	0 – 20%
2425	Amb. 8°C to 60°C	±0.25°C @ 37°C	0 – 20%
2475	Amb. 8°C to 60°C	±0.25°C @ 37°C	0 – 20%

Model	Amp	Voltage	Hertz
2325	5 amp	110-120V	50/60 Hz
2325-2	3.5 amp	208-240V	50/60 Hz
2375	5 amp	110-120V	50/60 Hz
2375-2	3.5 amp	208-240V	50/60 Hz
2425	5 amp	110-120V	50/60 Hz
2425-2	3.5 amp	208-240V	50/60 Hz
2475	5 amp	110-120V	50/60 Hz
2475-2	3.5 amp	208-240V	50/60 Hz