

Product Specifications

Power Source™ 300V Power Supply Specifications

Input Power (switchable)	110 VAC, 50-60 Hz 230 VAC, 50-60 Hz
Fuses	One 4A/250V, one extra fuse is provided
Output power in watts	90 watts
Output voltage range	2~300V
Output current range	4~500 mA
Duration Limits Time	~1min - 99hr 59min
Terminal pairs	4 (4 positive voltage and 4 negative voltage)
Operating Modes	
Constant Voltage	1V increment settings
Constant Current	1 mA increment settings
Crossover	Automatic
Display type	Backlit LCD Graphic type
Display size	53.64 x 15.64 mm (W x H)
Pause function	Yes
Safety features	No Load Detection Load Change Detection Overload Detection Ground Leak Detection Auto Restart
Stackable	Yes
Housing material	Flame retardant ABS
Housing size	200 x 290 x 70 mm (W x D x H)
Operating temp.	0°C-40°C
Environmental condition	85% RH, 75 KPa-106 Kpa, Altitude not to exceed 2000 meters
Weight	1.2 kg
Certifications	CE, RoHS, ETL, CSA
Warranty	2 years