

## Rotavac Valve Tec

P/N: 036304780



- Recommended for solvents with low or medium boiling points
- Vacuum can be controlled via valve operated or manual vacuum controllers
- A Vacuum valve is necessary, if you want to combine this vacuum pump with the Hei-VAP Expert Control / Hei-VAP Ultimate Control rotary evaporators

Condenser is not in scope of delivery.

## Rotovac Valve Tec - Technical Data

Weight	6 kg
Diaphragm pump	Two-stage
Suction capacity	0.75 m <sup>3</sup> /h
Ultimate vacuum	12 mbar
Combinable with a condenser	yes
Power input	80 W
Dimensions without condenser (w/d/h)	145 x 315 x 185 mm