

---

# Performance and Physical Data

Model No. . . . .	285A
Catalog No. . . . .	13-262-285A (120V), 13-262-286A (240V)
Vacuum Range . . . . .	Atmosphere to 30" Hg.* with adequate pump
Vacuum Leak Rate. . . . .	< .2" Hg per 24 hours
Temperature Range . . . . .	Ambient +10°C to 200°C
Shelf Capacity . . . . .	4.1 ft <sup>2</sup>
Inside Dimensions . . . . .	19"D x 12"W x 12"H
Outside Dimensions . . . . .	25"D x 22.5"W x 23.5"H
Standard Electrical Service . . . . .	1700 Watts, 120 VAC, 50/60 Hz, 14.2 Amps
Optional: . . . . .	1700 Watts, 240 VAC, 50/60 HZ, 7.1 Amps
Net Weight (Pounds) . . . . .	130 lb.

\* *within 10 Microns*