



Overview

PTFE-coated collector coil and chamber.

Catalog Number: 7934021

Specifications

- **Weight:** 270.0 lbs
- **Weight metric:** 122.0 kg
- **Dimensions:** 32.2" w x 28.0" d x 36.3" h
- **Dimensions metric:** 81.8 cm w x 71.1 cm d x 92.2 cm h
- **Electrical:** 115 volts, 60 Hz, 16.0 A
- **Region:** U.S. and Canada
- **Collector Temperature:** -84Â° C, -119Â° F
- **Conformance:** CAN/CSA C22.2, UL
- **Ice Holding Capacity:** 6L
- **Options Included:** Corrosion-Resistant Collector
- **Plug Type:** North America, 115V
- **Style:** Console
- **Vacuum Pump Displacement Required:** 98 L/min or larger displacement

For Product Assistance



[Andy Miller](#)
Senior Sales Representative



[Vivian Goldmon](#)
Customer Service Representative

