

## TX-750 Swinging Bucket Aluminum Rotor with Rectangular Buckets



| Rotor Specifications |           |             |
|----------------------|-----------|-------------|
| Capacity (ml)        | 4 x 250   |             |
| Tube Dimensions (mm) | 63 x 130  |             |
| Angle °              | 90        |             |
| Maximum Speed (rpm)  | 4.500     |             |
| K Factor             | 8983      |             |
| Net Weight (kg)      | 6.8       |             |
| RCF/Radius           | RCF (x g) | Radius (cm) |
| Maximum              | 4.415     | 19.5        |
| Minimum              | 2.015     | 8.9         |



## TX-750 Rotor Packages

| Cat. No. | Description  |
|----------|--|
| 75003607 | TX-750 Rotor Body  |
| 75003614 | Set of 4 250 ml Buckets without Lids                     |
| 75003615 | Set of 4 ClickSeal Lids                                  |
| 75003616 | Replacement O-rings for 75003615 Sealing Caps (set of 4) |

## Accessories

| Cat. No. | Description  |
|----------|--|
| 75003617 | Set of 4 Non-Sealed Microplate Carriers for Microplates or T-75 Nunc Easy Flasks |

## Tube Ordering Information

| Cat. No.  | Tube Vol. (ml) | Vol. Fill (ml) | Description                   | Qty. | Max. Speed (rpm) | Dims ØxL (mm) | Sealing Assembly Required: |          |             | Adapters Required: |          |               | Accessories Required: |          |               |
|-----------|----------------|----------------|-------------------------------|------|------------------|---------------|----------------------------|----------|-------------|--------------------|----------|---------------|-----------------------|----------|---------------|
|           |                |                |                               |      |                  |               | Cat. No.                   | Qty /Pkg | Description | Cat. No.           | Qty /Pkg | No. of Places | Cat. No.              | Qty /Pkg | No. of Places |
| 3141-0250 | 250            | 250            | PP Oak Ridge Bottle           | 4    | 4.700            | 63x130        | Included                   | 4        | PP Sealing  | 75003737           | 4        | 1             | -                     | -        |               |
| 3140-0250 | 250            | 250            | PC Oak Ridge Bottle           | 4    | 4.700            | 63x130        | Included                   | 4        | PP Sealing  | 75003737           | 4        | 1             | -                     | -        |               |
| -         | 150            | -              | Round Bottom Open Topped Tube | -    | -                | 56x124        | -                          | -        | -           | 75003738           | 4        | 1             | -                     | -        |               |
| 76009007  | 100            | 75             | Glass Bottle                  | 10   | 4.700            | 44x100        | -                          | -        | -           | 75003742           | 4        | 1             | -                     | -        |               |
| 76009084  | 100            | 75             | PP Bottle                     | 1    | 4.700            | 45x98         | -                          | -        | -           | 75003742           | 4        | 1             | -                     | -        |               |
| 76009095  | 100            | 75             | PC Bottle                     | 1    | 4.700            | 45x98         | -                          | -        | -           | 75003742           | 4        | 1             | -                     | -        |               |
| 334959    | 50             | 50             | Nunc Conical                  | 25   | 4.700            | 30x121        | Included                   | 25       | PP Sealing  | 75003685           | 4        | 5             | -                     | -        |               |
| -         | 50             | -              | Round Bottom Open Topped Tube | -    | -                | 35x121        | -                          | -        | -           | 75003749           | 4        | 3             | -                     | -        |               |
| 3139-0050 | 50             | 50             | PP Oak Ridge Tube             | 10   | 4.700            | 29x114        | Included                   | 10       | PP Sealing  | 75003750           | 4        | 4             | -                     | -        |               |
| 3138-0050 | 50             | 50             | PC Oak Ridge Tube             | 10   | 4.700            | 29x114        | Included                   | 10       | PP Sealing  | 75003750           | 4        | 4             | -                     | -        |               |
| 3114-0050 | 50             | 50             | Teflon Oak Ridge Tube         | 2    | 4.700            | 29x114        | Included                   | 2        | PP Sealing  | 75003750           | 4        | 4             | -                     | -        |               |
| -         | 45             | -              | Flat or Round Bottom Tube     | -    | -                | 30x121        | -                          | -        | -           | 75003750           | 4        | 4             | -                     | -        |               |
| 45500-30  | 30             | 25             | KIMAX Glass Tube              | 6    | 4.700            | 24x106        | -                          | -        | -           | 75003756           | 4        | 6             | -                     | -        |               |
| -         | 25             | -              | Universal                     | -    | -                | 26x121        | -                          | -        | -           | 75003755           | 4        | 4             | -                     | -        |               |