

# 12-Port Drying Chamber with PTFE



View online: <https://www.labconco.com/product/12-port-drying-chamber-ptfe/6573>

## Overview



Drying Chambers provide a means to attach flasks to your freeze dryer or contain samples in bulk.

The 12-Port Drying Chamber with PTFE is made of type 304 stainless steel with a corrosion-resistant PTFE-coated interior. Includes 12 neoprene valves with plastic knobs that accommodate both 1/2" and 3/4" adapters for connection of flasks, a 1/2" thick acrylic lid with neoprene gasket, and one latching pin. It is compatible with FreeZone 1, 2.5, 4.5, 6, 8, 12 or 18 Liter Freeze Dryers.

Catalog Number: 7522801

## Specifications

- **Weight:** 11.0 lbs
- **Weight metric:** 5.0 kg
- **Dimensions:** 8.0" diameter x 9.75" h
- **Dimensions metric:** 20.3 x 24.8 cm
- **Applicable Products:** Freeze Dryers

## For Product Assistance



[Andy Miller](#)  
Senior Sales Representative



[Vivian Goldmon](#)  
Customer Service Representative



Labconco Corporation  
8811 Prospect Avenue  
Kansas City, MO 64132-2696  
Phone: (816) 333-8811

