

Thermo Scientific General Purpose Rotors

BIOLiner High Capacity Rotor

BIOLiner High Capacity Rotor



Rotor Specifications	
Capacity (ml)	4 x 750
Tube Dimensions (mm)	98x124
Angle °	90
Maximum Speed (rpm)	3.500
K Factor	13.307
Net Weight (kg)	7.3

RCF/Radius	RCF (x g)	Radius (cm)
Maximum	2739	20
Minimum	1436	10.5



BIOLiner Rotor Packages

Cat. No.	Description
75003667	BIOLiner High Capacity Rotor Body
75003668	Set of 4 BIOLiner Buckets without sealing lids
75003669	Set of 4 ClickSeal Lids for BIOLiner Buckets
75003688	Replacement O-rings for BIOLiner Sealing Lids 75003669 (Set of 4)
75003670	Set of 4 Non-Sealed Microplate Carriers for BIOLiner Rotor
75003689	Set of 4 Sealing Lids for Microplate Carriers for BIOLiner Rotor
75003626	Replacement O-rings for BIOLiner Microplate Lids 75003689 (Set of 4)

Tube Ordering Information

Cat. No.	Tube Vol. (ml)	Vol. Fill (ml)	Description	Qty.	Max. Speed (rpm)	Dims ØxL (mm)	Sealing Assembly Required:			Adapters Required:		Accessories Required:		
							Cat. No.	Qty /Pkg	Description	Cat. No.	Qty /Pkg	No. of Places		
75003699	750	750	PP Bio-Bottle	12	3.500	98x135	Included	12	PP Sealing	75003692	4	1	-	-
-	500	500	Corning Conical	-	3.500	98x124	-	-	-	75003692+Corning Cushion#431124	4	1	-	-
-	250	250	Corning Conical	-	3.500	62x145	-	-	-	75003691	4	1	-	-
3141-0250	250	250	PP Oak Ridge Bottle	4	3.500	62x130	Included	4	PP Sealing	75003675	4	2	-	-
3140-0250	250	250	PC OakRidge Bottle	4	3.500	62x130	Included	4	PP Sealing	75003675	4	2	-	-
-	225	225	BD Conical	-	3.500	62x130	-	-	-	75003675+BD Cusion# 352090	4	2	-	-
376813	200	200	Nunc Conical	4	3.500	62x130	Included	4	PP Sealing	75003675+ Nunc Cushion #377585	4	2	-	-
-	175	175	BD Conical	-	3.500	62x130	-	-	-	75003675+BD Cusion# 352090	4	2	-	-
334959	50	50	Nunc Conical	25	3.500	30x120	Included	25	PP Sealing	75003674	4	10	*Max Capacity may vary by tube brand	
366036	15	15	Nunc Conical	50	3.500	17.5x121	Included	50	PP Sealing	75003673	4	22	-	-
-	10	-	Blood Collection	-	-	17x110	-	-	-	75003672	4	37	-	-
-	7	-	Blood Collection	-	-	14x110	-	-	-	75003671	4	49	-	-
-	5	-	Blood Collection	-	-	13x75	-	-	-	75003671	4	49	-	-
-	MP	-	Standard Microplates (up to 20mm)	-	-	Max Height 60mm	-	-	-	75003670	4	4	-	-
-	MP	-	Deepwell Plates (up to 30mm)	-	-	Max Height 60mm	-	-	-	75003670	4	2	-	-