

Technical Specification

	00077-00	00077-02
Speed	Fixed - 20 rpm	Variable - 2 to 40 rpm
Speed display	None	Digital
Timer	No	Yes – up to 9999 mins
Dimensions (w x d x h)	200 x 270 x 240 mm (7.6 x 10.6 x 9.5 in.)	200 x 270 x 240 mm (7.6 x 10.6 x 9.5 in.)
Weight	3.2 Kg (7 lb)	3.2 Kg (7 lb)
Electrical supply	230V, 50Hz, 50W	230V, 50Hz, 50W

	00077-01	00077-03
Speed	Fixed - 20 rpm	Variable - 2 to 40 rpm
Speed display	None	Digital
Timer	No	Yes – up to 9999 mins
Dimensions (w x d x h)	200 x 270 x 240 mm (7.6 x 10.6 x 9.5 in.)	200 x 270 x 240 mm (7.6 x 10.6 x 9.5 in.)
Weight	3.2 Kg (7 lb)	3.2 Kg (7 lb)
Electrical supply	120V, 60Hz, 50W	120V, 60Hz, 50W