

Cat. No.	Description
6795KIT	<p><b>Corning® Digital Display Hot Plate and Stirring Hot Plate Accessory Kit</b></p> <p>This accessory kit is designed for all Corning digital display hot plates and stirring hot plates (Models PC-400D, PC-420D, PC-600D, PC-620D).</p> <p>This kit includes: external temperature controller (Cat. No. 6795PR), vertical support rod (Cat. No 440129), Supplied as two 9 inch (22.9 cm) support rods that can be screwed together (18 inches), stir bar retriever (Cat. No. 6970SR).</p>
6795-420KIT	<p><b>Corning 5 x 7 in. Top PC-420D Stirring Hot Plate Kit with Digital Displays, 120V/60Hz</b></p> <p>This kit includes: One stirring hot plate (Cat. No. 6795-420D), external temperature controller Cat. No. 6795PR, vertical support rod (Cat. No 440129). Supplied as two 9 inch (22.9 cm) support rods that can be screwed together (18 inches), stir bar retriever (Cat. No. 6970SR).</p>
6795-620KIT	<p><b>6795-620KIT Corning 10 x 10 in. Top PC-620D Stirring Hot Plate Kit with Digital Displays, 120V/60Hz</b></p> <p>This kit includes: One stirring hot plate (Cat. No. 6795-620D), external temperature controller (Cat. No. 6795PR), vertical support rod (Cat. No 440129). Supplied as two 9 inch (22.9 cm) support rods that can be screwed together (18 inches), stir bar retriever (Cat. No. 6970SR).</p>

**Available for purchase separately and not included in the kits**

440140	Boss head clamp. Use for connecting the holding rod (Cat. No. 440141) at a 90° angle to the support rod (Cat. No. 440129). Aluminum.
440141	Holding rod. Use for holding the external temperature controller (Cat. No. 6795PR) in position. Aluminum.

## 11.0 Technical Specifications

Model No.	Power				Temp. Range <sup>2</sup>	Stir Range (rpm)	Weight
	120V/60Hz (6795 Cat. Nos.)	100V/60Hz (6797 Cat. Nos.)	230V/50Hz (6796 Cat. Nos.)	230V/50Hz (6798 Cat. Nos.) <sup>1</sup>			
PC-400D	665W/5.5A	685W/6.9A	635W/2.7A	635W/2.7A	5°C-550°C (41°F-1,022°F)	—	5.5 lbs. (2.5 kg)
PC-410D	15W/0.625A	15W/0.625A	15W/0.625A	15W/0.625A	—	60-1,150	6.2 lbs. (2.8 kg)
PC-420D	675W/5.6A	695W/6.9A	645W/2.8A	645W/2.8A	5°C-550°C (41°F-1,022°F)	60-1,150	6.2 lbs. (2.8 kg)
PC-600D	1215W/10.13A	1015W/10.15A	1215W/5.28A	1215W/5.28A	5°C-550°C (41°F-1,022°F)	—	10.0 lbs. (4.5 kg)
PC-610D	15W/0.625A	15W/0.625A	15W/0.625A	15W/0.625A	—	60-1,150	11.5 lbs. (5.2 kg)
PC-620D	1230W/10.3A	1030W/10.3A	1230W/5.4A	230W/5.4A	5°C-550°C (41°F-1,022°F)	60-1,150	11.5 lbs. (5.2 kg)

<sup>1</sup>Cat. Nos. beginning with 6798 in 230V are non-European configuration models.

<sup>2</sup>The Temperature Range using the External Temperature Controller is 5°C to 200°C (41°F to 392°F).

Model No.	Top Plate Size	Dimensions
PC-400D/410D/420D	5 x 7 in. (12.7 x 17.8 cm)	10.8 x 7.8 x 4.3 in. (27.5 x 19.7 x 10.9 cm)
PC-600D/610D/620D	10 x 10 in. (25.4 x 25.4 cm)	4.625 x 11 x 15.375 in. (11.75 x 19.7 x 39.05 cm)