

---

## ***Product Description***

<b>Models</b>	<b>Top Style</b>	<b>Top Size</b>	<b>Dimensions</b>
PC-200/210/220	Ceramic	4" x 5"	4.4H x 5.8W x 7.3D (inches) 11.2H x 14.7W x 18.5D (cm)
PC-400/410/420	Ceramic	5" x 7"	4.4H x 7.8W x 9.6D (inches) 11.2H x 19.8W x 24.4D (cm)
PC-600/610/620	Ceramic	10" x 10"	4.8H x 10.5W x 14.4D (inches) 12.2H x 26.7W x 36.6D (cm)

---

## ***Operating Conditions***

Corning PC-200, PC-400 and PC-600 Series Hotware is designed to be safe under the following conditions:

- Indoor use.
- Altitude up to 2000 m.
- Temperature of 0°C to 40°C
- Maximum relative humidity 80% for temperatures up to 31°C decreasing linearly to 50% relative humidity at 40°C.
- Main supply voltage fluctuations not to exceed  $\pm 10\%$  of the nominal voltage.
- Installation Category II.
- Pollution Degree 2.
- Position at least 12" from walls, 48" from ceilings and 12" between hot plates if using multiple units.

---

## **Technical Specifications - 120V**

<b>Model</b>	<b>Catalog Number</b>	<b>Power (Volts/Hz/Watts/Amp)</b>	<b>Temp. Range</b>	<b>Stir Range</b>	<b>Weight</b>
--------------	-----------------------	-----------------------------------	--------------------	-------------------	---------------

---

### **Hot Plates**

---

PG-200	5795-200	120V/60 Hz/253W/2.2A	77-1022°F 25-550°C	-	4.0 lbs. 1.8 kg.
--------	----------	----------------------	-----------------------	---	---------------------

---

PG-400	5795-400	120V/60 Hz/528W/5.3A	77-1022°F 25-550°C	-	6.0 lbs. 2.7 kg.
--------	----------	----------------------	-----------------------	---	---------------------

---

PG-600	5795-600	120V/60 Hz/1043W/8.7A	77-1022°F 25-550°C	-	10.0 lbs. 4.5 kg.
--------	----------	-----------------------	-----------------------	---	----------------------

---

### **Stirrers**

---

PG-210	5795-210	120V/60 Hz/33W/0.3A	-	60-1100	4.5 lbs. 2.0 kg.
--------	----------	---------------------	---	---------	---------------------

---

PG-410	5795-410	120V/60 Hz/73W/0.7A	-	60-1100	6.5 lbs. 2.9 kg.
--------	----------	---------------------	---	---------	---------------------

---

PG-610	5795-610	120V/60 Hz/73W/0.7A	-	60-1100	10.0 lbs. 4.5 kg.
--------	----------	---------------------	---	---------	----------------------

---

### **Stirrer/Hot Plates**

---

PG-220	5795-220	120V/60 Hz/283W/2.4A	77-1022°F 25-550°C	60-1100	5.0 lbs. 2.3 kg.
--------	----------	----------------------	-----------------------	---------	---------------------

---

PG-420	5795-420	120V/60 Hz/598W/5.9A	77-1022°F 25-550°C	60-1100	7.0 lbs. 3.2 kg.
--------	----------	----------------------	-----------------------	---------	---------------------

---

PG-620	5795-620	120V/60 Hz/1113W/9.3A	77-1022°F 25-550°C	60-1100	11.5 lbs. 5.2 kg.
--------	----------	-----------------------	-----------------------	---------	----------------------

---