



designed for scientists

Technical Data

| | |
|--|----------------------|
| Type of movement | overhead |
| Permissible shaking weight (incl. attachment) [kg] | 2 |
| Permissible ON time [%] | 100 |
| Speed min. [rpm] | 0 |
| Speed fix [rpm] | 20 |
| Speed display | none |
| Timer display | none |
| Operating mode | continuous operation |
| Dimensions (W x H x D) [mm] | 460 x 180 x 140 |
| Weight [kg] | 1.6 |
| Permissible ambient temperature [°C] | 4 - 50 |
| Permissible relative humidity [%] | 80 |
| Protection class according to DIN EN 60529 | IP 21 |
| Voltage [V] | 100 - 240 |
| Frequency [Hz] | 50/60 |
| DC Voltage [V=] | 24 |
| Current consumption [mA] | 1000 |