

Technical Specifications

EC-I05 Minicell Power Supply (115V/60Hz or 230V/50Hz)

Output mode: Constant Voltage

Voltage Steps:

Lo Range: 20,30,40,50,60,70,80,90, 100, 110,120, 130

Hi Range: 140,150,160,170,180,190,200,210,220,230,240,250

Maximum Voltage: 250 Volts

Maximum Current: 300 milliAmps @ 150-250 Volts

500 milliAmps below 150 Volts

Maximum Power: 75 Watts

Number of Outputs: 2

Fuse Replacement: 1 Amp slo blow

Temperature Range: 0°C -40°C (non-condensing atmosphere)

Dimensions: 8" W x 6" D x 4.5" H (20.3 cm x 13.3 cm x 11.4 cm)

Weight: (net) 5 lbs (2.3 kg)

Weight: (shipping) 7 lbs. (3.2 kg)