

7 Technical data

Power supply:	110 V AC (Nominal) 230 V AC (Nominal)
Input frequency:	47 – 63 Hz
Fuses:	T 5 A (T) (110 V) T 3.15 A (T) (230 V)
Power consumption:	Approx. 500 W (110 V) Approx. 400 W (230 V)
Internal heater: Heater plate:	Ambient to 220 °C max. Approx. 170 °C max.
Dimensions:	Width: 14 cm Depth: 15 cm Height: 33 cm
Weight:	3.4 kg
Instrument housing:	Sheet metal
Operating conditions:	10 to 35 °C, 10 % to 80 %, non-condensing
Storage conditions:	– 20 to 40 °C, 10 % to 80 %, non-condensing

Technical specifications subject to change!