











Rotors for the Centrifuge 5427 R

Centrifuge 5427 R

English (EN)

10.4 Rotor FA-45-30-11 and F-45-30-11

	Rotor FA-45-30-11 Aerosol-tight fixed-angle rotor for 30 micro test tubes Rotor F-45-30-11 Fixed-angle rotor for 30 micro test tubes	Max. <i>g</i> -force:	20 817 × <i>g</i>
		Max. speed:	14 000 rpm
		Max. load (adapter, tube and contents):	30 × 3.75 g

Tubes	Tube Capacity Tubes per adapter/ rotor	Adapter Order no. (international)	Adapter base Diameter	Max. <i>g</i> -force: Max. speed Centrifugation radius
	Micro test tube 1.5 mL/2.0 mL –/30	–	– Ø 11 mm	20 817 × <i>g</i> 14 000 rpm 9.5 cm
	PCR tube 0.2 mL 1/30	 5425 715.005	Conical Ø 6 mm	16 215 × <i>g</i> 14 000 rpm 7.4 cm
	Micro test tube 0.4 mL 1/30	 5425 717.008	Conical Ø 6 mm	20 817 × <i>g</i> 14 000 rpm 9.5 cm
	Micro test tube 0.5 mL 1/30	 5425 716.001	Open Ø 8 mm	18 407 × <i>g</i> 14 000 rpm 8.4 cm
	Microtainers 0.6 mL 1/30	 5425 716.001	Open Ø 8 mm	20 817 × <i>g</i> 14 000 rpm 9.5 cm