

Instructions for use

Mobile table for large benchtop centrifuges

English (EN)

7.2 Ambient conditions

Ambience	For indoor use only.
Ambient temperature	2 °C – 35 °C
Relative humidity	10 % – 75 %, non-condensing.

8 Transport, storage and disposal

8.1 Transport

- ▶ Only transport the mobile table in its original packing.
- ▶ Use a transport aid for transportation over a longer distance.

	Air temperature	Relative humidity	Atmospheric pressure
General transport	-25 °C – 60 °C	10 % – 75 %	30 kPa – 106 kPa
Air freight	-20 °C – 55 °C	10 % – 75 %	30 kPa – 106 kPa

8.2 Storage

	Air temperature	Relative humidity	Atmospheric pressure
In transport packing	-25 °C – 55 °C	10 % – 75 %	70 kPa – 106 kPa
Without transport packing	-5 °C – 45 °C	10 % – 75 %	70 kPa – 106 kPa

9 Ordering information

Order no. (International)	Description
	Mobile table for large benchtop centrifuges
5921 000.010	low version
5921 000.012	high version

Your local distributor: www.eppendorf.com/contact

Eppendorf AG · Barkhausenweg 1 · 22339 Hamburg · Germany

eppendorf@eppendorf.com · www.eppendorf.com