
Performance and Physical Data

Model No.	280A
Catalog No.	13-262-280A (120V), 13-262-281A (240V)
Vacuum Range	Atmosphere to 30" Hg.* with adequate pump
Vacuum Leak Rate	< .2" Hg per 24 hours
Temperature Range	Ambient to 200°C
Shelf Capacity	2.0 ft ²
Inside Dimensions	11.125"D x 9.75"W x 9.75"H
Outside Dimensions	16"D x 15"W x 20"H
Standard Electrical Service	1000 Watts, 120 VAC 50/60 Hz. 8.4 Amps Optional: 1000 Watts, 240 VAC, 50/60 Hz , 4.2 Amps
Net Weight (Pounds)	65 lb.

**Within 10 Microns*