

FIXED SPEED VORTEXER MIXER SPECIFICATIONS

Overall dimensions (L x W x H): 8 x 4.25 x 5"
(20.3 x 10.8 x 12.7cm)

Electrical (50/60 Hz): 120V, 1.2 amps, 150 watts
230V, 0.6 amps, 150 watts

Fuses: 5mm x 20mm, 5 amp quick acting

Speed range*: 3200rpm **

Orbit: 0.194" (4.9mm)

Duty rating: intermittent duty

Controls: none

Ship weight: 11.7lbs (5.3kg)

ANALOG VORTEXER MIXER SPECIFICATIONS

Overall dimensions (L x W x H): 8 x 4.25 x 5"
(20.3 x 10.8 x 12.7cm)

Electrical (50/60 Hz): 120V, 1.2 amps, 150 watts
230V, 0.6 amps, 150 watts

Fuses: 5mm x 20mm, 5 amp quick acting

Speed range*: 300 to 3200rpm **

Orbit: 0.194" (4.9mm)

Duty rating: intermittent duty

Controls: auto/off/on rocker switch,
speed knob: variable 1 to 10 dial marks

Ship weight: 11.7lbs (5.3kg)

DIGITAL & PULSING VORTEXER MIXER SPECIFICATIONS

Overall dimensions (L x W x H): 8 x 4.25 x 5"
(20.3 x 10.8 x 12.7cm)

Electrical (50/60 Hz): 120V, 1.2 amps, 150 watts
230V, 0.6 amps, 150 watts

Fuses: 5mm x 20mm, 5 amp quick acting

Speed range*, on mode: 500 to 3000rpm **
auto mode: 1000 to 3000rpm **
pulse mode: 3000rpm **

Timer: 1 second to 9999 minutes
(increased in 1 second increments)

Orbit, Digital unit: 0.194" (4.9mm)
Pulsing unit: 0.098" (2.5mm)

Duty rating: intermittent duty

Controls: see page 5

Ship weight: 11.7lbs (5.3kg)

* Maximum speed will vary depending on accessory being used.

** 230V units maximum speed: 2500rpm