

Thermo Scientific General Purpose Rotors

F20-6x100 Carbon Fiber Rotor

F15-6x100 Carbon Fiber Fixed Angle Rotor

| Rotor Specifications | | |
|----------------------|-----------|-------------|
| Capacity (ml) | 6 x 100 | |
| Tube Dimensions (mm) | 38 x 105 | |
| Angle ° | 25 | |
| Maximum Speed (rpm) | 15.000 | |
| K Factor | 1536 | |
| Net Weight (kg) | 4,5 | |
| RCF/Radius | RCF (x g) | Radius (cm) |
| Maximum | 24.652 | 9,8 |
| Minimum | 6.283 | 2,5 |



| F20-6x100 | |
|-----------|------------------------------|
| Cat. No. | Description |
| 75003698 | F20-6x100 Carbon Fiber Rotor |

| Accessories | |
|-------------|--------------------------|
| Cat. No. | Description |
| 099-069031 | Replacement Lid Assembly |
| 021-069031 | O-ring Refresher Kit |

| Tube Ordering Information | | | | | | | | | | | | | | |
|---------------------------|----------------|----------------|-----------------------|------|------------------|---------------|----------------------------|----------|-------------|--------------------|----------|---------------|-----------------------|---|
| Cat. No. | Tube Vol. (ml) | Vol. Fill (ml) | Description | Qty. | Max. Speed (rpm) | Dims ØxL (mm) | Sealing Assembly Required: | | | Adapters Required: | | | Accessories Required: | |
| | | | | | | | Cat. No. | Qty /Pkg | Description | Cat. No. | Qty /Pkg | No. of Places | | |
| 3118-0085 | 85 | 70 | PC Oak Ridge Tube | 10 | 15.000 | 38x105 | Included | 10 | PP Sealing | - | - | - | - | - |
| 334959 | 50 | 50 | Nunc Conical | 25 | - | 30x121 | Included | 25 | PP Sealing | 75003103 | 1 | 1 | - | - |
| 3139-0050 | 50 | 50 | PP Oak Ridge Tube | 10 | 15.000 | 29x114 | Included | 10 | PP Sealing | 75003102 | 1 | 1 | - | - |
| 3138-0050 | 50 | 50 | PC Oak Ridge Tube | 10 | 15.000 | 29x114 | Included | 10 | PP Sealing | 75003102 | 1 | 1 | - | - |
| 3114-0050 | 50 | 50 | Teflon Oak Ridge Tube | 2 | 15.000 | 29x114 | Included | 2 | PP Sealing | 75003102 | 1 | 1 | - | - |
| - | 30 | - | Glass Tube | - | - | 24x105 | - | - | - | 75003094 | 1 | 1 | - | - |
| 366036 | 15 | 15 | Nunc Conical | 50 | - | 17x120 | Included | 50 | PP Sealing | 75003095 | 1 | 1 | - | - |
| - | 10 | - | Blood Collection | - | - | 16X100 | - | - | - | 75002906 | 1 | 1 | - | - |
| 3139-0010 | 10 | 8 | PP Oak Ridge Tube | 10 | 15.000 | 16x82 | Included | 10 | PP Sealing | 75003093 | 1 | 2 | - | - |
| 3138-0010 | 10 | 8 | PC Oak Ridge Tube | 10 | 15.000 | 16x82 | Included | 10 | PP Sealing | 75003093 | 1 | 2 | - | - |
| - | 7 | - | Blood Collection | - | - | 13x100 | - | - | - | 75003092 | 1 | 2 | - | - |
| - | 5 | - | Blood Collection | - | - | 13x75 | - | - | - | 75003092 | 1 | 2 | - | - |
| - | 5 | - | RIA Tubes | - | - | 13x75 | - | - | - | 75003092 | 1 | 2 | - | - |
| - | 3,5 | - | Round Bottom Tube | - | - | 11x100 | - | - | - | 75003091 | 1 | 4 | - | - |
| - | 1,5/2 | - | Microtube | 1 | - | 11 x 42 | - | - | - | 75002905 | 1 | 4 | - | - |