

Thermo Scientific General Purpose Rotors

F15-8x50c Carbon Fiber Rotor

F15-8x50c Fixed Angle Carbon Fiber Rotor

Rotor Specifications	
Capacity (ml)	8 x 50 conical
Tube Dimensions (mm)	30 x 121
Angle °	25
Maximum Speed (rpm)	14,500
K Factor	1,062
Net Weight (kg)	4.4
RCF/Radius	RCF (x g)
Maximum	24,446
Minimum	10,098
RCF/Radius	Radius (cm)
Maximum	10.4
Minimum	4.3



F15-8x50c Carbon Fiber Rotor

Cat. No.	Description
75003663	F15-8x50c Carbon Fiber Rotor

Accessories

Cat. No.	Description
099-085077	Replacement Lid Assembly
021-085077	O-ring Refresher Kit

Tube Ordering Information

Cat. No.	Tube Vol. (ml)	Vol. Fill (ml)	Description	Qty.	Max. Speed (rpm)	Dims ØxL (mm)	Sealing Assembly Required:			Adapters Required:			Accessories Required:		
							Cat. No.	Qty /Pkg	Description	Cat. No.	Qty /Pkg	No. of Places	Cat. No.	Qty /Pkg	No. of Places
334959	50	50	Nunc Conical	25	14,500	30x115	Included	25	PP Sealing	-	-	-	-	-	-
366036	15	15	Nunc Conical	50	-	17x120	Included	50	PP Sealing	75100378	2	1	-	-	-