

## **Specifications**

### **Electrical Requirements**

Model 625	Cat. No. 13-247-625G	120 V, 50/60 Hz
	Cat. No. 13-247-626G	240 V, 50/60 Hz
Model 637	Cat. No. 13-247-637G	120 V, 50/60 Hz
	Cat. No. 13-247-638G	240 V, 50/60 Hz
Model 650	Cat. No. 13-247-650G	120 V, 50/60 Hz
	Cat. No. 13-247-651G	240 V, 50/60 Hz

### **Power Requirements**

Model 625	1300 W
Model 637	1300 W
Model 650	1300 W

### **Chamber Volumes**

Model 625	2.5 cu ft.
Model 637	3.8 cu ft.
Model 650	5.0 cu ft.

### **Chamber Dimensions (W x D x H)**

Model 625	18 x 18 x 13.5 in
Model 637	18 x 18 x 20 in
Model 650	18 x 18 x 26.5 in

---

---