

<i>Table 2:</i> (differences of 230V instruments are shown in parentheses)		
Rotor designation	Micro Liter Rotor 24 x 2.0 ml	
Order no.	7500 3337	
Maximum permissible load [g]	24 x 4	
Maximum speed n_{max} [rpm]	15 000	
Maximum RCF value at n_{max}	21 885	
Radius (max./min.) [cm]	8.7 / 5.9	
Tube angle [°]	45	
Acceleration / deceleration time [s]	30 / 35	
Min temperature at n_{max} [°C] relative to room temperature 23°C	< 0	
Speed at 4°C [rpm]	15 000	
Aerosol-tight *	yes	
Autoclavable	121°C	

* Checked by CAMR, Porton-Down, UK