

FIXED SPEED VORTEXER MIXER SPECIFICATIONS

Overall dimensions (L x W x H):	8 x 4.75 x 5.3" (20.3 x 12.1 x 13.5cm)
Electrical (50/60 Hz):	120V, 1.2 amps, 135 watts 230V, 0.6 amps, 135 watts
Speed range*:	3200rpm **
Orbit:	0.194" (4.9mm)
Duty rating:	intermittent duty
Controls:	none
Ship weight:	10lbs (4.5kg)

ANALOG VORTEXER MIXER SPECIFICATIONS

Overall dimensions (L x W x H):	8 x 4.75 x 5.3" (20.3 x 12.1 x 13.5cm)
Electrical (50/60 Hz):	120V, 1.2 amps, 135 watts 230V, 0.6 amps, 135 watts
Speed range*:	300 to 3200rpm **
Orbit:	0.194" (4.9mm)
Duty rating:	intermittent duty
Controls:	touch/off/on rocker switch, speed knob: variable 1 to 10 dial marks
Ship weight:	10lbs (4.5kg)

DIGITAL & PULSING VORTEXER MIXER SPECIFICATIONS

Overall dimensions (L x W x H):	8 x 4.75 x 5.3" (20.3 x 12.1 x 13.5cm)
Electrical (50/60 Hz):	120V, 1.2 amps, 135 watts 230V, 0.6 amps, 135 watts
Speed range*:	On mode: 500 to 3000rpm ** Touch mode: 1000 to 3000rpm ** Pulse mode: 3000rpm **
Timer:	1 second to 9999 minutes (increased in 1 second increments)
Orbit:	Digital unit: 0.194" (4.9mm) Pulsing unit: 0.098" (2.5mm)
Duty rating:	intermittent duty
Controls:	see page 5
Ship weight:	10lbs (4.5kg)

* *Maximum speed will vary depending on accessory being used.*

** *230V units maximum speed: 2500rpm*