

Section 3 Specifications

These are nominal specifications. Fisher Scientific reserves the right to change specifications or designs at any time without incurring obligation.

Temperature range
 Cover Open: Ambient to 75°C (model 102 to 85°C, model 102S to 70°C)
 Cover closed: Ambient to 100°C
Uniformity ±0.2°C at 37°C
 ±0.5°C at 37°C (shallow pan)
Stability ±0.2°C at 37°C (with cover)

Safety Specifications

Indoor use only
Altitude up to 2000 meters
Ambient Temperature 0°C to 40°C
Ambient Humidity 0 - 95% non-condensing
Voltage Requirements 120VAC or 230VAC nominal
Voltage Fluctuations ±10% of nominal voltage
Installation Category II
Pollution Degree 2

Certification

Declaration of Conformity is available from the factory

Section 3
Specifications

Catalog Number	Model Number	Specific Information	Voltage
15-460-2Q (2228)	102	2L	120V
15-460-3Q (2235)		2L	230V
15-460-2SQ (2230)		2L	120V
15-460-3SQ (2236)		2L	230V
15-460-5Q (2239)	105	5L	120V
15-460-6Q (2240)		5L	230V
15-460-10Q (2223)	110	10L	120V
15-460-11Q (2224)		10L	230V
15-460-15Q (2226)	115	5L/10L, dual chamber	120V
15-460-16Q (2227)		5L/10L, dual chamber	230V
15-460-20Q (2231)	120	20L	120V
15-460-21Q (2232)		20L	230V
15-460-28Q (2233)	128	28L	120V
15-460-29Q (2234)		28L	230V